

 Manufactured for:
Dynarex Corporation
10 Glenshaw Street
Orangeburg, NY 10962
USA · www.dynarex.com
Made in China

Reorder No. 30401

dynarex *Resp-O₂*™

(HEPA) Filter

Bacterial: 99.9999% • Viral: 99.999%
CO₂ Port



CAUTION: Federal (USA)
Law restricts this device to sale
by or on the order of a physician



SYMBOL GLOSSARY
For an explanation of symbols
used in Dynarex packaging, visit
dynarex.com/symbols.php

DESCRIPTION

The HEPA filter is a disposable, single-patient use device indicated for patient during the delivery of ventilator gases. The HEPA filter provides filtration for reducing possible cross-contamination between patient and equipment. The HEPA filter is ideal for use in-hospital, ICU, anesthesia, respiratory therapy, during transport and with resuscitators.

INSTRUCTIONS FOR USE

- Read these instructions carefully before using the product.
- Place the HEPA filter between the breathing circuit and the equipment being used.
- Always replace the HEPA Filter after each patient.
- When used continuously on a single patient, change the HEPA Filter every 24 hours or more frequently as required.

PRECAUTIONS

- All tubing and connections to the HEPA Filter shall be properly attached and checked for leakage prior to use.
- Before use, verify that the HEPA Filter has no occlusions and that air will flow through it.
- Compensation of ventilation may be necessary when using the HEPA Filter since dead space will be added to the system.
- During the use of the HEPA Filter, the patient shall be closely monitored and proper airway care administered if complications arise.
- If a sampling tube is not connected, make sure that the sampling port cap on the HEPA Filter is properly secured.
- The HEPA Filter must be changed between patients.

WARNINGS

- The HEPA Filter is designed for single-use only and should not be cleaned and reused.

TECHNICAL DATA

Tidal Volume Range	250-1500 mL
Resistance	at flow 30 L/min 1.2 cm H ₂ O at flow 60 L/min 2.6 cm H ₂ O
Dead Space	68 mL
Filtration Efficiency	Bacterial: 99.9999% Viral: 99.999%
Weight	39 g
Connections	22M/15F-22F/15M

Note: Measured by the method of Nelson Laboratories Inc. USA.

SYMBOLS

SYMBOL	SYMBOL NOTES
	Consult instructions for use
	Cautions and Warnings
	Keep dry
	Keep away from heat and sunlight
	Do not re-use
	Do not use if package is damaged
	Not made with natural rubber latex
	Not Manufactured with di(2-ethylhexyl)phthalate
	Batch code
	Prescription use only