
LAB CHOICE

ELECTRONIC PIPETTES PIPETAS ELECTRÓNICAS

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USER'S MANUAL
MANUAL DEL USUARIO

PRODUCT OVERVIEW

Dynarex Electronic Pipettes are designed to deliver accurate and precise performance in liquid handling. Its design and operation are based on the principle of air displacement and the use of disposable pipette tips. Frequency, Voltage: 5W/4.2V

These products have been tested according to ISO8655 quality management standards. In compliance with the quality control requirements of ISO8655-6 / DIN 12650, each electronic pipette is fitted with the original tip from the manufacturer and has been tested by gravimetric method with distilled water (DIN / ISO3696, grade at 22°C.)

Dynarex Electronic Pipettes are compatible with most universal tip brands covering full volume range.

CAUTION: Never use the pipette without a tip attached.




CONTENTS

- Electronic Pipette [single-channel (40100 & 40101) or eight-channel (40110, 40111, 40112)]
- AC Adapter
- USB Cable
- Charging Dock
- Calibration Card (ISO8655-6 / DIN12650)
- Silicone Grease
- Filter (40111 & 40112)


CAUTIONS FOR FILTER USE

- If the filter accidentally comes into contact with unintended liquids or solutions, discard immediately or replace.
- Do not use a filter tip with a pipette that is already equipped with a filter.

COMPONENTS AND DISPLAY

Ejector Assembly

GENERAL SET-UP/FUNCTION OVERVIEW


1A. Selecting a Function

1. To select a function, quickly rotate the Function Wheel left and right to activate the function setting.
2. Rotate the Function Wheel to select the desired setting. Settings include:
Aspiration Speed (three speeds)
Dispensing Speed (three speeds)
Sound Alerts (ON/OFF)
Power Off
3. To confirm a setting, press the Parameter Knob once.
4. To exit a function at any time, double-press the Parameter Knob.


1B. Adjustment of Aspiration/Dispensing Speed

1. Rotate Function Wheel to select aspiration or dispensing speed when corresponding symbols start to blink.
2. Rotate Parameter Knob to select among three speeds (1,2,3). [Speed 2 is shown in right illustration.]
3. Press Parameter Knob to confirm.

1C. Turning Sound Alerts On and Off


1. Use the Function Wheel to select the Sound Alert Icon when the corresponding sign starts to blink displaying status as ON or OFF.
2. Rotate Parameter Knob to select between ON and OFF.
3. Press Parameter Knob to confirm.

1D. Powering OFF

1. Rotate Function Wheel to select Power when the corresponding sign starts to blink.
2. Press Parameter Knob to switch the pipette off.


NOTE: The pipette will automatically switch off after idling for 8 minutes.

OPERATION INSTRUCTIONS


2A. Power On

1. Press the Parameter Knob for approximately 2 seconds.


2B. Function Selection

1. Turn the Function Wheel to select the desired function:
 - Pipetting (PI)
 - Stepper/Multi-Dispensing (ST)
 - Diluting (DI) (40100 & 40101 only)


2. After selecting the desired function, press Parameter Knob to confirm.


2C. Using the Pipetting Function: Basic Pipetting (PI)


1. Rotate the Parameter Knob to select "PI" for pipetting function.

2. Status Indicator upward arrow will begin to blink to indicate aspiration mode.


3. "PI" and Volume Range Parameter will blink. You can now set the volume range.


OPERATION INSTRUCTIONS (CONTINUED)

- To increase the volume, rotate the Parameter Knob counter-clockwise. To decrease the volume, rotate the Parameter Knob clockwise. Press Parameter Knob to confirm volume setting.

- Once the volume has been set, the **downward arrow** will begin to blink to indicate dispensing mode. Press the parameter knob to begin dispensing liquid.

2D. Using the Mixing Function

- If mixing is desired while in the pipetting mode, press Parameter Knob and keep holding it down until mixing activity begins. Release the Parameter Knob at any time during the process to stop mixing.

Three possible mixing volume formulas [Pipette volume (10-1000 µL)]

1. Pipetting Volume Greater Than 30%

If the pipetting volume is set to more than 30% of the total volume range of the pipette, 30% of the maximum volume of the pipette will be mixed. (Example: If the set volume is 700 µL, the actual mixing volume will be 300 µL.)

2. Pipetting Volume Less Than 30%

If the pipetting volume is set to less than 30% of the total volume range of the pipette (but more than the pipette's minimum volume range) 50% of the set pipetting volume will be mixed. (Example: If the set volume is 270 µL, the actual mixing volume will be 135 µL.)

3. Minimum Pipetting Volume

If 50% of the set pipetting volume is less than the minimum volume range of the pipette, the minimum volume will be mixed. (Example: If the set volume is 150 μL , the actual mixing volume will be 100 μL .)


2E. Using the Stepper/Multi-Dispensing Function (ST): Dispensing Settings

1. Rotate the Parameter Knob to select "ST" for Stepper/Multi-Dispensing function. The downward arrow will begin to blink.
2. Rotate the Parameter Knob to set the dispensing volume.
3. Press Parameter Knob to confirm.
4. Rotate the Parameter Knob to set the dispensing frequency. (Example: The image shows dispensing frequency is 5 times, and each volume is 200 μL .)

5. Press Parameter Knob to perform liquid dispensing. Each press of the Parameter Knob will dispense a liquid volume according to the set volume.

6. After dispensing the last batch of liquid, press the Parameter Knob to empty any residual liquid.
7. The pipette will then return to aspiration mode.

2F. Using the Diluting Function (DI) (available with 40100 & 40101 only)

A. Volume Setting


1. Rotate the Parameter Knob to set the volume of target A. Press Parameter Knob to confirm.
2. Rotate the Parameter Knob to set the volume of target B. Press Parameter Knob to confirm.

B. Diluting Instructions

1. Press Parameter Knob to aspirate target A.


2. Remove the pipette tip from the liquid and press Parameter Knob to aspirate a bit of air.


OPERATION INSTRUCTIONS (CONTINUED)

3. Press Parameter Knob to aspirate target B.

4. Press Parameter Knob to dispense targets A and B.

STORAGE GUIDELINES

- Keep the electronic pipette on the Charging Dock in a vertical position when not in use. [See page 11 for Charging Dock set-up instructions]
- The electronic pipette should be charged at least once a month, even when it is not in use for extended periods of time. Make sure there is always at least 50% power in the battery.
- During storage periods, the recommended temperature range is from 0-40°C with a humidity of no more than 80%.

CARE AND MAINTENANCE

For optimal functionality and performance, clean the electronic pipette before and after each use.

To clean, spray the outer surface with an appropriate pipette detergent or 75% ethanol. Wipe with lint-free cloth.


Note: Check the performance of the pipette every three months, and after any third-party repair and maintenance.

ANY WARRANTY WILL BE DEEMED AS VOID IF THE FAULT IS FOUND TO HAVE BEEN CAUSED BY UNAUTHORIZED CLEANING.

STERILIZATION INSTRUCTIONS FOR SINGLE-CHANNEL PIPETTES (40100-40101)

CAUTION: This pipette cannot be autoclaved as a whole piece. Any parts or components not identified in the illustration below must not be autoclaved or disassembled.

The autoclavable parts for Single-Channel Pipettes are labeled below.

3A. Disassembly and Cleaning Instructions


1. Rotate the Tip Ejector Collar counter-clockwise to remove it. **NOTE: The Tip Ejector Collar cannot be autoclaved.**


2. Rotate the Piston Cover counter-clockwise to remove the Ejector Spring.

3. Separate the Seal Gasket and O-Ring from the Piston Cover. **NOTE! Keep all parts in a safe place. Do not discard them, as they will be needed for reassembly.**

4. Separate the Piston from the Piston Spring.

5. Clean the disassembled autoclavable parts with 75% ethanol, or sterilize at 121°C, 1 bar over pressure for 20 minutes.
6. Apply an appropriate layer of silicone grease to the O-Ring, Seal Gasket, and Piston.
7. Reassemble.
8. After autoclaving and reassembling, let the pipette rest for 12 hours.


3B. Calibration

1. **Calibration should be performed after each sterilization.**
2. Calibration should take place at 20-25°C, kept constant within $\pm 0.5^\circ\text{C}$.
3. If the pipette does not function properly after calibration, please contact your distributor for assistance.

STERILIZATION INSTRUCTIONS FOR EIGHT-CHANNEL PIPETTES (40110-40112)

CAUTION: This pipette cannot be autoclaved as a whole piece. Any parts or components not identified in the illustration below must not be autoclaved or disassembled.

The autoclavable parts for Eight-Channel Pipettes are labeled below.

4A. Disassembly and Cleaning Instructions


1. Rotate the Connector counter-clockwise and pull out the Ejector Assembly in a downward motion.
2. Remove the Connector and Ejector Assembly shell. **NOTE: THESE PARTS CANNOT BE AUTOCLAVED!**

3. Open the front and back covers of the white Ejector Assembly shell along its middle seam and remove the Ejector Spring.



4. Remove the black Piston Assembly from inside.

5. Unscrew the screws at the four corners of the Piston Assembly.

6. Remove the Piston Rod from inside each Piston. **NOTE! Keep all parts in a safe place. Do not discard them, as they will be needed for reassembly.**


7. Clean the disassembled autoclavable parts with 75% ethanol, or sterilize at 121°C, 1 bar over pressure for 20 minutes.
8. Apply an appropriate layer of silicone grease to the O-Ring, Seal Gasket, and each Piston.
9. Reassemble.

4B. Calibration:


1. Calibration should be performed after each sterilization.
2. Calibration should take place at 20-25°C, kept constant within $\pm 0.5^\circ\text{C}$.
3. If the pipette does not function properly after calibration, please contact your distributor for assistance.

SETTING UP THE CHARGING DOCK

1. There are two strips of adhesive tape on the surface of the Charging Dock. Remove the backing paper from the adhesive you want to use and affix the Charging Dock onto a clean, flat surface.

2. Insert the Charging Port into the Charging Dock.
3. Insert USB Cable into USB Charging Port.

TROUBLESHOOTING GUIDE

Problem	Possible Cause	Solution
Remaining droplets	Incompatible tip	Use original tip.
	Plastic self-lubrication of the tip is not even.	
Leakage or inadequate aspiration	Tip is loose.	Attach the tip firmly.
	Foreign objects between tip and connector	Clean the Tip Cone and replace the tip.
	Pipette is contaminated.	Clean and grease O-Ring, Piston, and clean Tip Cone.
	Inadequate application of silicone grease to Piston or O-Ring	Apply silicone grease.
	O-Ring and Piston are unlocked or broken.	Replace O-Ring.
	Inappropriate operation	Follow the instructions in this User's Manual.
	Calibration is needed or liquid density is significantly different from that of water.	Re-calibrate according to the instructions in this User's Manual.
Pipette is broken.	Send the pipette for service.	
Inadequate aspiration due to jammed pipette	Liquid has leaked into Tip Connector and then dried, causing a blockage.	Clean and grease O-Ring and Piston; clean Tip Cone.
Ejector is jammed or cannot work properly.	Tip Connector and Ejector are contaminated.	Clean Tip Connector and Ejector.

LIMITED WARRANTY

Your Dynarex Product is warrantied to be free of defects in materials and workmanship for twelve (12) months from original date of purchase for the original purchaser. This item was built to exacting standards and carefully inspected prior to shipment. This Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

The Warranty shall not apply under the following conditions:

- Problems arising from normal wear
- Problems arising from failure to adhere to the product instructions
- Problems arising from misuse, negligence, accident or improper operation, maintenance or storage
- Problems arising from modifications or unauthorized repairs, parts or attachments
- Products where the serial number has been removed or defaced
- Problems with non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement

Dynarex shall not be liable for any consequential or incidental damages whatsoever. Dynarex shall repair or replace defective products at its option. The foregoing warranty is exclusive and in lieu of other express warranties, if any, including the implied warranties of merchantability and fitness of a particular purpose. The remedy for any violation of the implied warranty shall be limited to repair or replacement of the defective product pursuant to the terms contained herein.

If you have a question about your Dynarex device or this warranty, please contact an authorized Dynarex dealer.

Each and every electronic pipette has been calibrated and tested in compliance with ISO8655-6/DIN12650 (calibration card included in the contents of delivery package) when manufactured.

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Manufactured for/ Fabricado para: Dynarex Corporation • 11 Dynarex Drive • Middletown, NY 10941 • USA • www.dynarex.com