

 Manufactured for:
Dynarex Corporation
10 Glenshaw Street
Orangeburg, NY 10962
USA · www.dynarex.com
Made in Turkey

Reorder No. 30309



Heat Moisture Exchanger (HME) Tracheostomy

with Side Oxygen Port



CAUTION: Federal (USA)
Law restricts this device to sale
by or on the order of a physician



SYMBOL GLOSSARY
For an explanation of symbols
used in Dynarex packaging, visit
dynarex.com/symbols.php

DESCRIPTION

The HME is disposable single-use device indicated for tracheostomy patients requiring humidification. The HME is for use in hospital, respiratory therapy, and homecare. The HME can be used on adult and pediatric patients. The HME is indicated for use only under the supervision of qualified medical personnel. The HME is designed for disposable use and should be changed at least every 24 hours.

INSTRUCTIONS FOR USE

- Read these instructions carefully before using the product.
- Place the HME onto the tracheostomy device connection piece.
- When used continuously on the patient, change the HME every 24 hours or more frequently as required.

CONTRAINDICATIONS

- The HME is contraindicated in patients producing fulmination, frothy secretions within their airway and lungs.
- The HME shall not be used on patients with tidal volumes of less than 150 mL, or who weigh less than 33 lb. (15 kg).

PRECAUTIONS

- The HME shall be properly attached and checked for leakage prior to use.
- Before use, verify that the HME has no occlusions and that air will flow through it.
- During the use of the HME, the patient shall be monitored and proper airway care administered if complications arise.

WARNINGS

- HME is designed for single use only and should not be cleaned and reused.

TECHNICAL DATA

Tidal Volume Range	300-1500 mL
Resistance	at flow 30 L/min 0.8 cm H ₂ O at flow 60 L/min 1.8 cm H ₂ O
Dead Space	15 mL
Moisture Output (per ISO 9360 test procedure)	Tidal Volume = 300 mL 32.7 mg H ₂ O/L 500 mL 28.6 mg H ₂ O/L 1500 mL 25 mg H ₂ O/L
Weight	8.4 g
Connections	15F

SYMBOLS

SYMBOL	SYMBOL NOTES
	Consult instructions for use
	Cautions and Warnings
	Keep dry
	Keep away from heat and sunlight
	Do not re-use
	Do not use if package is damaged
	Not made with natural rubber latex
	Not Manufactured with di(2-ethylhexyl)phthalate
	Batch code
	Prescription use only