

HILLYARD

Ready-To-Use

CITRUS DISINFECTANT

ACTIVE INGREDIENT:

Citric Acid.....0.60%

OTHER INGREDIENTS:99.40%

TOTAL:100.00%

CAUTION:

KEEP OUT OF REACH OF CHILDREN

FIRST AID: IF IN EYES: Hold eyelids open and flush with a steady, gentle stream of water for 15 minutes. Remove contact lenses, if present, after the first five minutes, then continue rinsing. Call a Poison Control Center or doctor for further treatment advice.

Follow the directions for use on the pesticide label when applying this product. Further sale or distribution of this product is strictly prohibited.

Customer Affixed Label

EPA Reg. No. 34810-35-1658



HILLYARD INDUSTRIES, INC.
St. Joseph, MO 64502 U.S.A.
Telephone: 816-233-1321

In the presence of 5% blood serum and a contact time of 5 minutes, is effective against:

BACTERICIDAL: Pseudomonas aeruginosa, Staphylococcus aureus, Salmonella enterica, Listeria monocytogenes, Escherichia Coli 0157:H7, Staphylococcus aureus- MRSA, Staphylococcus epidermidis - MRSE, Enterococcus faecalis Vancomycin Resistant (VRE), Carbapenem-resistant Klebsiella pneumoniae.

VIRUCIDAL*: Respiratory Syncytial Virus, Rotavirus, Human Immunodeficiency Virus Type 1 (Strain HTLV-1118), Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2, Vaccinia Virus, Influenza A Virus (Strain Hong Kong), Adenovirus Type 2, Feline Calicivirus (surrogate for Norovirus). Contact time of 3 minutes is effective against: Hepatitis B Virus, Hepatitis C Virus, Rhinovirus.

FUNGICIDAL: Trichophyton interdigitale, Candida albicans.

TUBERCULOCIDAL: Effective against Mycobacterium bovis in 5 minutes at 20°C/68°F on previously cleaned hard non-porous inanimate surfaces.

SANITIZER: Staphylococcus aureus and Escherichia coli in 60 seconds.

PRECAUTIONARY STATEMENTS – HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Product is ready to use. DO NOT DILUTE. Preclean areas with visible soil prior to disinfection.

TO CLEAN/DEODORIZE: Apply undiluted to hard, non-porous surfaces by spraying, wiping with a cloth or towel until clean.

TO SANITIZE: Food Contact: Preclean. Wet the surface. Let stand for 60 seconds on hard, non-porous surfaces.

Wipe or let air dry. No rinsing or wiping is required, even on food contact surfaces.

TO DISINFECT: Apply Citrus Disinfectant undiluted to hard, non-porous surfaces with a sprayer, cloth, mop or brush, or by immersing equipment to be disinfected. Wet all surfaces thoroughly. The surface must remain visibly wet for 5 minutes. Wipe or let air dry. No rinsing or wiping is required.

FUNGICIDAL DIRECTIONS: This product is effective against Trichophyton interdigitale (athlete's foot fungus and a cause of ringworm) on inanimate surfaces in locker rooms, dressing rooms, exercise facilities, and shower and bath areas. Follow disinfection directions.

AREAS OF USE: Hospitals, Kitchens, Bathrooms, Gyms, Health Care Facilities, Medical Offices, Dental Offices, Day Care Centers, Nurseries, Schools, Nursing Homes. **TYPES OF SURFACES:** On washable hard non-porous surfaces of floors, walls, fixtures, furniture, metal, stainless steel, chrome, CuVerro® antimicrobial copper, glazed porcelain and glazed ceramic tile, vinyl, plastic, baked enamel, and glass. This includes infant care equipment¹, kitchen surfaces, wrestling mats¹, Exercise equipment¹, radiology equipment surfaces, computer surfaces and keyboards.

DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/ HIGH LEVEL DISINFECTION:

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or, (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Citrus Disinfectant KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES / OBJECTS PREVIOUSLY SOILED WITH BLOOD / BODY

FLUIDS in all settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS virus).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES / OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1.

PERSONAL PROTECTION: When handling objects soiled with blood or body fluids use disposable latex gloves, gowns, mask and eye coverings.

CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

CONTACT TIME: DO NOT DILUTE. Leave precleaned surfaces wet for five minutes using mop, cloth or sponge to mitigate HIV-1 virus. Use a 5-minute contact time to kill all bacteria, viruses* and fungi listed on the label.

DISPOSAL OF INFECTIOUS MATERIAL: Blood and other fluids should be autoclaved and disposed of according to Federal, State and Local regulations for infectious waste disposal.

STORAGE AND DISPOSAL: Store airtight at room temperature. Nonrefillable container. Do not reuse or refill this empty container. Offer empty container for recycling. If recycling is not available, discard in trash.

*Viruses

¹on hard, non-porous surfaces