

Hemostat Gauze Dressings

Soluble Regenerated Cellulose
For Bleeding Control



Instructions for Use

Description

Hemostat Gauze Dressing is a white or light yellow sterile fabric prepared by the chemical treatment of regenerated cellulose. It is an effective and innovative soluble hemostat widely used in first-aid and surgery. When the fabric comes in contact with blood or exudates, it can quickly transform into a clear viscous gel, which fills wound voids, seals capillary ends, activates the clotting system and thereby helps to stop bleeding.

Indications

Hemostat Gauze Dressing is indicated for the management of topical wounds and to temporarily control external surface bleeding, external topical bleeding of minor cuts and lacerations.

For External Use Only.

Directions

1. Clean the wound.
2. Apply to cover the bleeding site. Dressings can be cut or folded as needed.
3. Gently press the wound site until bleeding stops.
4. Moisten with clean water after the bleeding stops.
5. Excess gel can be rinsed away with clean water.
6. Repeat the steps above with additional dressings as needed. If bleeding continues, contact a doctor or emergency room.

Warnings

In case of wound swelling or redness for more than 48 hours after use, contact a doctor. For external topical and temporary use only.

Cautions

Hemostat Gauze Dressing is supplied sterile. Do not use if package is open or damaged. Hemostat Gauze Dressing should not be used on dirty, contaminated or infected wounds. Consult your doctor if you notice any itching, burning, rash or signs of infection such as fever, pain, redness or swelling.

Specifications

Reorder No.	Size
3686	2 in. × 4 in. (5 cm × 10 cm)
3687	4 in. × 4 in. (10 cm × 10 cm)

Storage

Keep in a cool dry place for maximum shelf life. The preferred storage temperature should be between 41°F and 104°F or 5°C and 40°C.

Hemostat Gauze Dressing is sterilized in the original packaging by gamma radiation.

Shelf Life

Three years.

Manufactured for:
Dynarex Corporation
10 Glenshaw Street
Orangeburg, NY 10962
USA • www.dynarex.com



SYMBOL GLOSSARY
For an explanation of symbols used in Dynarex packaging, visit dynarex.com/symbols.php

STERILE

