



## **Evenflo® QuickAction™ Electric Nasal Aspirator**

Read all instructions before  
using this product.

**KEEP INSTRUCTIONS FOR  
FUTURE USE.**



## Table of Contents

Introduction.....	2
Warnings .....	3
Important Information Before Use .....	3
Product Identification .....	4
Recharging the Battery.....	4
Instructions for Use.....	5
Cleaning and Disinfecting.....	6
Technical Specifications.....	7
Troubleshooting.....	7
EMC Tables.....	8
FCC Interference Statement .....	11
Limited Warranty .....	12

## Introduction

1. Only use the **Evenflo QuickAction™ Nasal Aspirator** to remove nasal secretions.
2. **Soft Silicone Tip Design** to reduce possibility of injury during operation. Safe and comfortable for children.
3. **Replaceable Tip:** Bulbous and Cone-shaped tips available.
4. **Compact and Lightweight.**
5. **Rechargeable Lithium Ion Battery**
6. **Rechargeable** via USB-C Cable.
7. **Easy Clean.** Silicone tips and collection cup can be removed and washed or boiled to disinfect.
8. **Quiet.** Low noise for nighttime use.
9. **Soft Amber Night Light / Low Battery Indicator** for use in low light situations. When light is flashing, Nasal Aspirator needs recharging.

**NOTE:** Use the **Evenflo QuickAction Nasal Aspirator** to remove nasal secretions under advice of your pediatrician.

## **WARNING!**

- Product is intended to be used by adults.
- Product is intended to be used **ONLY** by adults.
- This is not a toy. **DO NOT** allow children to play with it.
- Keep away from children and pets. Small parts may be inhaled or swallowed.
- Never immerse Nasal Aspirator main body in water.
- Never push the silicone tip too far into the child's nostril.
- Always wash and dry the tip before and after each use.
- **NEVER** use in the child's mouth or ears.

## **Important Information Before Use**

The **Evenflo QuickAction Nasal Aspirator** is designed to provide intermittent suction for removal of nasal secretions and mucus in children (age 0-12 years old) in a home environment.




### **WARNING!:**

**DO NOT** modify this device. **DO NOT** expose the device to extreme temperatures, very high humidity, or direct sunlight.

1. Avoid dropping the device.
2. Charge battery before first use.
3. Low battery is indicated by the nightlight flashing (see pg. 4).
4. 3 green LED lights display when battery is fully charged.
5. The Evenflo QuickAction Nasal Aspirator is for removing nasal secretion only. Do not use this device other than for the intended purposes.
6. For assistance with setup, use, or maintenance of this device, please contact ParentLink at 1-800-233-5921.
7. Report any serious incidents involving this Nasal Aspirator to Evenflo by calling ParentLink at 1-800-233-5921.

## Product Identification


1. Silicone Tip
2. Collection Cup
3. Night Light & Low Battery Indicator
4. 3-Level Lighted Suction Indicator
5. Power Button
6. Main Body
7. USB-C Charging Port
8. USB-C Cable



## Recharging the Battery



1. To recharge battery pack, use the included USB-C Cable.
2. Do not use a Nasal Aspirator that is damaged or modified.
3. Do not expose a battery pack or Nasal Aspirator to fire or excessive temperature. Fire or temperature above 265°F (130°C) may cause an explosion.
4. Charging improperly at temperatures outside of the specified range may damage the battery and increase the risk of fire. -25°C~+55°C (-13°F~+131°F) with relative humidity up to 85% (non condensing).

## Instructions for Use

1. Choose the correct size included tip based on child's age:
  - **Bulbous-shaped tip (1A)** with the narrow opening: recommended for children 0-2 years.
  - **Cone-shaped tip (1B)** with the wider opening: recommended for children 2+ years.
2. Attach the selected tip firmly to the Nasal Aspirator body; it should fit flush against the main body.
3. Press and release the power button  to turn on. Test the Nasal Aspirator on the back of your hand to ensure it is activated and suctioning. The aspirator will continue to suction for one minute and then will shut off automatically to conserve battery life.
4. If selected to desired suction level, hold the child upright and hold the Nasal Aspirator gently inside the child's nostril to begin suctioning to clear the mucus. Suction power is mild and should only take a few seconds; consult your pediatrician if child is uncomfortable.
5. The Nasal Aspirator has a soft amber night light for low light situations that also acts as the low battery indicator. When the light is flashing, the Nasal Aspirator needs to be recharged.

**NOTE:** Mucus may not immediately be visible in the collection cup. It is recommended that the silicone tip and collection cup be removed and rinsed between each use to keep any mucus from drying and clogging the device.

### Suction Level Selection

1. The Nasal Aspirator is designed with 3 suction levels. The suction level defaults to the last level used.
2. To change suction level, press and release the power button  to turn Nasal Aspirator on. Then press and hold power button  until desired suction level is reached (see pg. 4). Repeat for next suction level.

**NOTE:** 1 light = low, 2 lights = medium, 3 lights = high


## Cleaning and Disinfecting

The following cleaning procedures should be performed daily or when the device is being used regularly. Failure to properly clean your unit may result in clogging and shortened device life.

1. Carefully separate the silicone tip from the collection cup.
2. Unscrew the collection cup from the device main body.
3. Silicone tip and collection cup should be rinsed well with hot water after each use to prevent mucus from drying and clogging unit.
4. **To Clean (at the end of the day / as needed):** Wash the collection cup and silicone tip with mild, non-abrasive soap and water.
5. **Disinfecting (after an illness / as needed):** After cleaning the pieces, rinse the tip and collection cup thoroughly with hot water.
6. Do not wash or immerse the main body of the aspirator. Use a soft, dry cloth to wipe Nasal Aspirator. Never use abrasive cleaning agents or any industrial strength cleaners.
7. Do not expose the device to extreme temperatures, very high humidity, or direct sunlight.
8. Avoid dropping the device.
9. Allow all parts to dry completely before reassembly.

### NOTE:

- Protect the Nasal Aspirator from drops or heavy impact.
- Store in an area away from sunlight, dust and humidity.
- If tips degrade or get damaged, please replace them.



## Technical Specifications

- **Effective Vacuum Pressure:** Up to 65kPa (487 mmHg)
- **Sound Pressure Level:** < 60dB (30 cm away from the device)
- **Power supply:** Rechargeable Lithium Ion Battery
- **Charging:** USB-C Cable
- **Aspiration Tip:** Soft Silicone
- **Operating environment:** 16°C~35°C (60.8°F~95°F) with relative humidity up to 85% (non condensing)
- **Storage/ Transportation environment:** -25°C~+55°C (-13°F~+131°F) with relative humidity up to 85% (non condensing)
- **Weight:** approx. 4.9 oz (140.8g) with battery
- **Product Dimensions with tip:** approx. 7.36 in x 2.16 in x 1.4 in (18.7 cm x 5.5 cm x 3.5 cm)

## Troubleshooting

Problem	Possible Cause	Solution
The suction is weak or there is no suction	Battery not charged	Charge battery
	The O-ring located where the collection cup screws in is missing.	Make sure there is an O-ring where the collection cup screws into the main body.
	The collection cup is not properly screwed in.	Make sure the three guides on each side of the collection cup line up with the guides on the main body.
	The Nasal Aspirator is clogged.	Wash collection cup and silicone tip with mild, non-abrasive soap and water.

## EMC Tables

60311000/NS15 is intended for use in the electromagnetic environment specified below. The customer or the user of 60311000/NS15 must make sure that it is used in such an environment.

Guidance and manufacturer's declaration - Electromagnetic emissions

Phenomenon	Professional healthcare facility environment a)	HOME HEALTHCARE ENVIRONMENT a)
Conducted and radiated RF EMISSIONS	a)	CISPR 11 Group 1 Class B
Harmonic distortion	Not applicable	
Voltage fluctuations and flickering	Not applicable	

a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive acceptance limits of Group 1 Class B (CISPR 11) have been considered and applied. The equipment is suitable for use in the mentioned environments when directly connected to the Public Mains Network.

b) The test is not applicable in this environment unless the ME EQUIPMENT and ME SYSTEM used will be connected to the PUBLIC MAINS NETWORK and the power input is otherwise within the scope of the Basic EMC standard.

Guidance and manufacturer's declaration - Electromagnetic immunity - Enclosure port

Phenomenon	Basic EMC standard or test method	Immunity test levels	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT



## EMC Tables

ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8kV contact ± 2 kV, ±4kV ±, ±8 kV, ±15 kV air	
Radiated RF EM fields	IEC 61000-4-3	a)	10 V/m b) 80MHz - 2.7 GHz 80% AM at 1kHz
Proximity fields from RF wireless	IEC 61000-4-3	COMPLIANT NOTE: Further information about distances to be maintained between portable and mobile RF communications equipment (transmitters) and the 60311000/NS15 can be requested from supplier using the contact information provided in this manual. However, it is advisable to keep the equipment at an adequate distance of, at least, 0.5 m from mobile phones or other RF communications transmitters to minimize possible interference.	
RATED power frequency magnetic fields	IEC 61000-4-8	30 A/m c) 50 Hz or 60 Hz	

a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive IMMUNITY acceptance limits have been considered and applied.

b) Before modulation is applied.

c) This test level assumes a minimum distance of at least 15 cm between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic fields.

## FCC Interference Statement

- This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:  
(1) this device may not cause harmful interference, and  
(2) this device must accept any interference received, including interference that may cause undesired operation.

**WARNING!** Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## Limited Warranty

For a period of 1 year from the original purchase of this Product, Evenflo warrants to the original end user (“Purchaser”) this Product (including any accessories) against defects in material or workmanship. Evenflo’s sole obligation under this express limited warranty shall be, at Evenflo’s option, to repair or replace any Product that is determined to be defective by Evenflo and determined to be covered by this warranty.

Repair or replacement as provided under this warranty is the exclusive remedy of the Purchaser. Proof of purchase in the form of a receipted invoice or bill of sale evidencing that the Product is within the warranty period must be presented to obtain warranty service. This express limited warranty is extended by Evenflo ONLY to the original Purchaser of the Product and is not assignable or transferable to subsequent purchasers or end users of the Product. For warranty service, contact Evenflo’s ParentLink Consumer Resource Center at 1-800-233-5921 or [www.evenflo.com](http://www.evenflo.com).

EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THIS PRODUCT IS HEREBY DISCLAIMED. EVENFLO SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES OF ANY KIND FOR BREACH OF THE EXPRESS LIMITED WARRANTY ON THIS PRODUCT, OR ANY WARRANTY IMPLIED BY OPERATION OF LAW. OTHER THAN THE EXPRESS LIMITED WARRANTY SET FORTH ABOVE, THERE ARE NO OTHER WARRANTIES THAT ACCOMPANY THIS PRODUCT AND ANY ORAL, WRITTEN OR ANY OTHER EXPRESS REPRESENTATION OF ANY KIND IS HEREBY DISCLAIMED.

# ParentLink®

answers. advice. affirmation.

**[www.evenflo.com](http://www.evenflo.com)**

1-800-233-5921 8 am - 5 pm EST

Manufactured in China by AViTA (WuJiang) Co., LTD.  
Distributed in the U.S.A. by Evenflo Company, Inc.

Model No.: 60311000/NS15

Manufacturer: AViTA (WuJiang) Co. LTD.

Tel: +86-512-8207-8288

Add.: No. 858, Jiao Tong Rd., WuJiang Economic  
Development Zone, Suzhou City, Jiangsu Province, P.R.C.  
Made in China