


# DynaSorb™

## Super Absorbent Dressing

### Self-Adherent

#### Description

DynaSorb™ Super Absorbent Dressing with Self-Adherent Contact Layer (Bordered) is a sterile dressing for use on heavily exuding wounds, consisting of breathable and waterproof nonwoven layer, a super absorbent pad layer, a fabric of wicking layer and soft silicone contact layer. The multi-layered construction facilitates fluid management to provide an optimal moist wound environment to help promote faster wound closure. The absorbent pad is designed to provide optimal absorption and to help lock exudates away and help minimize the risk of maceration and damage to peri-wound skin.

#### Indications

DynaSorb™ Super Absorbent Dressing with Self-Adherent Contact Layer (Bordered) is indicated for management of moderate to high exuding wounds, including: leg ulcers, pressure ulcers, diabetic foot ulcers and dehisced surgical wounds.

#### Contraindications

1. Do not use this product on dry or low exuding wounds, eyes, mucous membranes or in wound cavities because the dressing swells during absorption. Do not use on third-degree burns.
2. DynaSorb™ Super Absorbent Dressing with Self-Adherent Contact Layer (Bordered) with silicone layer should not be used on individuals with a known allergy or hypersensitivity to the dressing or its components.

#### Instructions for Use

1. After wound is cleaned and sterilized, the appropriate size dressing is selected to allow for sufficient space to the wound edge.
2. DynaSorb™ Super Absorbent Dressing with Self-Adherent Contact Layer (Bordered) can be used both as a primary dressing or as a secondary dressing. When used as a primary dressing, DynaSorb™ Super Absorbent Dressing with Self-Adherent Contact Layer (Bordered) with silicone layer must be in close contact with the wound bed to ensure wound fluid is absorbed. Dynarex Super Absorbent Bordered Dressing with silicone layer is self adherent but may require additional support with a mild-to-moderate compression bandaging.
3. The frequency of changing dressing depends on the amount of exudates.
4. Dressing removal can be facilitated by moistening the dressing with saline solution.
5. Dispose of as normal waste. If contaminated, follow local regulations.

#### Warnings

Do not cut the dressing, the cover must remain intact.

#### Shelf Life

Dressings expires 2 years from date of manufacture.

#### Specifications

Reorder No.	Size	Dressing Core
3089	4 in. × 4 in. (10 cm × 10 cm)	3 in. × 3 in. (7.5 cm × 7.5 cm)
3090	6 in. × 6 in. (15 cm × 15 cm)	4.72 in. × 4.72 in. (12 cm × 12 cm)

Manufactured for:  
**Dynarex Corporation**  
 10 Glenshaw Street  
 Orangeburg, NY 10962  
 USA • [www.dynarex.com](http://www.dynarex.com)  
 Made in China



**SYMBOL GLOSSARY**  
 For an explanation of symbols  
 used in Dynarex packaging, visit  
[dynarex.com/symbols.php](http://dynarex.com/symbols.php)