



Directions for Use

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Open lid. Pull out wipe from center of roll and through lid. Close lid to retain moisture. **TO CLEAN:** Unfold premoistened cloth and wipe area to be cleaned for 30 seconds or until clean. **TO DISINFECT:** Product is ready to use. Preclean areas with gross filth and visible soil prior to disinfection. Wipe hard non-porous area or item to be disinfected. Wet all surfaces thoroughly. The surface must remain visibly wet for at least 5 minutes. Wipe or let air dry. No rinsing or wiping is required except on food contact surfaces which require a potable water rinse after use. **FUNGICIDAL DIRECTIONS:** This product is effective against *Trichophyton interdigitale* (athlete's foot fungus and a cause of ringworm) on inanimate surfaces. Follow disinfection directions. **AREAS OF USE:** On washable hard, non-porous surfaces of Hospitals, Kitchens, Bathrooms, Health Care Facilities, Medical Offices, Day Care Centers, Schools, Nursing Homes, and Grocery Stores. **TYPES OF SURFACES:** On hard non-porous surfaces of floors, walls, fixtures, metal, stainless steel, glazed porcelain and glazed ceramic tile, vinyl, plastic, and glass, infant care equipment, kitchen surfaces, wrestling mats, computer surfaces, keyboards, and PPE (Personal Protective Equipment). For plastic and painted surfaces, spot test on an inconspicuous area before use. Do not use on unpainted wood or natural marble or brass. **DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/HIGH LEVEL DISINFECTION:** This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or, (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical

Disinfection Through The Power of Citrus. An alternative to quat, alcohol, or bleach-based disinfectants. EPA Category IV (lowest level of toxicity) requiring no precautionary statements or PPE - will not irritate skin, eyes or lungs. Highly biodegradable ingredients, leaving no harmful or irritating residues.

Features & Benefits

Powerful Disinfection - Broad Efficacy

Ideal For Hard, Non-Porous Surfaces - Ready-To-Use

One-Step, No-Rinse Cleaning & Disinfection

Item Number & Unit of Measure

Item Number	HIL01173
Unit of Measure	160 Wipes/Container

Specifications

Color	White
Fragrance	Citrus
Formula Type	Wipes
PH	2 - 2.4
EPA Registration #	34810-36-1658

Safety

See safety datasheet (if applicable) and product label for safety information, handling and proper use.

devices prior to sterilization or high-level disinfection. ***KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in all settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS virus). **SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS VIRUS) PERSONAL PROTECTION:** When handling objects soiled with blood or body fluids use disposable latex gloves, gowns, mask and eye coverings. **CLEANING PROCEDURE:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant. **CONTACT TIME: Do Not Dilute.** Precleaned surfaces must remain visibly wet for one minute using mop, cloth or sponge to mitigate HIV-1 virus. Use a 5-minute contact time to kill all bacteria, viruses* and fungi listed on the label. **DISPOSAL OF INFECTIOUS MATERIAL:** Blood and other fluids should be autoclaved and disposed of according to Federal, State and Local regulations for infectious waste disposal.

Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2, Vaccinia Virus. **FUNGICIDAL: 2-minute contact time:** Trichophyton interdigitale. **5-minute contact time:** Candida albicans. **TUBERCULOCIDAL: 5-minute contact time:** Mycobacterium bovis

Flammability RTU	0
Health RTU	0
Reactivity RTU	0

Efficacy

In the presence of 5% blood serum, is effective against: **BACTERICIDAL: 2-minute contact time:** Pseudomonas aeruginosa, Salmonella enterica, Staphylococcus aureus. **5-minute contact time:** Carbapenem-resistant Klebsiella pneumoniae, Enterococcus faecalis Vancomycin Resistant (VRE), Staphylococcus aureus - MRSA, Escherichia coli O157:H7, Listeria monocytogenes, Staphylococcus epidermidis - MRSE. **VIRUCIDAL*: 30-second contact time:** SARS-CoV-2 (cause of COVID-19) (Strain USA-WA 1/2020). **3-minute contact time:** Rhinovirus, Hepatitis B Virus (HBV), Hepatitis C Virus (HCV). **5-minute contact time:** Human Immunodeficiency Virus Type 1 (HIV-1), Influenza A Virus, Respiratory Syncytial Virus (RSV), Rotavirus, Feline Calicivirus (surrogate for Norovirus), Adenovirus Type 2,