

D-Lactate Assay Kit (Colorimetric)

LS-K236-100 (100 Tests) • Store at -20°C



Introduction

Lactate is generated by lactate dehydrogenase (LDH) under hypoxic or anaerobic conditions. Monitoring lactate levels is, therefore, a good indicator of the balance between tissue oxygen demand and utilization and is useful when studying cellular and animal physiology. D-lactate is produced in only minor quantities in animals and measuring for D-lactate in animal samples is a means to determine the presence of bacterial infection.

Simple, direct and automation-ready procedures for measuring lactate concentration are very desirable. LSPBio's lactate assay kit is based on lactate dehydrogenase catalyzed oxidation of lactate, in which the formed NADH reduces a formazan (MTT) Reagent. The intensity of the product color, measured at 565 nm, is proportionate to the lactate concentration in the sample.

Key Features

- Sensitive and accurate. Detection limit of 0.05 mM and linearity up to 2 mM D-lactate in 96-well plate assay. For cell culture samples containing phenol red: detection limit of 0.1 mM and linearity up to 1 mM D-lactate in 96-well plate assay.
- Convenient. The procedure involves adding a single working reagent, and reading the optical density at time zero and at 20 min. Room temperature assay. No 37°C heater is needed.
- High-throughput. Can be readily automated as a high-throughput 96-well plate assay for thousands of samples per day.

Applications

- Direct Assays: D-lactate in serum, plasma, and cell media samples.

Components

Component	K236-100
	100 Tests
Assay Buffer	10 mL
Enzyme A	120 µL
Enzyme B	120 µL
NAD Solution	1 mL
MTT Solution	1.5 mL
Standard (20 mM D-lactate)	1.0 mL

Materials Not Supplied

Pipetting (multi-channel) devices. Clear-bottom 96-well plates (e.g. Corning Costar) and plate reader.

Storage

The kit is shipped on ice. Store all kit components at -20 °C.

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121
www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com

D-Lactate Assay Kit (Colorimetric)

LS-K236-100 (100 Tests) • Store at -20°C



Assay Procedure

Important: this assay is based on an enzyme-catalyzed kinetic reaction. Addition of Working Reagent should be quick and mixing should be brief but thorough. Use of a multi-channel pipette is recommended.

1. Standard Curve. Prepare 1000 μL 2.0 mM D-lactate Premix by mixing 100 μL 20 mM Standard and 900 μL distilled water. For cell culture samples containing phenol red, prepare 1000 μL 1.0 mM lactate Premix by mixing 50 μL 20 mM Standard and 950 μL culture medium without serum. Dilute standard as follows. Transfer 20 μL standards into wells of a clear bottom 96-well plate.

No	Premix + H ₂ O or Medium	Vol (μL)	D-lactate (mM)
1	100 μL + 0 μL	100	2.0 or 1.0
2	80 μL + 20 μL	100	1.6 or 0.8
3	60 μL + 40 μL	100	1.2 or 0.6
4	40 μL + 60 μL	100	0.8 or 0.4
5	30 μL + 70 μL	100	0.6 or 0.3
6	20 μL + 80 μL	100	0.4 or 0.2
7	10 μL + 90 μL	100	0.2 or 0.1
8	0 μL + 100 μL	100	0

Samples. Add 20 μL sample per well in separate wells. For samples with potential endogenous enzyme activity (i.e. serum, plasma, tissue extracts), two reactions should be run: one with added Enzyme A and a No Enzyme A control. Serum and Plasma should be diluted at least 2 $\times$ with dH₂O prior to assay.

2. Reagent Preparation. Spin the Enzyme tubes briefly before pipetting. For each Sample and Standard well, prepare Working Reagent by mixing 60 μL Assay Buffer, 1 μL Enzyme A, 1 μL Enzyme B, 10 μL NAD and 14 μL MTT. Fresh reconstitution is recommended. For the No Enzyme A control, the Working Reagent includes 60 μL Assay Buffer, 1 μL Enzyme B, 10 μL NAD and 14 μL MTT.
3. Reaction. Add 80 μL Working Reagent per reaction well quickly. Tap plate to mix briefly and thoroughly.
4. Read optical density (OD₀) for time “zero” at 565 nm (520-600nm) and OD₂₀ after a 20-min incubation at room temperature.
5. Calculation. Subtract OD₀ from OD₂₀ for the standard and sample wells. Use the D_{OD} values to determine the sample D-lactate concentration from the standard curve. For samples requiring a No Enzyme A control, subtract the $\Delta\text{OD}_{\text{NoEnz}}$ value from the $\Delta\text{OD}_{\text{Sample}}$ and use this $\Delta\Delta\text{OD}$ value to determine the sample D-lactate concentration from the standard curve.

Note: if the sample OD value is higher than OD for 2 mM D-lactate standard, dilute sample in water and repeat the assay. Multiply the results by the dilution factor.

GENERAL CONSIDERATIONS

The following substances interfere and should be avoided in sample preparation: EDTA (>0.5 mM), ascorbic acid, SDS (>0.2%), sodium azide, NP-40 (>1%) and Tween-20 (>1%).

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121


www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com

D-Lactate Assay Kit (Colorimetric)

LS-K236-100 (100 Tests) • Store at -20°C



Sample Data



Standard Curve in 96-well plate assay in water.

Version: V.08.09.2018

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121
www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com