# EN FINGER TIP PULSE OXIMETER

**General description**  
Oxygen binds to hemoglobin in red blood cells when moving through the lungs. It is transported throughout the body as arterial blood. A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO2. A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO2 level.

**Measurement principle**  
Principle of the oximeter is as follows: The pulse oximeter works by applying a sensor to a pulsating arterial vascular bed. The sensor contains a dual light source and photo detector. The one wavelength of light source is 660nm, which is red light; the other is 950nm, which is infrared-red light. Skin, bone, tissue and venous vessels normally absorb a constant amount of light over time. The photo detector in finger sensor collects and converts the light into electronic signal which is proportional to the light intensity. The arterial bed normally pulsates and absorbs variable amounts of light during systole and diastole, as blood volume increases and decreases. The ratio of light absorbed at systole and diastole is translated into an oxygen saturation measurement. This measurement is referred to as SpO2.

**Diagram of operation principle**  
1. Red and Infrared-ray Emission Tube  
2. Red and Infrared-ray Recept Tube


- Precautions for use**
- Before use, carefully read the manual.
  - Operation of the fingertip pulse oximeter may be affected by the use of an electrical unit (ESU).
  - The fingertip pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO2 measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO2 measurement.
  - Do not use the fingertip pulse oximeter in an MRI or CT environment.
  - Do not use the fingertip pulse oximeter in situations where alarms are required. The device has no alarms. It is not for continuous monitoring.
  - Do not use the fingertip pulse oximeter in an explosive atmosphere.
  - The fingertip pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
  - In order to ensure correct sensor alignment and skin integrity, the maximum application time at a single site for our device should be less than half an hour.
  - Do not sterilize the device using autoclaving, ethylene oxide sterilizing, or immersing the device in liquid. The device is not intended for sterilization.

- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- This equipment complies with IEC 60601-1-2:2014 for electromagnetic compatibility for medical electrical equipment and/or systems. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device.

- Portable and mobile RF communications equipment can affect medical electrical equipment
- This equipment is not intended for use during patient transport outside the healthcare facility.
- This equipment should not be used adjacent to or stacked with other equipment.
- It may be unsafe to:
  - use accessories, detachable parts and materials not described in the instructions for use
  - interconnect this equipment with other equipment not described in the instructions for use
  - disassemble, repair or modify the equipment
- These materials that contact with the patient's skin contain medical silicone and ABS plastic enclosure are all pass the ISO10993-5 Tests for in vitro cytotoxicity and ISO10993-10 Tests for irritation and delayed-type hypersensitivity.
- When the signal is not stable, the reading may be inaccurate. Please do not rely on it.

- Contraindications**
- It is not for continuous monitoring.
  - Inaccurate measurements may be caused by**
    - Significant levels of dysfunctional hemoglobin (such as carboxy- hemoglobin or methemoglobin).
    - Intravascular dyes such as indocyanine green or methylene blue.
    - High ambient light. Shield the sensor area if necessary.
    - Excessive patient movement.
    - High-frequency electrocautery and defibrillators.
    - Venous pulsations.
    - Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
    - The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia.
    - The patient is in cardiac arrest or is in shock.
    - Fingernail polish or false fingernails.
    - Weak pulse quality low perfusion.
    - Low hemoglobin.

- Product features**
- Simple to operate and convenient to carry.
  - Small volume, light weight and low consumption index.
  - Dual color LED displays SpO2, PR, PI (Perfusion Index), Pulse bar, and waveform.
  - 7 display modes.
  - Level 1-10 adjustable brightness.
  - 2pcs AAA-size alkaline batteries; real-time battery status indication.
  - Weak or unstable signal prompt provides more accurate measurements.
  - The device automatically shuts off after no operation in 8 seconds when "Tinger out" displays.
  - Multiple-patient reusability.
- Intended use**  
The Fingertip Pulse Oximeter is a handheld non-invasive device intended for spot-checking of oxygen saturation of arterial hemoglobin (SpO2) and Pulse Rate of adult, adolescent and child patients in hospitals, hospital-type facilities and homecare.
- Operations**
- Install two AAA batteries according to the battery installation instructions.
  - Place one of your fingers into the rubber opening of the Pulse Oximeter.
  - Press the switch button one time on front panel to turn the pulse oximeter on.
  - Keep your hands still for the reading. Do not shake your finger during the test. It is recommended that you do not move your body while taking a reading.
  - Read the data from the display screen.
  - The display modes are as follows:

- Notes:**
- Short press the power button to switch the display modes.
  - Long press the power button to adjust the brightness of the oximeter. There are 10 levels of brightness. The default is level four.
  - Take out your finger, the screen displays "FingerOut". It means the measurement ending.

**Reading your results**

Battery low indicator  
Signal indication  
Oxygen saturation: 98.34% (SpO2), 197.72 (PR), 96.72 (PI)

- Notes:**
- The pulse bar less than 30% indicates signal inadequacy and the displayed SpO2 or pulse rate value is potentially incorrect.
  - If the screen display "7", it means the signal is unstable, please keep your hands still and retry.
  - PI means Perfusion Index.

**Battery installation**

- Slide the battery door cover horizontally along the arrow shown on the picture.
- Install two AAA batteries into the battery compartment. Match the plus (+) and minus (-) signs in the compartment. If the polarities are not matched, damage may be caused to the oximeter.
- Close the battery door.

**Notes:** Please remove the batteries if the pulse oximeter will not be used for long periods of time. Please replace the battery when the power indicator starting flickering.

**Using the lanyard**

- Thread thinner end of the lanyard through the hanging hole.
- Thread thicker end of the lanyard through the threaded end before pulling it tightly.

**Warnings!**

- Keep the oximeter away from young children. Small items such as the battery door, battery, and lanyard are choking hazards.
- Do not hang the lanyard from the device's electrical wire.
- Please notice that the lanyard which is tied to the oximeter may cause strangulation due to excessive length.

**Maintenance and storage**

- Replace the batteries in a timely manner when low voltage lamp is lighted.
- Clean surface of the fingertip oximeter before it is used as diagnosis for patient.
- Remove the batteries if the oximeter is not operated for a long time.
- It is best to store the product in -25°C ~ +70°C and ±93% humidity.
- Keep a dry place. Extreme moisture may affect oximeter lifetime and may cause damage.
- Dispose of battery properly, follow any applicable local battery disposal laws.

**Cleaning and disinfecting the fingertip pulse oximeter**

- Please use medical alcohol to clean the silicone tubing the finger inside of oximeter with a soft cloth dampened with 70% isopropyl alcohol. Also clean the biting tested finger using alcohol before and after each test.
- Do not use any spray liquids or other cleaning products, and do not allow any liquid to enter any openings in the device. Allow the oximeter to dry thoroughly before reuse.
- The fingertip pulse oximeter requires no routine calibration or maintenance other than replacement of batteries.
- The applied parts touching the patients' body are required to be disinfected once after each use. The recommended disinfectants include ethanol 70%, isopropanol 70%, glutaraldehyde-type 2% liquid disinfectants.
- Disinfection may cause damage to the equipment and is therefore not recommended for this pulse oximeter unless otherwise indicated in your hospital's servicing schedule. Clean the pulse oximeter before disinfecting it.

**CAUTION:** Never use EO or formaldehyde for disinfection.  
The use life of the device is five years when it is used for 15 measurements every day and 10 minutes per one measurement. Stop using and contact your local service center if one of the following cases occurs:

- An error in the Possible Problems and solutions is displayed on screen.
- The oximeter cannot be powered on in any case and not the reasons of battery.
- There is a crack on the oximeter or damage on the display resulting readings cannot be identified; the spring is invalid; or the key is unresponsive or unavaliable.

**Specifications**

- Display Type: OLED display
- SpO2
  - Display range: 0%-100%
  - Measurement range: 70%-100%
  - Accuracy: 70% - 100% ±2%, 90% - 69% no definition.
  - Resolution: 1%

**Pulse rate**

- Display range: 30bpm-250bpm
- Measure range: 30bpm-250bpm
- Accuracy: 30bpm-79bpm, ±2bpm; 100bpm-250bpm, ±2%
- Resolution: 1bpm

**Perfusion index**

- Display range: 0.1%-20%
- Measure range: 0.2-20.0%
- Accuracy: 0.2%-1.5%, ±0.26bpm; 1.1%-20.0%, ±20%
- Resolution: 0.1%

**Probe LED specifications**

RED	Wavelength	Radiant Power
IR	660±2nm	3.2mw
	905±10nm	2.4mw

**Power requirements**


- Two AAA alkaline Batteries
- Power consumption: Less than 0.6mA
- Battery Life: Two AAA 1.5, 1200mAh alkaline batteries could be continuously operated as long as 24 hours.

**Environment requirements**

- Operation Temperature: 5°C ~ 40°C
- Storage Temperature: -25°C ~ 70°C
- Relative Humidity: 15% ~ 93% (non-condensing in operation; ±93% no condensation in storage/transport)
- Atmosphere pressure: 70kPa ~ 106kPa

**Equipment response time**


- As shown in the figure
- Response time of slower average is 8s.

**Classification**

- According to the type of protection against electric shock: INTERNALLY POWERED EQUIPMENT.
- According to the degree of protection against electric shock: TYPE BF APPLIED PART. (The application part is rubber instead of the Pulse Diometer).
- According to the degree of protection against ingress of dust and water: IP22
- According to the mode of operation: CONTINUOUS OPERATION

**Clinical study summary**  
The following tables are provided to disclose actual operation observed in the clinical validation study of healthy adult volunteers. The ARMS value analysis system and Band-Altman plot are shown as follows:

**Possible problems and solutions**

Problems	Possible reason	Solution
SpO2 or PR can not be shown normally	1. Finger is not inserted correctly 2. Patient's SpO2 value is too low to be measured	1. Retry by inserting the finger 2. There is excessive illumination 3. Try some more times. If you can make sure no problem exists in the product, please go to a hospital timely for exact diagnosis.
SpO2 or PR is shown unstable	1. Finger might not be inserted deep enough. 2. Excessive patient movement	1. Retry by inserting the finger 2. Be calmness
The oximeter cannot be powered on	1. No battery or low power of battery 2. Batteries might be installed incorrectly 3. The oximeter might be damaged	1. Please replace batteries 2. Please reinstall the batteries 3. Please contact with local customer service centre
Indication lamps are suddenly off	1. The product is automatically powered off when no signal is detected longer than 8 seconds 2. The battery power is too low to work	1. Normal 2. Replace the batteries
"Error" or "Error" is displayed on screen	1. Err 3 means red emission LED is damaged 2. Err 4 means infra-red emission LED is damaged	1. Check the red emission LED 2. Check the infra-red emission LED
"Error"	Err 6 means screen failure	Please contact the local customer service centre
"Error" is displayed on screen	Err 7 means all the emission LED or reception diode is damaged.	Please contact with local customer service centre

**Symbol definitions**

Symbol	Definition	Symbol	Definition
	Type BF applied part		Attention
	Follow instructions for use	%SpO2 symbol"/>	Oxygen saturation
	Pulse rate (BPM)		Lower power indication
	No SpO2 Alarm		Serial No.
	Storage temperature and relative humidity		The degree of protection against ingress of dust and water
	Authorised representative in the European community		European Union approval
	Conformity to WEEE Directive		Indicates the signal is not stable

**Box contents**

- Fingertip pulse oximeter
- Lanyard
- Two AAA batteries
- One instruction manual

**Notes:** The illustrations used in this manual may differ slightly from the appearance of the actual product. The specifications are subject to change without prior notice.

**2 YEAR WARRANTY**  
FKA Brands Ltd guarantees this product from defect in material and workmanship for a period of 2 years from the date of purchase, except as noted below. This FKA Brands Ltd product guarantee does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorised accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of FKA Brands Ltd. This guarantee is effective only if the product is purchased and operated in the UK/ EU. A product that requires modification or adaptation to enable it to operate in any country other than the country for which it was designed, manufactured, approved (A) or authorised (or repair of products damaged by these modifications is not covered under this guarantee. FKA Brands Ltd shall not be responsible for any type of incidental, consequential or special damages.  
To obtain guarantee service on your product, return the product post-paid to your local service centre along with your dated sales receipt (proof of purchase). Upon receipt, FKA Brands Ltd will repair or replace, as appropriate, your product and return it to you, post-paid. Guarantee is solely through the HomeDics Service Centre. Service of this product by anyone other than the HomeDics Service Centre voids the guarantee.  
This guarantee does not affect your statutory rights. For your local HomeDics Service Centre, go to [www.homedics.co.uk/servicecentres](http://www.homedics.co.uk/servicecentres)

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