EN

FINGERTIP PULSE OXIMETER

General description

Oxygen binds to hemoglobin in red blood cells when moving through the lungs. It is transported throughout the body as arterial blood. A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO2. A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO2 level.

Measurement principle

Principle of the device is as follows: The pulse oximeter works by applying a sensor to a pulsating arterial vascular bed. The sensor contains a dual light source and photo detector. The true wavelength of light source is 660nm, which is red light; the other is 950nm, which is infrared-red light. Skin, bone, tissue and venous vessels normally absorb a constant amount of light over time. The photo detector in finger sensor collects and converts the light into electronic signal which is proportional to the light intensity. The arterial bed normally pulsates and absorbs variable amounts of light during systole and diastole, as blood volume increases and decreases. The ratio of light absorbed at systole and diastole is translated into an oxygen saturation measurement. This measurement is referred to as SpO2.

Diagram of operation principle

1. Red and Infrared-rd Emission Tube

2. Red and Infrared-rd Receipt Tube

Precautions for use

1. Before use, carefully read the manual.

2. Operation of the fingertip pulse oximeter may be affected by the use of an electrocussorial unit (ESU).

3. The fingertip pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO2 measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO2 measurement.

4. Do not use the fingertip pulse oximeter in an MRI or CT environment.

5. Do not use the fingertip pulse oximeter in situations where alarms are required. The device has no alarms. It is not for continuous monitoring.

6. Do not use the fingertip pulse oximeter in an explosive atmosphere.

7. The fingertip pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.

8. In order to ensure correct sensor alignment and skin integrity, the maximum application time in a single site for our device should be less than half an hour.

9. Do not sterilize the device using autoclaving, ethylene oxide sterilization, or immersing the device in liquid. The device is not intended for sterilization.

10. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.

11. This equipment complies with IEC 60601-1-2:2014 for electromagnetic compatibility for medical electrical equipment and/or systems. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device.

12. Portable and mobile RF communications equipment can affect medical electrical equipment.

13. This equipment is not intended for use during patient medical electrical healthcare facility.

14. This equipment should not be used adjacent to or stacked with other devices.

15. It may be unsafe to use accessories, detachable parts and materials not described in the instructions for use.

16. Use accessories, detachable parts and materials not described in the instructions for use.

17. When the device is not in use, carefully read the manual.

Specifications

• Display Type: LED display

• SpO2

• Display range: 0%-100%

• Measurement range: 70%-100%

• Accuracy: 70%-100%±1.5%, 0%-65% no definition.

• Resolution: 1%

Pulse rate

• Display range: 30bpm-250bpm

• Measure range: 30bpm-250bpm

• Accuracy: 30bpm-199bpm: ±2bpm; 100bpm-250bpm: ±2%

• Resolution: 1bpm

Perfusion index

• Display range: 0.1%-20%

• Measure range: 0.1%-20.0%

• Accuracy: 0.2%-1.5%, ±0.2bpm; 1.1%-20.0%, ±2.0%

• Resolution: 0.1%

Probe LED specifications

	Wavelength	Radiant Power
RED	660±3nm	2.2mw
IR	950±10nm	2.4mw

Power requirements

• Two AAA alkaline Batteries

• Power consumption: Less than 40mA

• Battery Life: Two AAA 1.5V 1200mAh alkaline batteries can be continuously operated as long as 24 hours.

Environment requirements

• Operation Temperature: -5°C ~40°C

• Storage Temperature: -25°C ~70°C

• Ambient Humidity: 15% ~93% no condensation in operation; >93% no condensation in storage/transport

• Atmosphere pressure: 70kPa~106kPa

Equipment response time

As shown in the figure

Response time of slower average is 8s.

Classification

• According to the type of protection against electric shock: INTERNALLY POWERED EQUIPMENT;

• According to the degree of protection against electric shock: TYPE BF APPLIED PART. (The application part is rubber insided of the Pulse Oximeter);

• According to the degree of protection against ingress of dust and water: IP22;

• According to the mode of operation: CONTINUOUS OPERATION

Clinical study summary

The following details are provided to disclose actual performance observed in the clinical validation study of healthy adult volunteers. The ARMS value analysis statement and Bland-Altman plot of data is shown as following:

ARMS value analysis statement

Item	90-100	-90--100	-40--80
Fpts	78	66	43

Bland-Altman

ARMS	1.62	0.40	-0.48
Bias	1.66	1.46	1.93

Product features

1. Simple to operate and convenient to carry.

2. Small volume, light weight and low power consumption.

3. Dual color OLED displays SpO2, PR (Perfusion Index), Pulse bar, and waveform.

4. 7 display mode.

5. Level 1-10 adjustable brightness.

6. 2pcs AAA-size alkaline batteries; real-time battery status indication.

7. Weak or unstable signal prompt provides more accurate measurement.

8. The device automatically shuts off after no operation in 3 seconds when "Tinger off" displays.

9. Multiple patient reusability.

Intended use

The Fingertip Pulse Oximeter is a handheld non-invasive device intended for spot-checking of oxygen saturation of arterial hemoglobin (SpO2) and Pulse Rate of adult, adolescent and child patients in hospitals, hospital-type facilities and homecare.

Operations

1. Install the AAA batteries according to the Battery Installation instructions.

2. Place one of your fingers into the rubber opening of the Pulse Oximeter.

3. Press the switch button one time to turn on /stand by the pulse oximeter.

4. Press your hands still for the reading. Do not shake your finger during the test.

5. It is recommended that you do not move your body while taking a reading.

6. Read the data from the display screen.

7. The display modes are as follows:

%SpO2

PR bpm

%SpO2

PR bpm

%SpO2

PR bpm

%SpO2

PR bpm

96 34

96 72

96 72

96 72

Notes:

1. The pulse bar is more than 30% indicates signal inadequacy and the displayed SpO2 or pulse rate value is potentially incorrect.

2. If the screen display "T", it means the signal is unstable, please keep your hands still and rely.

3. PR means Perfusion Index.

Battery installation

1. Slide the battery door cover horizontally along the arrow shown as the picture.

2. Install two AAA batteries into the battery compartment. Match the plus (+) and minus (-) signs in the compartment. If the polarities are not matched, damage may be caused to the oximeter.

3. Close the battery door.

Notes:

1. Please remove the batteries of the pulse oximeter after use for long periods of time.

2. Please replace the battery when the power indicator starting flickering.

Using the lanyard

1. Thread thinner end of the lanyard through the hanging hole.

2. Thread thicker end of the lanyard through the threaded end before pulling it tightly.

Box contents

1. Fingertip pulse oximeter

2. One lanyard

3. Two AAA batteries

4. One instruction manual

Notes:

The illustrations used in this manual may differ slightly from the appearance of the actual product. The specifications are subject to change without prior notice.

Warnings

1. Keep the oximeter away from young children. Small items such as the battery door, battery, and lanyard are choking hazards.

2. Do not hang the lanyard from the device's electrical wire.

3. Please notice that the lanyard which is tied to the oximeter can cause strangulation due to excessive length.

Maintenance and storage

1. Replace the batteries in a timely manner when low voltage lamp is lighted.

2. Clean surface of the fingertip oximeter before it is used in diagnosis for patients.

3. Remove the batteries if the oximeter is not operated for a long time.

4. It is best to store the product in -25°C ~+70°C and 93% humidity.

5. Keep in a dry place. Extreme moisture may cause internal corrosion and may cause damage.

6. Dispose of battery properly, follow any applicable local battery disposal laws.

Cleaning and disinfecting the fingertip pulse oximeter

• Please use medical alcohol to clean the oximeter touching the finger inside of oximeter with a soft cloth dampened with 70% isopropyl alcohol. Also clean the long tested finger using alcohol before and after each test.

• Do not put or spray liquids into the device, and do not allow any liquid to enter any openings in the device. Allow the oximeter to dry thoroughly before reuse.

• The fingertip pulse oximeter requires no routine calibration or maintenance other than replacement of batteries.

• The applied parts touching the patient's body are required to be disinfected once after each use. The recommended disinfectants include: ethanol 70%, isopropanol 70%, glutaraldehyde-type 2 liquid disinfectants.

• Disinfection may cause damage to the equipment and it is therefore not recommended for this pulse oximeter otherwise indicated in your hospital's sterilizing schedule. Clean the pulse oximeter before disinfecting it.

• CAUTION: Never use EO or formaldehyde for disinfection.

The use life of the device is five years when it is used for 15 measurements every day for 10 minutes per use measurement. Stop using and contact your local service center if one of the following cases occurs:

• An error in the Possible Problems and Solutions is displayed on screen.

• The oximeter cannot be powered on in any case and not the reason of the battery.

• There is a crack on the oximeter or damage on the display resulting readings cannot be identified; the spring is invalid; or the key is unresponsive or unavailable.

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