

WARRANTY

Your Dynarex Product is warrantied to be free of defects in materials and workmanship for (2) two years for the original purchaser. This item was built to exacting standards and carefully inspected prior to shipment. This Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

The Warranty shall not apply under the following conditions:

- Problems arising from normal wear
- Problems arising from failure to adhere to the product instructions
- Problems arising from misuse, negligence, accident or improper operation, maintenance or storage
- Problems arising from modifications or unauthorized repairs, parts or attachments
- Products where the serial number has been removed or defaced
- Problems with non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement

WARRANTY CONTINUED

Dynarex shall not be liable for any consequential or incidental damages whatsoever. Dynarex shall repair or replace defective products at its option. The foregoing warranty is exclusive and in lieu of other express warranties, if any, including the implied warranties of merchantability and fitness of a particular purpose. The remedy for any violation of the implied warranty shall be limited to repair or replacement of the defective product pursuant to the terms contained herein.

If you have a question about your Dynarex device or this warranty, please contact an authorized Dynarex dealer.

 **Manufactured for:**
Dynarex Corporation
10 Glenshaw Street
Orangeburg, NY 10962
USA • www.dynarex.com
Made in China

R210104



SYMBOL GLOSSARY
For an explanation of symbols used in Dynarex packaging, visit dynarex.com/symbols.php



Oxygen Regulators CGA 870 Style

User's Manual



Reorder Nos. **36350-36353**

Oxygen Regulators CGA 870 Style

Barb Outlet • Brass-Sleeved

Reorder Nos. **36354-36357**

Mini Oxygen Regulators CGA 870 Style

Barb Outlet • Brass-Sleeved

Reorder Nos. **36358-36359**

Extended Oxygen Regulators CGA 870 Style

DISS Outlet • Brass-Sleeved

Reorder Nos. **36360**

Oxygen Regulators CGA 870 Style

Barb Outlet • All Brass

Reorder Nos. **36361**

Oxygen Regulators CGA 870 Style

Barb Outlet • 2 DISS Outlets • All Brass


Please read the user's manual closely before using!

OXYGEN FLOW CONTROL REGULATORS

Oxygen Flow Control Regulators are pressure-reducing devices which lower oxygen pressure from a cylinder to levels which can be safely used by the user. Oxygen Flow Control Regulators are intended for the administration of oxygen to patients who need increased oxygen levels to improve their breathing conditions.

WARNINGS

- This unit should only be used on patients who are breathing; do not use to resuscitate non-breathing patients.
- Do not use with oil, grease or any other flammable material. Keep out of reach of children.
- Do not smoke in an area where oxygen is in use.
- If the pressure relief valve on the underside of the regulator begins to vent oxygen, turn off the system and remove the unit until it can be serviced.
- NEVER attempt to repair a regulator or cylinder yourself.
- (36361 only) This device should only be used for emergency purposes.
- NEVER try to remove the regulator unless the cylinder is off and pressure has been relieved.



SETUP

- Inspect the regulator for oil or grease. If these are visible, DO NOT use the regulator.
- Be sure the sealing washer is in place on the yoke inlet. Inspect the washer. It should be free of cracks and debris. If cracks and debris exist, the washer should be replaced immediately.
- Slip the yoke over the cylinder post and fit the pins into the holes on the cylinder valve. The pins should fit easily into these holes.
- Turn the T-Handle clockwise until the screw point is seated in the dimple on the cylinder valve. Be sure the regulator is tightened snugly onto the cylinder.
- Attach the oxygen supply tubing to the barbed outlet connection.

USE

- Be sure the regulator is in the “off” position by rotating the knob counterclockwise to the “0” position.
- For safety, be sure you are not directly in front of or behind the regulator when opening the cylinder valve.
- SLOWLY turn the cylinder valve on (counterclockwise), about one full turn.
- If you can hear a hissing sound, there is a leak in the systematizer. Turn the cylinder off and turn the regulator on to relieve any built-up pressure. Try tightening the regulator to the cylinder and opening the cylinder again. If this does not help, shut the system down again and contact a service representative.
- If you suspect leakage, inspect and/or replace the yoke washer.
- Adjust the flow rate setting by turning the knob clockwise until the desired setting shows through the window. For emergency situations, this should be at least 6 LPM.
- If the cylinder pressure falls below 300 PSI (500 PSI for 5227 only), you should exchange the cylinder for a full one.
- To turn off the system, close the cylinder valve, wait for oxygen to flow out of the regulator and turn the regulator to the “off” position.