

Quick Read Forehead

User's Guide TH068, TH073

Please be sure to read all warnings, instructions and notations before using this product. Please keep these instructions as they contain important information.

Basics of Your Thermometer

Care and Maintenance

Cleaning Temperature Sensor Lens:

- Clean very gently with a soft, alcohol-moistened cloth.
- Do not use an excessive amount of alcohol.
- Do not pour alcohol directly on sensor.

Cleaning Thermometer Body:

- Use a soft, dry cloth to clean thermometer body.

Notes:


- **DO** remove battery if stored for a long period of time.
- **DO** store thermometer in a dry place at room temperature.
- **DO** take special care not to scratch the sensor or display window.
- **DO NOT** use abrasive cleaning agents, thinners or benzene for cleaning.
- **DO NOT** expose thermometer to extreme temperatures, humidity or direct sunlight.
- **DO NOT** drop thermometer.
- **DO NOT** submerge any part of the product in water.

Normal Body Temperature

Temperature readings are different from person to person and also vary depending upon the body location from which the temperature is measured. Forehead temperatures tend to be comparable to oral temperatures and approximately 1° lower than rectal temperatures for the same individual. Temperature readings are affected by many factors including exercise, drinking hot or cold beverages prior to taking a reading as well as technique.

"Normal" temperature refers to each individual's average temperature while healthy. You may want to measure and record your child's temperature a few times when he or she is well, so that you will know what is normal for your own child. On average, temperature ranges are as follows:

TEMPERATURE METHOD	NORMAL TEMPERATURE	FEVER RANGE
Forehead	97°F to 100.3°F	100.3°F +

WARNING

- This is **NOT** an ear thermometer. Take temperature measurements from the center of forehead.
- Use caution when taking temperature readings of bath water and food. This product is not intended to touch the person or item whose temperature is being measured.
- **ALWAYS** supervise your child when using the thermometer. **NEVER** allow your child to use it alone.
- **NEVER** allow children to play with this thermometer. It is not a toy. This is a sensitive electronic device and may be damaged if not used properly under parental supervision.
- **DO NOT** expose the thermometer to extreme temperatures, humidity, direct sunlight or shock, and/or dust.
- Keep battery out of child's reach. Swallowing a battery could be fatal. Battery should not be charged or placed into extreme heat as they may explode.
- Used batteries should be disposed of properly.

Other Battery Information:

We recommend alkaline batteries for longer battery life. **DO NOT** mix old and new batteries. **DO NOT** mix alkaline, standard (carbon zinc) or rechargeable (nickel-cadmium) batteries. Non-rechargeable batteries are not to be recharged. Rechargeable batteries are to be removed from the unit before being charged. Rechargeable batteries are only to be charged under adult supervision. Only batteries of the same or equivalent type as recommended are to be used. Batteries are to be inserted with the correct polarity. The supply terminals are not to be short-circuited. Check that all contact surfaces are clean and bright before installing batteries.


How to Measure Temperature Make sure the thermometer has been at room temperature for 30 minutes before use.

1. Turn Thermometer On

Press and release the ON/MEASURE button. Thermometer runs a brief self-check, then beeps and displays the last reading with the memory icon.

2. Position Thermometer Near Forehead


Aim the temperature sensor lens at center of the forehead. Position lens approx. 1/2 inch from forehead. Do not touch forehead. **Hint:** If a child is uncomfortable seeing the thermometer coming toward his face, try turning the thermometer upside-down. Just be sure the lens is aimed directly at the forehead.

3. Measure Temperature

Hold thermometer steady. Do not move thermometer while measuring temperature. Press and release the ON/MEASURE button. The "Measuring in Progress" arrow displays and thermometer beeps once when measurement is complete.

4. Temp Readings with Fever/No Fever Icons

The "No Fever" icon displays when temp is in normal range. Thermometer gives additional beeps and the "Fever" icon displays when fever temp is detected. **NOTE:** High, prolonged fever requires medical attention, especially for young children. Please contact your physician. Fever chart included.

5. Auto Shut-Off


Thermometer beeps and shuts off automatically after 20-30 seconds. Wait 1 minute before taking another measurement, so that thermometer has time to re-acclimate to the room's temperature.

Troubleshooting When a malfunction or incorrect temperature measurement occurs, an error message will appear.


Display	Meaning	Recommended Action	Display	Recommended Action
	• Measured temperature is below this thermometer's range. See Technical Specifications section for Measurement Temperature Range.	• Ensure thermometer is being positioned correctly (See How to Measure Temperature). • Ensure icon does not appear. If it does, replace batteries. • Consistent readings may indicate that the sensor is dirty (see Care and Maintenance).		• Batteries are low. Replace them. (See Battery Information section).
	• Measured temperature is above this thermometer's range. See Technical Specifications section for Measurement Temperature Range.	• Ensure thermometer is being positioned correctly (See How to Measure Temperature). • Allow thermometer to reach room temperature, then re-take temperature measurement.		• For best accuracy, wait one minute between temperature measurements so that thermometer has time to re-acclimate to the room's temperature.
	• Thermometer may have been stored in a cold or hot environment. • This thermometer cannot function if the environment's temperature is outside the permissible range. See Operation Environment for range.	• Allow thermometer to reach room temperature, then try again.		• Do not return this product to the place of purchase. If any parts are missing, email consumer@digusa.com, call Consumer Relations at (800) 544-1108, or fax at (800) 207-8182. You can also visit our website at www.safety1st.com. Have ready the model number (TH068, TH073) and date code (manufacture date located on package). • Install new batteries and re-try function before calling Consumer Relations on electronics.

How to Switch Between Fahrenheit & Celsius

- Use a Phillips head screwdriver to loosen battery cover screw. Remove cover.

- While thermometer is ON, gently press and release the °C/°F button.


- Place the battery cover on the thermometer and tighten the screw to secure it in place.

How to Replace the Batteries


Low Battery Warning:

- When the icon appears on the display screen when powering up the thermometer, the batteries are low but may still work for a very short time.
- When the icon appears with the **Lo** message the batteries are exhausted and need to be replaced.

To Replace Batteries:

- 1 Use a Phillips head screwdriver to loosen battery cover screw. Remove cover.
- 2 Remove old batteries
- 3 Insert two AAA (or equivalent) batteries, correctly positioned with positive "+" and negative "-" as shown in illustration.
- 4 Place the battery cover on the thermometer and tighten the screw to secure it in place.

Technical Specifications

ASTM laboratory accuracy requirements in the display range of 37 to 39 °C (98 to 102 °F) for IR thermometers is +/- 0.2 °C (+/- 0.4 °F), whereas for electronic thermometers, the requirement per ASTM Standard E 1112-86 is +/- 0.1 °C (+/- 0.2 °F).

Measurement Temperature Range: 71.6°F to 104.0°F

Measurement Accuracy: +/- 0.4°F Temps 96.8°F to 102.2°F
+/- 0.5°F All Other Temps

Operation Environment: 60.8°F to 104.0°F
Humidity up to 85% (non condensing)

Storage Environment: -4°F to 122°F
Humidity up to 85% (non condensing)

Test Protocol:

- 1 Place the IR thermometer into water bath chamber for at least 2 hours.
- 2 After 2 hours, press MEASURE button without releasing key until the decimal point disappears.
- 3 Put the thermometer into the fixture of black body.
- 4 Press and release MEASURE button, the device will beep once and measurement will be shown on the display.

One (1) Year Limited Warranty

Dorel Juvenile Group, Inc. warrants to the original purchaser that this product (Thermometer), is free from material and workmanship defects when used under normal conditions for a period of one (1) year from the date of purchase. Should the product contain defects in material or workmanship Dorel Juvenile Group, Inc., will repair or replace the product, at our option, free of charge. Purchaser will be responsible for all costs associated with packaging and shipping the product to Dorel Juvenile Group Consumer Relations Department at the address noted on the front of this document and all other freight or insurance costs associated with the return. Dorel Juvenile Group will bear the cost of shipping the repaired or replaced product to the purchaser. Product should be returned in its original package accompanied by a proof of purchase, either a sales receipt or other proof that the product is within the warranty period. This warranty is void if the owner repairs or modifies the product or the product has been damaged as a result of misuse.

This warranty excludes any liability other than that expressly stated previously, including but not limited to, any incidental or consequential damages.

SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS THAT MAY VARY FROM STATE TO STATE.

©2014 Dorel Juvenile.
www.digusa.com (800) 544-1108
www.safety1st.com
Made in CHINA.
Styles and colors may vary.
Distributed by Dorel Juvenile Group, Inc.,
2525 State St., Columbus, IN 47201-7494
Dorel Distribution Canada, 873 Hodge, St. Laurent, QC H4N 2B1
03/12/14 4358-6388A

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Distributed by (distribuido por) Dorel Juvenile Group, Inc.,