Homedics

body analyser scale

✓ Register your product today at
www.homedics.co.uk/product-registration

3 YEAR GUARANTEE

SCL-BC145-BK3R

PRODUCT FEATURES:

PREPARING YOUR SCALE:

1. Open the battery compartment on the scale underside.
2. Insert the supplied batteries observing the polarity signs (+ and -) inside the battery compartment.
3. Close the battery compartment.
4. Position scale on a firm flat surface.
5. Set up the scale by briefly stepping on and then off the platform, allowing the display to set to zero.
6. With the display showing zero, select KG, ST/LB or LB weight units by pressing either arrow key until your preferred units are displayed.

Your scale is now ready to use.

If you move the scale between weighings, repeat the setup process in step 5 before next use.

SETTING UP YOUR USER PROFILE:

For the scale to calculate your body composition, you need to enter some basic details which the scale will store in your user profile. The scale can store details for up to 9 users:

1. Activate the scale by pressing firmly on the centre of the platform. The display will show zero.
2. Press the 'SET' button and the user number will be displayed (flashing).
3. Use the arrow keys to select your preferred user number, then press 'SET'
4. Use the arrow keys to select male or female, then press 'SET'
5. Use the arrow keys to enter your age, then press 'SET'
6. Use the arrow keys to enter your height, then press 'SET'
7. Now step on the scale to weigh yourself for the first time and to complete the user setup.

USING YOUR SCALE:

1. Step on the scale and stand still.
2. The scale will show your weight reading. It will lock on to your weight reading and then compute for a few seconds before sequencing twice through your readings of weight, body fat %, body water %, bone mass, muscle %, BMI and BMR readings before switching off.

NOTE: If the scale is unsure which user you are, it will display two or more options for you to select from. Simply use the arrow key to select the correct user number.

For more information about your readings please scan the QR code:

WARNING INDICATORS:

- 'Lo' is displayed when batteries are low
- 'Err' is displayed when the weight limit is exceeded
- 'Err2' is displayed if the scale is unable to take a body fat reading. This may happen if you have not set up a user profile or if you are not in bare feet

ADVICE FOR USE AND CARE:

- Always weigh yourself on the same scale placed on the same floor surface. Do not compare weight readings from one scale to another as some differences will exist due to manufacturing tolerances.
- Placing your scale on a hard, even floor will ensure the greatest accuracy and repeatability.
- Weigh yourself at the same time each day, before meals and without footwear. First thing in the morning is a good time.
- Your scale rounds up or down to the nearest increment. If you weigh yourself twice and get two different readings, your weight lies between the two.
- Clean your scale with a damp cloth. Do not use chemical cleaning agents.
- Do not allow your scale to become saturated with water as this can damage the electronics.

TROUBLESHOOTING GUIDE:

If you experience any difficulty in the use of your scale:

- Check the battery is correctly fitted.
- Check that you have selected your preference of stones/pounds, kilograms or pounds.
- Check that the scale is on a flat, level floor and not touching against a wall.
- If when you use the scale there is no display at all, replace the batteries.
- Avoid use on carpeted or uneven surfaces.
- If you have particularly dry skin, you may need to moisturize your feet to get a reading through the metal contacts. Alternatively, use after a bath or shower, but always dry feet first to avoid a slip hazard.

READ ALL INSTRUCTIONS BEFORE USE. SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE.

- DO NOT use if you have either implanted or external medical devices such as a cardiac pacemaker, artificial lung, or wearable electrocardiograph.
- DO NOT use while pregnant.
- Keep batteries away from children.
- Keep scale dry and away from water. Do not step onto the glass with wet feet as it will be a slip hazard.
- DO NOT jump on the scale, or stand on the outer edge of the platform which could cause it to tip.
- DO NOT use the scale if the glass appears cracked or broken.
- DO NOT overload the scale beyond its stated capacity.
- This product is intended for personal monitoring only and should not be used for treating or diagnosing any medical conditions.

3 YEAR GUARANTEE

FKA Brands Ltd guarantees this product from defect in material and workmanship for a period of 3 years from the date of purchase, except as noted below. This FKA Brands Ltd product guarantee does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorised accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of FKA Brands Ltd. This guarantee is effective only if the product is purchased and operated in the UK / EU. A product that requires modification or adaptation to enable it to operate in any country other than the country for which it was designed, manufactured, approved and / or authorised, or repair of products damaged by these modifications is not covered under this guarantee. FKA Brands Ltd shall not be responsible for any type of incidental, consequential or special damages.

To obtain guarantee service on your product, return the product post-paid to your local service centre along with your dated sales receipt (as proof of purchase). Upon receipt, FKA Brands Ltd will repair or replace, as appropriate, your product and return it to you, post-paid. Guarantee is solely through HoMedics Service Centre. Service of this product by anyone other than the HoMedics Service Centre voids the guarantee. This guarantee does not affect your statutory rights.

For your local HoMedics Service Centre, go to www.homedics.co.uk/servicecentres

WEEE explanation

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

Battery directive

This symbol indicates that batteries must not be disposed of in the domestic waste as they contain substances which can be damaging to the environment and health. Please dispose of batteries in designated collection points.

CAUTION

- Replacement of the battery with an incorrect type can defeat a safeguard.
- Do not dispose of a battery into fire or a hot oven, or mechanically crushing or cutting of a battery, that can result in an explosion;
- Do not leave a battery in an extremely high temperature surrounding environment that can result in an explosion or the leakage of flammable liquid or gas
- Do not expose the battery to extremely low air pressure that may result in an explosion or the leakage of flammable liquid or gas

ELECTROMAGNETIC COMPATIBILITY

This product is intended for domestic use.

ADVICE TO THE USER

This product is a sensitive electronic instrument and as such may temporarily be affected by Radio Transmitting Devices being used in close proximity to it (such as Mobile Phones, Walkie Talkies, CB Radios, Radio Model Controllers, and some Microwave Ovens, etc.) Additionally, extreme levels of Electrostatic Interference may cause this product to temporarily malfunction. In such case it may be necessary to remove and re-install the battery to re-establish normal working. If the problem persists then contact your local agent.

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FKA Brands Ltd, Somerhill Business Park, Tonbridge, Kent TN11 0GP, UK
FKA Brands Ltd, 29 Earlsfort Terrace, Dublin D02AY28, Ireland