# SALTER®

SINCE 1760

## Body Analyser Scale

### Instruction Manual

EN | Please retain instructions for future reference.

#### Safety Instructions

**BATTERY WARNING: KEEP OUT OF THE REACH OF CHILDREN.**

This product is not a toy.

This product contains no user serviceable parts. Only a qualified electrician should attempt repairs.

If any part of the appliance is malfunctioning or if it has been dropped or damaged, cease using the product immediately to avoid potential injury.

Do not place, store or use the product near direct sunlight or direct heat sources or in a damp/high humidity environment.

This product is a high precision device; handle it with care.

This product is intended for domestic use only.

**CAUTION: DO NOT USE THE ANALYSER FUNCTION IF YOU HAVE A PACEMAKER OR OTHER MEDICAL DEVICE FITTED. IF IN DOUBT, CONSULT A DOCTOR.**

#### Battery Safety

This product is supplied with 1 x CR2032 battery. Only batteries of the same or equivalent type as those recommended should be used with this product.

Batteries are to be inserted with the correct polarity.

Exhausted batteries should be removed from the product to avoid leakage.

#### WARNING

##### THIS PRODUCT CONTAINS A BUTTON/COIN BATTERY

Swallowing can cause severe or fatal injuries within 2 hours.

Keep batteries out of reach of children (whether new or used).

If you think batteries may have been swallowed or placed inside any part of the body, seek immediate medical attention.

**WARNING: IF THE BATTERY COMPARTMENT DOES NOT CLOSE SECURELY, STOP USING THE PRODUCT AND KEEP IT AWAY FROM CHILDREN.**

**CAUTION: VIGILANCE MUST BE MAINTAINED WITH "FLAT" OR SPARE BUTTON OR COIN BATTERIES IN THE HOME AND THE PRODUCTS THAT CONTAIN THEM.**

There are no obvious specific symptoms when a button or coin battery is stuck in a child's oesophagus (food pipe). However, the child might exhibit the following symptoms:

- Coughing, gagging or drooling
- Appear to have a stomach upset or virus
- Vomiting
- Pointing to their throat or stomach
- Pain in their abdomen, chest or throat
- Tiredness and lethargy
- Be quieter or more clingy than usual
- Loss or reduction of appetite
- Refusal or inability to eat solid foods

These sorts of symptoms may fluctuate.

**A SPECIFIC SYMPTOM OF BUTTON OR COIN BATTERY INGESTION IS VOMITING FRESH (BRIGHT RED) BLOOD. SEEK IMMEDIATE MEDICAL ATTENTION IF THIS OCCURS.**

#### Care and Maintenance

Wipe the exterior of the scale with a soft, damp cloth and allow it to dry thoroughly.

Never immerse the scale in water or any other liquid.

Do not use harsh or abrasive cleaning detergents or scourers to clean the scale, as this could cause damage.

#### Before First Use

Loosen the battery safety screw and open the battery compartment on the underside of the scale. Remove the plastic isolator tab from the battery compartment (if present) prior to use. Check that the battery is inserted following the correct polarity. Close the battery compartment and tighten the battery compartment screw. Position the scale on a firm flat surface.

#### Changing the Unit

**STEP 1:** Press the centre of the glass platform to turn on the scale.

**STEP 2:** Select kg, st or lbs weight mode by pressing [UNIT] on the base of the scale. The scale will then turn off.

#### Setting the Scale to Zero

**STEP 1:** Press the centre of the glass platform and make sure no weight is on the scale. [0.00/ZERO] will be displayed.

**STEP 2:** The scale will switch off and is now ready for use.

**NOTE: THE ZERO-SETTING PROCESS MUST BE REPEATED IF THE SCALE IS MOVED. AT ALL OTHER TIMES STEP STRAIGHT ONTO THE SCALE.**

#### Using the Body Analyser for a Weight Reading Only

**STEP 1:** Step onto the platform and stand still.

**STEP 2:** After 2–3 seconds the weight reading will be displayed.

**STEP 3:** Step off the platform; the weight reading will be displayed for several seconds before the scale turns off.

#### Setting User Memory

**STEP 1:** Press the centre of the platform and make sure no weight is on the scale. [0.00/ZERO] will be displayed

**STEP 2:** Press and hold [O] for 2 seconds.

**STEP 3:** Once the user number is flashing, select by pressing [>] or [<]. Press [O] to confirm the selection.

**STEP 4:** The male [M] or female [F] symbol will flash. Press [>] or [<] to choose between 'Male' [M], 'Female' [F], 'Male Athlete' [M+\*], 'Female Athlete' [F+\*] then press [O].

**STEP 5:** The height display will flash next. Use [>] or [<] as necessary to select user height, then press [O].

**STEP 6:** To set user age, wait until the age display flashes. Press [>] and [<] as necessary, then press [O] to confirm.

**STEP 7:** The scale will switch off and is now ready for use.

**STEP 8:** To add another user or update the memorised data, repeat steps 1–7.

**ATHLETE MODE: AN ATHLETE IS DEFINED AS A PERSON WHO IS INVOLVED IN INTENSE PHYSICAL ACTIVITY OF APPROXIMATELY 12 HOURS PER WEEK AND WHO HAS A RESTING HEARTBEAT OF APPROXIMATELY 60 BEATS PER MINUTE OR LESS. WHEN ATHLETE MODE IS SELECTED, ONLY AGES 15–85 ARE AVAILABLE. WHEN ATHLETE MODE IS NOT SELECTED, AGES 10–85 ARE AVAILABLE.**

## Using the Body Analyser Scale

This scale uses BIA (Bio Impedance Analysis) technology which passes a small electrical impulse through the body to distinguish fat from lean tissue. The electrical impulse cannot be felt and is perfectly safe. Contact with the body is made via stainless steel pads on the glass platform.

**NOTE: DURING PREGNANCY A WOMAN'S BODY COMPOSITION CHANGES CONSIDERABLY TO SUPPORT THE DEVELOPING CHILD. UNDER THESE CIRCUMSTANCES, BODY FAT PERCENTAGE READINGS COULD BE INACCURATE AND MISLEADING. PREGNANT WOMEN SHOULD THEREFORE ONLY USE THE WEIGHT FUNCTION. FOR MORE INFORMATION, VISIT [SALTER.COM/UNDERSTANDING-YOUR-MEASUREMENTS](http://SALTER.COM/UNDERSTANDING-YOUR-MEASUREMENTS).**

### First Weight and Body Fat Readings

**STEP 1:** Position the scale on a flat, level surface. Press on the platform centre and remove weight from the scale.

**STEP 2:** Select user number by using [➤] or [◀]. Press [O] again to confirm the selection.

**STEP 3:** When [0.00/ZERO] is displayed, with bare feet, step onto the glass platform and stand still. Ensure feet are evenly spaced on the platform.

**STEP 4:** After 2–3 seconds, the weight will be displayed.

**STEP 5:** Remain standing on the scale while other readings are taken.

**STEP 6:** The BMI is displayed followed by body fat %, total body water %, muscle mass and bone mass.

### WEIGHT AND BODY FAT AUTO USER DETECTION

After first use this scale will automatically detect the user for weight and body fat readings.

**STEP 1:** Position the scale on a flat, level surface.

**STEP 2:** With bare feet, step onto the glass platform and stand still. Ensure feet are evenly spaced on the platform.

**STEP 3:** After 2–3 seconds, the weight will be displayed.

**STEP 4:** Remain standing on the scale while other readings are taken.

**STEP 5:** If the user number can be identified, the weight is displayed followed by BMI, body fat %, total body water %, muscle mass and bone mass.

**STEP 6:** If more than one possible previous user number can be identified, the display will show all possible previous user numbers e.g. [P1P4]. Press [◀] to select [P1]. Press [➤] to select [P4].

**STEP 7:** The weight is displayed followed by BMI, body fat %, total body water %, muscle mass and bone mass.

### PREVIOUS WEIGHT & BODY FAT READINGS

Previous weight and body fat readings can be recalled.

**STEP 1:** Position the scale on a flat, level surface.

**STEP 2:** Press the platform centre and remove weight from the scale.

**STEP 3:** [0.00/ZERO] will be displayed.

**STEP 4:** Select a user number by pressing the [➤] or [◀] buttons. Press [O] to confirm the selection.

**STEP 5:** The previous weight and body fat readings will be displayed.

**STEP 6:** If there are no previous weight and body fat readings for the user number selected [----] will be displayed.

### Warning indicators

**Lo** – Replace battery.

---- – Weight exceeds maximum capacity.

**C** – Incorrect zero reading. Repeat setting the scale to zero.

### Specifications

Product code: 9182

Display: LCD

Battery: 1 x 3 V CR2032 coin cell battery (included)

Unit of measurement: kg, st, lbs

### Disposal of Waste Batteries and Electrical and Electronic Equipment

This symbol on the product, its batteries or its packaging means that this product and any batteries it contains must not be disposed of with household waste. Instead, it is the user's responsibility to hand this over to

an applicable collection point for the recycling of batteries and electrical and electronic equipment. This separate collection and recycling will help to conserve natural resources and prevent potential negative consequences for human health and the environment due to the possible presence of hazardous substances in batteries and electrical and electronic equipment, which could be caused by inappropriate disposal. Some retailers provide take-back services which allow the user to return exhausted equipment for appropriate disposal. **It is the user's responsibility to delete any data on electrical and electronic equipment prior to disposal.** For more information about where to drop batteries, electrical and electronic waste off, please contact the local city/municipality office, household waste disposal service, or the retailer.

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**EN** If this product does not reach you in an acceptable condition, please contact our Customer Services Department at **salter.com**. Please have your delivery note to hand as details from it will be required. If you wish to return this product, please return it to the retailer from where it was purchased with your receipt (subject to their terms and conditions).

\*To be eligible for the extended guarantee, go to **guarantee.upgs.com/salter** and register your product within 30 days of purchase.

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©Salter trademark. All rights reserved.
UP Global Sourcing UK Ltd., Victoria Street, Manchester OL9 0DD. UK.