


EN I			
1. Mini Massage Gun	7. Round massage head		
2. Touch button (O)	8. Conical massage head		
3. Charging port	9. U-shaped massage head		
4. Type-C charging cable	10. Flat massage head		
5. LCD screen	11. Carry bag		
6. Charge level indicator/Speed level indicator			

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EN I Please retain instructions for future reference.

Safety Instructions

When using electrical products, basic safety precautions should always be followed. Check that the voltage indicated on the rating plate corresponds with that of the local network before connecting the product to the mains power supply. Children aged from 8 years and people with reduced physical, sensory or mental capabilities or lack of experience and knowledge can use this product, only if they have been given supervision or

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instruction concerning use of the product in a safe way and understand the hazards involved. Children should be supervised to ensure that they do not play with the product. Unless they are older than 8 and supervised, children should not perform cleaning or user maintenance. If the power supply cord, plug or any part of the product is malfunctioning or if it has been dropped or damaged, cease using the product immediately to avoid potential injury. This product contains no user serviceable parts, only a qualified electrician should carry out repairs. Improper repairs may place the user at risk of harm. Do not use the product for anything other than its intended use. Use of an extension cord with the appliance is not recommended. Only use the product with the power supply unit provided.

This appliance should not be operated by means of an external timer or separate remote control system. Do not use this product outdoors, it is intended for indoor use only. Ensure that all hair, clothing and jewellery are kept clear of the massaging mechanisms or any other moving parts at all times. This product is intended for domestic use only. It should not be used for commercial purposes. Do not immerse in any liquid to clean. Never use abrasive cleaners to clean this product. Do not operate in wet or damp conditions. If there are any concerns regarding health or illnesses which may cause concern, consult a doctor before using this product. It is recommended to consult a doctor before using this product if the user:

- Has osteoporosis
- Has heart problems
- Is fitted with a Pacemaker or other electronic medical device
- Is pregnant
- Is diabetic
- Has suffered trauma to the skin

This product is not recommended for use by people with sensory deficiencies. Never cover the product when it is in use. Do not use the product directly on skin. Use over a piece of clothing or towel. Do not use the product for more than 15 minutes at a time and no more than twice a day. Never use the product on swollen or inflamed areas of the skin. This product is a non-professional product designed to provide a soothing massage only. Do not use this product as a substitute for medical attention. Do not use this product before going to bed as the massager has a stimulating effect and may delay sleep. Do not use the product while in bed.

Battery Safety

The supply terminals are not to be short-circuited. In the event of a battery and cell leaking, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash the affected area with copious amounts of water and seek medical advice. Secondary cells and batteries need to be charged before use. After extended periods of storage, it may be necessary to charge and discharge the cells or batteries several times to obtain maximum performance. Keep batteries out of the reach of children. Do not crush, puncture, dismantle or otherwise damage batteries. Do not store above 60 °C. Avoid storage in direct sunlight. Do not dispose of exhausted batteries in household waste. Do not expose the batteries to heat or fire as they may explode. Once the batteries have been exhausted, dispose of the product at a local recycling station. Contact the local government offices for further details. This product is fitted with a 7.4 V 2000 mAh Li-Ion battery (non-removable).

Disposal of Waste Batteries and Electrical and Electronic Equipment

This symbol on the product, its batteries or its packaging means that this product and any batteries it contains must not be disposed of with household waste. Instead, it is the user's responsibility to hand this over to an applicable collection point for the recycling of batteries and electrical and electronic equipment. This separate collection and recycling will help to conserve natural resources and prevent potential hazardous consequences for human health and the environment due to the possible presence of hazardous substances in batteries and electrical and electronic equipment, which could be caused by inappropriate disposal. Some retailers provide take-back services which allow the user to return exhausted equipment for appropriate disposal. **It is the user's responsibility to delete any data on electrical and electronic equipment prior to disposal.** For more information about where to drop batteries, electrical and electronic waste off, please contact the local city/municipality office, household waste disposal service, or the retailer.

Care and Maintenance

Wipe the external surface of the massage gun with a soft, damp cloth and allow to dry thoroughly. Do not use any harsh or abrasive detergents or cleaning agents to clean any part of the massage gun.

DO NOT:

Use any abrasive cleaning solvents to clean the massage gun. Follow the instructions in the section entitled 'Charging the Mini Massage Gun'. Drop the massage gun.

Attempt to disassemble the massage gun.

Dispose of the massage gun in a fire as it will result in an explosion.

Massage Head Attachments

Round Massage Head

For large muscle areas.

Conical Massage Head

For large tight muscles, deep tissue and target areas.

U-shaped Massage Head

For spine, neck and Achilles tendon.

Flat Massage Head

For all areas.

Instructions for Use

Using the Mini Massage Gun

STEP 1: Push the required massage head into the front of the massage gun.
STEP 2: To turn the massage gun on, press 'O' for approx. 3 seconds on the LCD touch screen.
STEP 3: Once the massage gun is switched on, press 'O' once to cycle through the 6 speed/intensity levels.
STEP 4: To turn the massage gun off, press and hold 'O' for 3 seconds on the LCD touch screen.

NOTE: The massage gun is fitted with an auto cut-off feature. After 15 minutes the massage gun will turn off.

WARNING: Long term use of the massage gun can cause the massage gun to overheat. If this occurs, stop using immediately and allow the massage gun to fully cool before using again.

Charging the Mini Massage Gun

STEP 1: Make sure the massage gun is off before charging it.
STEP 2: Insert the Type-C end of the charging cable into the charging port located at the base of the massage gun
STEP 3: Insert the USB end of the charging cable into a suitable USB port to commence charging.

NOTE: The massage gun will take 3–5 hours to charge depending on the current charge. Once charging is complete, please remove the cable from the

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port and store the cable away safely. One full charge will provide 2–4 hours of use dependent on the speed settings used for the massage gun. Press 'O' once and the LCD screen will display the battery percentage; the max. charge the massage gun is able to show is 99 %.

Storage

Check that the massage gun is cool, clean and dry before storing in a cool, dry place.

Troubleshooting

The massage gun will not turn on.

The massage gun is not charged. Follow the instructions in the section entitled 'Charging the Mini Massage Gun'.

The massage gun has turned off.

The massage gun is fitted with an auto-cut off feature. The massage gun will stop after 15 minutes to prevent over-use or injury. Follow the instructions in the section entitled 'Using the Mini Massage Gun' to switch the massage gun back on.

The massage gun will not charge.

The Type-C charging cable hasn't been fully inserted. Follow the instructions in the section entitled 'Charging the Mini Massage Gun'. The massage gun is not compatible with the selected adaptor. Check that the adaptor is compatible and ensure it has a power supply of 5 V === 1 A. The massage gun has overheated. Wait until the massage gun has fully cooled before attempting to charge it.

The massage gun is not pulsing.

The massage gun has low battery.

Follow the instructions in the section entitled 'Charging the Mini Massage Gun'. Too much pressure is being applied with the massage gun. Reduce the amount of pressure being applied.

Specifications

Product code: EE7427	Operating temperature: -10 °C–40 °C
Input: 5 V === 1 A	Noise level < 60 dB
Output: 25 W	Product dimensions (without attachment): approx.13.5 x 13.5 x 4.5 cm

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