SALTER®

SINCE 1760

Ultra Slim Glass Analyser Scale

Instruction Manual

EN | Please retain instructions for future reference.

Safety Instructions

BATTERY WARNING: KEEP OUT OF THE REACH OF CHILDREN.

This product is not a toy.

This product contains no user serviceable parts. Only a qualified electrician should attempt repairs.

If any part of the appliance is malfunctioning or if it has been dropped or damaged, cease using the product immediately to avoid potential injury.

Do not place, store or use the product near direct sunlight or direct heat sources or in a damp/high humidity environment.

This product is a high precision device; handle it with care.

This product is intended for domestic use only.

CAUTION: DO NOT USE THE ANALYSER FUNCTION IF YOU HAVE A PACEMAKER OR OTHER MEDICAL DEVICE FITTED. IF IN DOUBT, CONSULT YOUR DOCTOR.

Battery Safety

This product is supplied with 1 x CR2032 batteries. Only batteries of the same or equivalent type as those recommended should be used with this product.

Batteries are to be inserted with the correct polarity.

Exhausted batteries should be removed from the product to avoid leakage.

WARNING

THIS PRODUCT CONTAINS A BUTTON/COIN BATTERY

Swallowing can cause severe or fatal injuries within 2 hours.

Keep batteries out of reach of children (whether new or used).

If you think batteries may have been swallowed or placed inside any part of the body, seek immediate medical attention.

WARNING: IF THE BATTERY COMPARTMENT DOES NOT CLOSE SECURELY, STOP USING THE PRODUCT AND KEEP IT AWAY FROM CHILDREN.

CAUTION: VIGILANCE MUST BE MAINTAINED WITH "FLAT" OR SPARE BUTTON OR COIN BATTERIES IN THE HOME AND THE PRODUCTS THAT CONTAIN THEM.

There are no obvious specific symptoms when a button or coin battery is stuck in a child's oesophagus (food pipe). However, the child might exhibit the following symptoms:

- Coughing, gagging or drooling
- Appear to have a stomach upset or virus
- Vomiting
- Pointing to their throat or stomach
- Pain in their abdomen, chest or throat
- Tiredness and lethargy
- Be quieter or more clingy than usual
- Loss or reduction of appetite
- Refusal or inability to eat solid foods

These sorts of symptoms may fluctuate.

A specific symptom of button or coin battery ingestion is vomiting fresh (bright red) blood. Seek immediate medical attention if this occurs.

Care and Maintenance

Wipe the exterior of the scale with a soft, damp cloth and allow to dry thoroughly. Never immerse the scale in water or any other liquid.

Do not use harsh or abrasive cleaning detergents or scourers to clean the scale, as this could cause damage.

Before First Use

Loosen the battery safety screw and open the battery compartment on the underside of the scale. Remove the plastic isolator tab from the battery compartment (if present) prior to use. Make sure the battery is inserted following the correct polarity. Close the battery compartment and tighten the battery compartment screw. Position the scale on a firm flat surface. For use on carpet, attach the enclosed carpet feet.

Changing the Units

STEP 1: Press the centre of the glass platform to turn on the scale.

STEP 2: Select kg, st or lb weight mode using the switch on the underside of the scale. The scale will then turn off.

Using the Body Analyser Scale to take Weight Readings

STEP 1: Position the scale on a flat, level surface. Press the centre of the glass platform to turn on the scale.

STEP 2: When [0.00/ZERO] is displayed, step on to the platform and stand still.

STEP 3: After 2–3 seconds the weight reading will be displayed.

STEP 4: Step off the platform; the weight reading will be displayed for several seconds before the scale turns off.

Setting User Memory

STEP 1: Position the scale on a flat, level surface. Press the centre of the glass platform to turn on the scale.

STEP 2: Press [SET].

STEP 3: Once the user number is flashing, select by pressing [^] or [v]. Press [SET] to confirm the selection.

STEP 4: The male or female symbol will flash. Press [^] until the correct symbol is flashing, then press [SET].

STEP 5: The height display will flash next. Press [^] or [v] as necessary to set user height, then press [SET].

STEP 6: To set user age, wait until the age display flashes. Press [^] and [v] as necessary, then press [SET] to confirm.

STEP 7: The scale will switch off and is now ready for use.

STEP 8: To add another user or update the memorised data, repeat steps 1–7.

Using the Body Analyser Scale

This scale uses BIA (Bio Impedance Analysis) technology which passes a small electrical impulse through the body to distinguish fat from lean tissue. The electrical impulse cannot be felt and is perfectly safe. Contact with the body is made via pads under the glass platform.

NOTE: DURING PREGNANCY A WOMAN'S BODY COMPOSITION CHANGES CONSIDERABLY TO SUPPORT THE DEVELOPING CHILD. UNDER THESE CIRCUMSTANCES, BODY FAT PERCENTAGE READINGS COULD BE INACCURATE AND MISLEADING. PREGNANT WOMEN SHOULD THEREFORE ONLY USE THE WEIGHT FUNCTION. FOR MORE INFORMATION, VISIT SALTER.COM/UNDERSTANDING-YOUR-MEASUREMENTS.

Weight and BIA Readings

STEP 1: Position the scale on a flat, level surface. Press on the platform centre and remove weight from the scale.

STEP 2: Press [SET]. Select user number by pressing [^] or [v].

OR

STEP 1: Step onto the platform and stand still.

STEP 2: The display will cycle through users 01 02 03 04. Step off the scale when the desired user number is displayed.

STEP 3: The LCD display will show the saved data for that user, followed by a zero reading.

STEP 4: When [0.00/ZERO] is displayed, with bare feet, step onto the glass platform and stand still. Ensure feet are evenly spaced on the platform.

STEP 5: After 2–3 seconds, the weight, BMI, body fat % and body water % will be displayed.

NOTE: FOR CONSISTENCY, ALWAYS WEIGH ON THE SAME SCALE PLACED ON THE SAME FLOOR SURFACE.

DO NOT COMPARE WEIGHT READINGS FROM ONE SCALE TO ANOTHER AS SOME DIFFERENCES WILL EXIST DUE TO MANUFACTURING TOLERANCES.

PLACING THE SCALE ON A HARD, EVEN FLOOR WILL ENSURE THE GREATEST ACCURACY AND CONSISTENCY.

IT IS RECOMMENDED THAT READINGS ARE TAKEN AT THE SAME TIME OF DAY, PREFERABLY EARLY EVENING AND BEFORE A MEAL FOR THE MOST CONSISTENT RESULTS.

THE SCALE ROUNDS UP OR DOWN TO THE NEAREST INCREMENT. IF WEIGHT IS MEASURED TWICE AND THE READINGS DIFFER BY ONE INCREMENT, THE WEIGHT LIES BETWEEN THE TWO READINGS.

Troubleshooting

Problem	Cause	Solution
"Err2" is showing on the LCD display when taking body fat and body water measurements.	The user does not have bare feet. The condition of the user's skin is affecting the reading.	Remove shoes and socks before using the scale. Wipe feet with a damp cloth, leaving them slightly moist before stepping onto the scale.

Warning indicators

Lo – Replace battery.

Err1 – Unstable weight. Stand still.

Err2 – Out of range, incorrect operation or poor feet contact.

Specifications

Product code: 9141

Display: LCD

Battery: 1 x 3V CR2032 coin cell

battery (included)

Unit of measurement: kg, st, lbs

Disposal of Waste Batteries and Electrical and Electronic Equipment

This symbol on the product, its batteries or its packaging means that this product and any batteries it contains must not be disposed of with household waste. Instead, it is the user's responsibility to hand this over to an applicable collection point for the recycling of batteries and electrical and electronic equipment. This separate collection and recycling will help to conserve natural

resources and prevent potential negative consequences for human health and the environment due to the possible presence of hazardous substances in batteries and electrical and electronic equipment, which could be caused by inappropriate disposal.

Some retailers provide take-back services which allow the user to return exhausted equipment for appropriate disposal. **It is the user's responsibility to delete any data on electrical and electronic equipment prior to disposal.** For more information about where to drop batteries, electrical and electronic waste off, please contact the local city/municipality office, household waste disposal service, or the retailer.

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©Salter trademark. All rights reserved.
UP Global Sourcing UK Ltd.,