(NPWT) Negative Pressure Wound Therapy Pump

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User's Manual
Manual del usuario

Table of Contents

1. Introduction	3
2. Symbol List	4-5
3. Device Specifications	6
4. Indications for Use	7
5. Contraindications for Use	7
6. Warnings	8
7. Precautions	9
7.1) Be aware for any of the following conditions	9
7.2) Prior to Therapy	10
7.3) Periwound Skin	11
8. Features	11
8.1) System overview	11
8.2) Control Panel User Interface	11
8.3) Control Panel User Interface consists of 4 buttons	12
8.4) Alarm Features / Troubleshooting	12
8.4.1 Alarm Table	12-13
8.4.2 LCD Screen in Alarm Condition	13
8.4.3 Mute the alarm	14
9. Instructions for Use	14
9.1) Canister Application	14
9.2) Operation the Device	14
9.2.1 Power the VCare1000-300S NPWT On or Off	14
9.2.2 Screen Lock	15
9.2.3 Battery Level Indicator	15
9.2.4 Charging Indicator	15
9.3) Clinical Settings	16
9.3.1 Clinical Mode Settings	16
9.3.2 Pressure and Time Settings	16
9.3.3 Therapy Device Mode	17
10. Disposal	17
11. Maintenance and Replacement Parts	17
12. Cleaning	18
13. Warranty Information	18
Appendix 1	19
Appendix 2	20-23

WARNING

Indications, Contraindications, Warnings, Precautions and other Safety Information are contained in this User Manual.

To reduce risk of serious or fatal injury, all caregivers and patients must carefully read and follow all user instructions and safety information.

1. Introduction

The Dynarex Negative Pressure Wound Therapy Pump is a portable, battery-powered system intended to generate negative pressure or suction to remove wound exudates, infectious material, and tissue debris from the wound bed to help promote wound healing. The system consists of one pump and one AC power cable. Wound dressings, additional canisters (with 140 mL, 400 mL and 600 mL options), carrying cases, and other accessories are sold separately. The Negative Pressure Wound Therapy Pump, when used with the proprietary wound dressing, creates a negative pressure environment.

The pump and the wound dressing are designed to produce a negative pressure environment in continuous or intermittent mode. This allows the user to program the specific pressure ranging from -20mmHg to -200mmHg. In continuous mode, the pressure is applied to the wound as long as the pump is powered on. In variable pressure mode, the pump provides gentle intermittent pressure, eliminating the painful high and low spikes characteristic of traditional intermittent pressure.

Clinically suggested wound types that can be treated using negative pressure wound therapy (NPWT) are:

- Full-thickness wounds
- Chronic wounds that have failed to heal despite good wound care, including venous and arterial ulcers, diabetic foot ulcers, neuropathic ulcers, pressure ulcers, and other chronic wounds
- Acute wounds such as dehiscent surgical wounds, traumatic wounds, and other open ulcers which would benefit from NPWT are also appropriate.

2. Symbol List

	General Mandatory Action Sign
	Class II Equipment
IP22	Ingress Protection: - Protected against access to hazardous parts with a finger - Protected against solid foreign objects of 12.5 mm (Z and greater) - Protected against vertically falling water drops when enclosure is tilted up to 15°
	Single Use Only
	Consult Instructions for Use
	Tripping Hazard
	Type BF Applied Part
	Mute the alarm (button in the front panel), or status of alarm being muted (symbol shows in LCD)
	Temperature Limits
	Humidity Limitation
	Keep Dry
	Caution: Federal law restricts this device to sale by or on the order of a licensed health care provider.

2. Symbol List (continued)

	Waste Electrical Goods Recycled
	No Open Flame, No Open Ignition Source, and Smoking Is Prohibited
	Unique Device Identification System
	Batch Code
	Date of Manufacture
	Serial Number
	Continuous Mode
	Intermittent Mode

3. Device Specifications

Operating Data	Weight (with empty canister attached)	11.3 oz. (320 g)
	Dimensions	Length: 3.58" (9.1 cm) Height: 5.04" (12.8 cm) Width: 2.24" (5.7 cm)
	Mobility	Portable. Carrying case available
	Canister	Volume: 140 mL, 400 mL and 600 mL (optional)
	Therapy Device Modes	Continuous mode Intermittent mode
	Therapy Pressure Options	-20mmHg to -200mmHg
Electrical Data	External Power Supply Input	AC100V-240V, 50-60HZ 0.5A/Voltage tolerance $\pm 10\%$
	External Power Supply Output	DC5V ± 0.25 , 1.5A
	Battery Type	Lithium (rechargeable), 1500mAh/3.7V
	Battery Run Life	72 hours, depending on settings
	Recharging Time	<5 hours from a fully discharged status
IEC Classification	Protection Against Electrical Shock	Class II
	Applied Part	Type BF
	Ingress Protection	IP22
	Equipment not suitable for use in presence of flammable anesthetic mixture with air, oxygen, or nitrous oxide	
Storage Conditions	Temperature Range	-10°C to 45°C (14°F to 113°F)
	Relative Humidity Range	15% to 93% non-condensing
Operating Conditions	Temperature Range	5°C to 40°C (41°F to 104°F)
	Relative Humidity Range	30% to 75% non-condensing
	Atmospheric Pressure Range	700hPa to 1060hPa

4. Indications for Use

The Negative Pressure Wound Therapy Pump is a portable, battery-powered system intended to generate negative pressure or suction to remove wound exudates, infectious material, and tissue debris from the wound bed to help promote wound healing.

5. Contraindications for Use

The pump should **NOT** be used in the following conditions:

- Exposed vessels, organs, or nerves
- Exposed anastomotic sites
- Untreated osteomyelitis
- Exposed vasculature, arteries, or veins in a wound. All exposed vessels and organs in and around the wound must be completely covered prior to initiating NPWT.
- **Note:** A thick layer of natural tissue is preferred. Several layers of fine meshed non-adherent material or bio-engineered tissue may be an alternative. Ensure that protective materials will maintain their position throughout therapy.
- Non-enteric and unexplored fistulas
- Actively bleeding wounds
- Malignancy in the wound. The pump may be used after excision of malignancy if all margins are clear.
- Excess amount of necrotic tissue with eschar
- **Note:** The pump may be used with a minimal amount of slough or after debridement of necrotic tissue.
- Poor blood flow to the wound
- Wounds that are too large or too deep to be accommodated by the dressing
- Patient is unable to follow up with a medical professional or keep scheduled appointments
- Allergy to urethane dressings and adhesives
- Use of topical products that require more frequent applications than the dressing change schedule allows

6. Warnings

- 1) **WARNING:** Review this User's Manual prior to using the pump to reduce risk of serious or fatal injury.
- 2) Do NOT use the pump around explosive or flammable material.
- 3) It is important to note that the pump cannot be taken into a magnetic resonance imaging (MRI) suite or hyperbaric oxygen therapy (HBO) chamber.
- 4) This device should be used only under the direction and/or supervision of a trained professional, such as a doctor or nurse. Performing therapy applications without the express direction and/or supervision of a trained professional may lead to improper product performance and the potential for serious or fatal injury.
- 5) Larger canister sizes (600 mL or larger) should only be used in a facility where drainage can be closely monitored.
- 6) Negative pressure wound therapy has not been cleared for use on children.
- 7) Use a properly rated charger to charge the lithium battery. Incorrect voltage and/or current can cause fire.
- 8) Do not place this device at temperatures greater than 170°F or 76.7°C for more than 2 hours, as this may cause a battery fire.
- 9) If battery swells, gets hot, or smokes while charging, disconnect the charger immediately. This may cause the battery to leak, and the reaction with air may cause the chemicals to ignite, resulting in fire.
- 10) This device shall not be serviced or maintained while in use with a patient.
- 11) Do not tilt the pump more than 90 degrees with respect to the vertical for more than 10 seconds.
- 12) Do not modify the device.
- 13) Operator should keep the device within reach when in use.

7. Precautions

7.1) Be aware of any of the following conditions:

There are additional conditions to take into account before using the pump, such as:

1. **BLEEDING:** There is a risk of bleeding/hemorrhaging with negative pressure wound therapy. If hemostasis cannot be achieved, if the patient is on anticoagulants or platelet aggregation factors, or if the patient has friable blood vessels or infected vascular anastomosis, he or she may have an increased risk of bleeding. Accordingly these patients should be treated in an inpatient care facility per their treating physician. If active bleeding develops suddenly or in large amounts during therapy, immediately disconnect the pump, leave the wound dressing in place, and take measures to stop bleeding. Seek medical attention immediately.
2. **VESSEL AND BONE PROTECTION:** Precautionary measures should be taken if any bones, vessels, ligaments, or tendons are exposed. Additionally, sharp edges (due to bone fragments) require special attention. These areas should be covered and smoothed wherever possible. These conditions should be factored into the therapy prescription as the attending clinician sees fit.
3. **ENVIRONMENT:**
 - a. Defibrillation: Remove the dressing if defibrillation is required in the area of dressing placement. Failure to remove the wound dressings may inhibit transmission of electrical energy and/or patient resuscitation.
 - b. Magnetic Resonance Imaging (MRI): The Negative Pressure Wound Therapy Pump is unsafe in the MR environment. Do not take the Negative Pressure Wound Therapy Pump into the MR environment. However, wound dressings can typically stay on the patient with minimal risk in an MR environment, assuming that the use of the pump is not interrupted for more than two hours.
 - c. Hyperbaric Oxygen Therapy (HBO): Do not take the Negative Pressure Wound Therapy Pump into a hyperbaric oxygen chamber. After disconnecting the device, either (i) replace the dressing with another HBO-compatible material during the hyperbaric treatment, or (ii) cover the unclamped end of the tubing. For HBO therapy, the tubing must not be clamped. Never leave a dressing in place without active negative pressure wound therapy for more than two hours.

4. **INFECTION:** Infected wounds and osteomyelitis pose significant risks for Negative Pressure Wound Therapy. If untreated osteomyelitis is present, therapy should not be initiated. If the wound is infected, it should be closely monitored and dressings should be changed frequently.

Additionally, to reduce the risk of transmission of infectious agents, standard precautions should be taken when handling or working with therapeutic parts or equipment.

5. **PATIENT SIZE AND WEIGHT:** Patient size and weight should be taken into account when prescribing therapy. In addition, small adults, young adults, and elderly patients should be closely monitored.

6. **SPINAL CORD INJURY:** If a patient experiences autonomic dysreflexia (sudden changes in blood pressure or heart rate because of sympathetic nervous system stimulation) discontinue therapy to minimize sensory stimulation and give immediate medical assistance.

7. **MODE:** In unstable anatomical structures, continuous rather than intermittent therapy is recommended to help minimize movement and instability. Continuous therapy is also recommended in patients with an increased bleed risk, profusely exuding wounds, fresh grafts and/or flaps, and wounds with acute enteric fistulae.

8. **ENTERIC FISTULAS:** Wounds with enteric fistulas require special consideration to be effective in negative pressure wound therapy. If enteric fistula effluent management or containment is the only goal of such therapy, negative pressure wound therapy is not recommended.

9. **CIRCUMFERENTIAL DRESSING:** Do not use circumferential dressings.

10. **BRADYCARDIA:** Avoid placement of the Negative Pressure Wound Therapy Pump next to the vagus nerve to minimize the risk of bradycardia.

If any of this information is not understood, contact the manufacturer before using this device.

7.2) Prior to Therapy


- Patient should be assessed, and measures should be taken to optimize and stabilize their medical condition. Nutrition, medication, blood glucose, blood pressure, and circulation as well as other medical issues should be addressed.
- The wound should be recently debrided by whatever measure is appropriate and the amount of necrotic tissue should be minimized.
- Issues of infection should be addressed.

7.3) Periwound Skin

- Ensure that the skin that will be under the dressing is clean, dry, and free of surfactants and oil. Any hair should be clipped.
- The periwound area should be cleaned and allowed to air-dry. The use of a skin preparation wipe is also recommended.
- A thin film dressing or hydrocolloid may be used as additional protection.
- Monitor skin for signs of irritation or breakdown. Treatment may be discontinued if this occurs and cannot be managed.

8. Features

8.1) System Overview

8.2) Control Panel User Interface

8.3) Control Panel User Interface consists of 4 buttons:

Button	Function
	• Turns the system on or off (extended press) • Confirms the programmed pressure or clinical mode • Turns the pump on or off • Changes the language (Press  and  simultaneously)
	• Enters into pump setting menu (extended press) • Increases the programmed pressure • Increases the programmed time in intermittent mode setting • Changes clinical mode in setting menu • Locks or unlocks the screen (Press and hold  and  simultaneously)
	• Reduces the programmed pressure • Reduces the programmed time in intermittent mode setting • Locks or unlocks the screen (Press and hold  and  simultaneously)
	• Mutes the alarm • Changes the language (Press  and  simultaneously)

8.4) Alarm Features/Troubleshooting:

8.4.1) Alarm Table

Alarm Type	Cause/Description	Audio Alarm Features	Visual Alarm Features	System Status	Suggested Mitigation
Low Battery Alarm	Less than 25% power remaining (3.7V) indicating the system will power off soon.	1 beep every 20 seconds	Alarm symbol flashes on LCD screen. Yellow LED remains illuminated.	Pump continues to function until Under Voltage Alarm sounds.	Plug the Negative Pressure Wound Therapy Pump in, allowing it to function and charge simultaneously.
Under Voltage Alarm	Less than 10% power remaining (3.5V). The pump will power off automatically.	1 beep every 20 seconds	Alarm symbol flashes on LCD screen. Yellow LED remains illuminated.	Pump shuts off in 10 seconds.	Plug the Negative Pressure Wound Therapy Pump in, allowing it to function and charge simultaneously.

8.4.1) Alarm Table (continued)

Alarm Type	Cause/Description	Audio Alarm Features	Visual Alarm Features	System Status	Suggested Mitigation
Air Leakage Alarm	Pump unable to reach 80% of the preset pressure after 2 minutes of pump effort	No beep when battery is full; 1 beep every 20 seconds when battery is less than 50% charged.	Alarm symbol flashes on LCD screen. Yellow LED remains illuminated.	Pump will stop running in 10 minutes.	Verify connection between 1) canister and wound dressing; 2) wound and wound dressing 3) canister and pump. If necessary, power off and back on to restart the pump after adjustment. Install a new canister.
Canister Full Alarm	Alarm is activated if the canister is almost full of exudates after 2 minutes of pump effort	1 beep every 20 seconds	Alarm symbol flashes on LCD screen. Yellow LED remains illuminated.	Pump will stop running immediately.	Install a new canister
Blockage Alarm	Alarm is activated if there is a pipe blockage after 5 minutes of pump effort.	1 beep every 20 seconds	Alarm symbol flashes on LCD screen. Yellow LED remains illuminated.	Pump remains in function.	Check tubing between wound dressing and canister to determine if tube is blocked by impurities or tube kinking.
Pump Idle Alarm	No action after 5 minutes.	1 beep every 20 seconds	Alarm symbol flashes on LCD screen. Yellow LED remains illuminated.	Pump remains in current status.	Press any button to lift alarm

8.4.2) LCD Screen in Alarm Condition

 Low Battery Alarm	 Under Voltage Alarm	 Air Leakage Alarm
 Canister Full Alarm	 Blockage Alarm	 Pump Idle Alarm

 An Alarm will be displayed on the LCD screen User Interface, with a delay of about 10 seconds, when the Negative Pressure Wound Therapy Pump detects a condition that requires immediate patient or caregiver attention in order to ensure the prescribed therapy is being delivered.


8.4.3) Mute the Alarm

1. Press  to mute the alarm.  Symbol will be shown on the LCD screen.
- 2 The alarm audio will resume in 5 minutes.

9. Instructions for Use

9.1) Canister Application

1. Install one canister
2. Connect the canister tubing

 Excess tubing may present a tripping hazard. Ensure that excess tubing is out of areas where people may walk.

9.2) Operating the Device

9.2.1) Power the system On or Off


1. Press and hold  to turn the system ON or OFF.

System On


1. The system will go through a self-check and then presents a “welcome” screen


The next interface shows the version number of the Pump

2. The "home screen" then will be shown (The setting mode is continuous mode)

The "home screen" then will be shown (The setting mode is intermittent mode)

3. Press  to enter Settings menu;

Press  to start therapy by continuing the previous clinical mode

9.2.2) Screen Lock

1. If there is no operation for 5 minutes, the screen will be locked.

2.  Symbol will be shown on the screen.

3. Unlock/lock the screen by pressing and holding  and  simultaneously.

9.2.3) Battery Level Indicator

Low Battery <25%
Flashing

In Use 25%~50%

In Use >50%

Fully Charged No flash

9.2.4) Charging Indicator

a) Battery charging, icons of four levels flash circularly, LED GREEN

b) Battery fully charged , no flashing

Power cords may present a tripping hazard. Ensure that all cords are out of areas where people may walk. Power cords may result in strangulation. Ensure that all cords are out of breathing system hoses areas.


9.3) Clinical Settings

CAUTION NOTICE TO PATIENTS: This section contains specific information intended for clinicians ONLY. Do not attempt to perform any applications or setting adjustments in this section without the express direction and supervision of your treating physician.

9.3.1) Clinical Mode Settings

When system is in the setting menu:

1. Program the clinical mode by pressing .
2. Press  to confirm the programmed Clinical Mode.


9.3.2) Pressure and Time Settings

WARNING: Negative pressure at the wound may be below set pressure, potentially compromising therapeutic benefits.



1. If Continuous Mode is confirmed (see 10.3.1), set the pressure and time:

- a) Program the pressure value by pressing  or .
- b) Confirm the programmed pressure value by pressing .



2. If Intermittent Mode is confirmed (see 10.3.1), set the pressure and time:

- a) Program the **pressure** value by pressing  or .
- b) Confirm the programmed **pressure** value by pressing .
- c) Program the **operating time** by pressing  or .
- d) Confirm the programmed **operating time** by pressing .


9.3.3) Therapy Device Mode


If Clinical mode setting, pressure, and time settings are confirmed, Press  to start therapy.

Below picture shows the operating interface:

Continuous Mode

Intermittent Mode

10.) Disposal

The Negative Pressure Wound Therapy Pump is powered electromechanically by a battery, which should be recycled according to the local regulations governing such products and Waste Electrical and Electronic Equipment (WEEE) Directive.

The wound dressings, tubing, and canister can be disposed of according to policy for wound care dressings after use.

11.) Maintenance and Replacement Parts

The Negative Pressure Wound Therapy Pump contains no user serviceable parts inside: **Opening or tampering with this device will void the warranty.** In the event the system requires repairs, it should be returned to the medical equipment company.

Power Adapter: the system should only be recharged using the AC/DC adapter provided or an equivalent IEC60601-1 compliant adapter with +5VDC±0.25/1.5A output.

Battery: Do not attempt to open, disassemble, or service the battery pack. Do not crush, puncture, short external contacts, or dispose of in fire or water. Use only a manufacturer-approved battery. If the device will not be in use for an extended period of time, the battery should be maintained by recharging regularly. Battery should be stored in a safe and dry place.

12. Cleaning

The pump can be wiped as necessary with a damp cloth using a mild household cleaner. If additional information is required, please contact your local representative or the manufacturer.

13. Warranty

LIMITED LIFETIME WARRANTY

Your Dynarex Product is warrantied to be free of defects in materials and workmanship for the lifetime of the product for the original purchaser. This item was built to exacting standards and carefully inspected prior to shipment. This Limited Lifetime Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

The Warranty shall not apply under the following conditions:

- Problems arising from normal wear
- Problems arising from failure to adhere to the product instructions
- Problems arising from misuse, negligence, accident or improper operation, maintenance or storage
- Problems arising from modifications or unauthorized repairs, parts or attachments
- Products where the serial number has been removed or defaced
- Problems with non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement

Dynarex shall not be liable for any consequential or incidental damages whatsoever. Dynarex shall repair or replace defective products at its option. The foregoing warranty is exclusive and in lieu of other express warranties, if any, including the implied warranties of merchantability and fitness of a particular purpose. The remedy for any violation of the implied warranty shall be limited to repair or replacement of the defective product pursuant to the terms contained herein.

If you have a question about your Dynarex device or this warranty, please contact an authorized Dynarex dealer.

Appendix 1

Product Classification (according to IEC 60601):

- According to the type of protection against electrical shock, this device is classified as class II Equipment, and Type BF Equipment that is powered by an external electrical power source.
- According to the degree of protection against harmful ingress of water this System is classified as water Protected Equipment (IP22: Protection against vertically falling water drops.)
- **CAUTION:** This device has been tested and confirmed to comply with the IEC60601-2:2020 and essential requirements of Medical Device Directive 93/42/EEC. However with the proliferation of radio-frequency transmitting equipment, and other sources of electrical noise in a healthcare environment, high levels of interference may induce an abnormal stoppage or other disruption of this device. This device may also cause adverse effects in other nearby equipment. It is strongly recommended that this device be isolated from other electromagnetic equipment when in use.
- This System is classified as Equipment not Suitable for use in the presence of a Flammable Anesthetic Mixture with Air or Oxygen or Nitrous Oxide.
- According to the mode of operation this System is classified as Equipment that can be used for Continuous Operation.
- Unit is packaged for transportation by common carrier.

Appendix 2

EMC Information

Guidance and Manufacturer's Declaration - Electromagnetic Emissions - for all ME EQUIPMENT and ME SYSTEM.

Table 1: Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The Negative Pressure Wound Therapy Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the pump should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Negative Pressure Wound Therapy Pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Negative Pressure Wound Therapy Pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Compliance	

Guidance and Manufacturer's Declaration - Electromagnetic Emissions - for all ME EQUIPMENT and ME SYSTEM.

Table 2: Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The Negative Pressure Wound Therapy Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Negative Pressure Wound Therapy Pump should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge(ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	Electrostatic Discharge(ESD) IEC 61000-4-2	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment

Table 2: Guidance and Manufacturer's Declaration – Electromagnetic Immunity (cont'd)

The Negative Pressure Wound Therapy Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Negative Pressure Wound Therapy Pump should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Surge IEC 61000-4-5	± 1 kV lines to lines ± 2 kV lines to earth	± 1 kV lines to lines ± 2 kV lines to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ ($> 95\%$ dip in U_T) For 0,5 cycle 40% U_T (60% dip in U_T) For 5 cycle 70% U_T (30% dip in U_T)	$< 5\% U_T$ ($> 95\%$ dip in U_T) For 0,5 cycle 40% U_T (60% dip in U_T) For 5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Negative Pressure Wound Therapy Pump requires continued operation during power mains interruptions, it is recommended that the Negative Pressure Wound Therapy Pump be powered from an uninterruptible power supply or a battery.
	For 25 cycle $< 5\% U_T$ ($> 95\%$ dip in U_T) For 5 cycle	70% UT (30% dip in U_T) For 25 cycle $< 5\% U_T$ ($> 95\%$ dip in U_T) For 5 cycle 3A/m	
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
Note: UT is the a.c. mains voltage prior to application of the test level.			

Table 3: Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The Negative Pressure Wound Therapy Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Negative Pressure Wound Therapy Pump should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80MHz 3 Vrms 80MHz to 2.5GHz	3V 3V/m	Portable and mobile RF communications equipment should be used no closer to any part of the system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the ME EQUIPMENT and ME SYSTEM – for ME EQUIPMENT and ME SYSTEM that Are Not LIFE-SUPPORTING.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Negative Pressure Wound Therapy Pump

The Negative Pressure Wound Therapy Pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Negative Pressure Wound Therapy Pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Negative Pressure Wound Therapy Pump as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80 MHz $d = [\frac{3.5}{V_1}] \sqrt{P}$	80MHz to 800MHz $d = [\frac{3.5}{E_1}] \sqrt{P}$	800MHz to 2.5GHz $d = [\frac{7}{E_1}] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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Dynarex Corporation