

Suggested Care and Cleaning

Caron II, Excel, Perval, Sonata, Beth & Staph-Chek Pillows

Pillows should be treated in the same manner as the mattress at patient changes.
LAUNDERING IS NOT RECOMMENDED!

- **ROUTINE CLEANING:** Soils or stains can be removed with neutral soap and lukewarm water. Rinse and dry after cleaning to remove residual soap is recommended.
- **DISINFECTION:** Do not exceed manufacturer's recommended dilution for disinfectants.
- **HARD TO CLEAN SPOTS:** Clean with soft sponge using liquid household type cleaners.
- **TO REMOVE BLOOD:** Use routine cleaning procedure. If needed, use cold water and hydrogen peroxide or cold water and ammonia.
- **DO NOT USE HARSH CLEANERS OR SOLVENTS:**

Disinfectant products including a two-step cleaner/disinfectant process can be used on pillows when following manufacturers' recommended concentration or instructions. When used properly, disinfectants typically pose no real problem to the expected use life of the pillow ticking. However, disinfectants applied at full concentration or in highly concentrated solutions may eventually cause problems and prove harmful to the fabric. It is recommended that the facility perform a test trial with their specific disinfectant to ensure damage does not occur to the pillows.

Note: Standard Textile does not endorse any specific chemical manufacturer or disinfectant.