



Oxivir® HC Disinfectant Cleaner

Oxivir® HC - Engineered for Compatibility

A ready-to use disinfectant cleaner, designed using the power of AHP® for protecting patients, staff and hospital assets.

- One-step cleaner disinfectant powered by AHP® technology
- *Virucide, Bactericide, Sanitizer, Oxidizing Agent
- Effective against key pathogens – MRSA, Norovirus and VRE
- Kills HBV, HCV and HIV in 1 minute
- No VOCs, no added fragrance and no NPEs. Environmentally responsible – active ingredient breaks down to oxygen and water
- Engineered for Compatibility and for use on hard, nonporous surfaces. Compatible with healthcare equipment, computers, screens, plastics, gloves and other surfaces commonly found in healthcare settings
- Non-irritating to eyes and skin

Applications

- Clean and disinfect high-touch surfaces in just ONE Minute
- Compatible with most hard, non-porous surfaces and is a 30-second soft-surface sanitizer





Oxivir[®] HC Disinfectant Cleaner

Oxivir[®] HC - Engineered for Compatibility

Use Instructions: As a One-Step Cleaner Disinfectant

- Pre-clean heavily soiled areas.
- Apply product by coarse trigger spray to hard, non-porous inanimate surfaces.
- All surfaces must remain visibly wet for 1 minute.
- Air Dry, wipe surfaces to dry and remove any residue, or rinse with potable water as necessary.

Note: Rinsing is not necessary unless floors are to be coated with finish or restorer. All food contact surfaces must be rinsed with potable water. Do not use on glassware, utensils, or dishes. Prior to first use on a particular surface, it is recommended to test in a small, discreet area if there is concern of incompatibility, or contact a Diversey representative.

To Spot Sanitize Soft Surfaces

- Remove gross filth prior to sanitization.
- Hold bottle upright 6" – 8" from surface.
- Spray until fabric is wet. **DO NOT SATURATE.** Allow surface to remain wet for 30 seconds.
- Gently blot with a clean, damp, color safe cloth and/or let air dry. For difficult odors or heavy fabrics, repeat application.

Technical data	Oxivir [®] HC Disinfectant Cleaner
EPA Reg. No.	70627-79
Color	Colorless, liquid
pH	3
Scent	No fragrance added
Shelf Life	2 Years

Product	Pack size	Dilution	Product code
Oxivir [®] HC Disinfectant Cleaner	1 x 5 gal / 18.9 L pail	Ready-to-use	101104467 
Oxivir [®] HC Disinfectant Cleaner	12 x 32 oz. / 946 mL	Ready-to-use	101103716 

Safe handling

Please make sure your employees read and understand the product label and Safety Data Sheet before using this product. The label contains directions for use; and both the label and SDS contain hazard warnings, precautionary statements and first aid procedures. SDS are available online at www.diversey.com or by calling 888.352.2249. Improper use or dilution may result in damage to surfaces and may result in health and physical hazards that match those of the concentrate.

AHP[®] and Design are trademarks of Diversey, Inc.