

INVACARE®
SEMI-ELECTRIC
BED



The **Invacare® Semi-Electric Bed** combines effortless positioning of the upper body and knees with the economy of manual bed height adjustment. The pendant control provides motorized positioning of the upper body and/or knees. A manual crank located at the foot end of the bed adjusts bed frame height for caregiver assistance.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

Features

High impact bed end panels are durable, impact and scratch resistant; washable for easy cleaning

Universal bed end can be used as either the headboard or footboard on the bed

Rubber caster tread for better tracking with non-marring material

DC low voltage, compact motor system

Quiet motor operation and fast cycle time on lowering height of bed for efficient bed set-up

Motor system with waterproof enclosures is rated IPX6. Washable with standard water pressure. The bed must not be power washed

Color-coded connections on junction box for easier maintenance

Double-insulated junction box helps provide extra protection against electrical shock

Ergonomic hand pendant comes with durable strain relief and heavy-duty cord. Waterproof casing for easy cleaning, conveniently attaches to bed sheet with clip or rail with hand pendant cradle (sold separately)

Hand pendant cradle available for convenient storage of hand pendant (model # 26975)

Interchangeability simplifies storage, reduces inventory complexity

HCPCS Code: E0295



1. Patented universal bed ends
2. Ergonomic hand pendant
3. Color coded connections
4. Emergency manual crank
5. Thermal overload protection
6. Double insulated electrical components/wiring

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2021
Invacare Corporation.
All rights reserved.
Trademarks are
identified by the symbols
™ and ®. All trademarks
are owned by or licensed
to Invacare Corporation
unless otherwise
noted. Specifications
are subject to change
without notification.
Form No. 10-124
Rev. 03/21
210351



Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
5310IVC	88" L 15"-23" H 36" W	80" L 36" W	Safe Working Load: 450 lb.* Max. User Weight: 350 lb.	N/A	N/A	UL962 CSA Certified	5 Years Welds 2 Years Mechanical/ Electrical
5000IVC Headspring				45 lb.	43" L 8" W 36" D		5 Years Welds 2 Years Mechanical/ Electrical
5890IVC Footspring				77 lb.	47" L 9" W 36" D		5 Years Welds 2 Years Mechanical/ Electrical
5301IVC Universal Bed End				55 lb.	39" L 7" W 24" D		

*Includes user, mattress and all accessories.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.