

INVACARE®
HOMEFILL®
CYLINDERS



Invacare offers a wide variety of **Invacare® HomeFill® Cylinders** to fit your client's lifestyle and daily oxygen needs. These cylinders range in size from as small as a water bottle to large enough to last over 15 hours*.

Now HomeFill cylinders are also available with CGA870 compatible post valves and integrated conservers. With a full line of cylinders, the Invacare HomeFill Oxygen System meets the needs of most oxygen patients.

*Duration Varies Depending on the User's Oxygen Needs.

Invacare® HomeFill® Oxygen Cylinders

Integrated Converter Cylinders

HCPSC Code: K0738*

Continuous Flow Cylinder

HCPSC Code: K0738*

M2 Cylinder

ML4 Cylinder

ML6 Cylinder

M9 Cylinder

D-Tank



Model #	HF2PCM2KIT	HF2PCL4	HF2PCL6	HF2PC9	IOH2MD
---------	------------	---------	---------	--------	--------



Height	8.25"	9"	11"	14"	19.7"
--------	-------	----	-----	-----	-------



Weight	2 lb	3.1 lb	3.5 lb	4.8 lb	7.0 lb
--------	------	--------	--------	--------	--------



Diameter	3"	4.38"	4.38"	4.38"	4.38"
----------	----	-------	-------	-------	-------



Capacity	57 Liters	113 Liters	170 Liters	255 Liters	425 Liters
----------	-----------	------------	------------	------------	------------

Conserved Duration (hrs. @ 20BPM)



Pulse Setting 1	3.5 Hrs	6.9 Hrs	10.4 Hrs	15.7 Hrs	N/A
-----------------	---------	---------	----------	----------	-----



Pulse Setting 2	2.2 Hrs	4.4 Hrs	6.7 Hrs	10.1 Hrs	N/A
-----------------	---------	---------	---------	----------	-----



Pulse Setting 3	1.6 Hrs	3.3 Hrs	5.0 Hrs	7.5 Hrs	N/A
-----------------	---------	---------	---------	---------	-----



Pulse Setting 4	1.3 Hrs	2.6 Hrs	4.0 Hrs	6.0 Hrs	N/A
-----------------	---------	---------	---------	---------	-----



Pulse Setting 5	1.1 Hrs	2.1 Hrs	3.2 Hrs	4.9 Hrs	N/A
-----------------	---------	---------	---------	---------	-----



Duration @ 2 LPM Flow	24 min.	55 min.	75 min.	125 min.	3.5 hrs.
-----------------------	---------	---------	---------	----------	----------



Fill Time	24 min.	55 min.	75 min.	125 min.	3.5 hrs.
-----------	---------	---------	---------	----------	----------

Conserved Duration (hrs. @ 20 BPM)

1142661 Conserver Pulse Setting	1	2	3	4	5	Cylinder Volume
Pulse Volume Range for Setting (mL)	10-17	18-24	25-31	32-38	39-47	N/A
Pulse Volume Average for Setting (mL)	13.5	21	28	35	43	N/A

Post Valve Cylinders

HCPCS Code: K0738*

Post Valve ML6 Post Valve M9 Post Valve D




Model # HF2POST6 HF2POST9 HF2POSTD

	Height	12.0"	15.0"	20.7"
	Weight	3.2 lb	4.1 lb	6.0 lb
	Diameter	4.38"	4.38"	4.38"
	Capacity	170 Liters	255 Liters	425 Liters
	Duration @ 2 LPM Flow	75 min.	125 min.	3.5 hrs.
	Fill Time	75 min.	125 min.	3.5 hrs.

Post Valve Replacement Kits

HCPCS Code: K0738*



Non-compatible
HomeFill Cylinder

HF2POST-R
HF2POSTD-R



*All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

Invacare Corporation

www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2018

Invacare Corporation.

All rights reserved.

Trademarks are identified

by the symbols [™] and

®. All trademarks are

owned by or licensed

to Invacare Corporation

unless otherwise noted.

Specifications are subject

to change without

notification.

Form No. 18-199

180719



Yes, you can.®