

ANALYSIS



Step	Procedure	Result
1	Sample	
2	Solid	
3	Liquid	
4	Dissolve in water	
5	Filter	
6	Filter residue	
7	Filterate	
8	Dry and weigh	
9	Ignite and weigh	
10	Residue	
11	Loss on ignition	
12	Evaporate to dryness	
13	Ignite and weigh	
14	Residue	
15	Loss on ignition	
16	Distill off water	
17	Residue	
18	Dry and weigh	
19	Ignite and weigh	
20	Residue	
21	Loss on ignition	