

INVACARE BEDS FOR POST-ACUTE CARE

A bed for every facility.

Quality and comfort for every resident.



Great looks and great performance? Get the best of both worlds with Invacare beds for Post-Acute Care. From the proven durability of Invacare DLX Series beds, to the flagship Carroll CS Series™ portfolio, to the stylish new Invacare® ETUDE™ HC Bed, Invacare is your first choice for the quality, versatility and durability needed in Post-Acute Care.



Yes, you can.

INVACARE ETUDE HC BED

European Style at Home in Post-Acute Care.

The new Invacare ETUDE HC Bed not only has the styling of a fine Scandinavian bed, but also packs impressive features into its attractive design. From Semi-Trendelenburg tilt to innovative Shear Management Technology, the ETUDE HC Bed is at home in post-acute care settings.



- Shear Management Technology is designed to help prevent harmful shear as the bed is articulated, aiding in the reduction of risk for pressure injuries
- Semi-Trendelenburg 11° tilt for respiratory, digestive or circulatory therapy
- Easy to use 8-function hand pendant with lock-out capability
- Dual height brackets allow bed to go as low as 7.8"



ETUDE BED ACCESSORIES



1. Verso Side Support Rails

Designed to provide a unique collapsable, folding solution, accommodating the needs of those at a higher risk of falls or requiring enhanced care.

2. Lifting Support Trapeze Kit

Provides safe and convenient overhead repositioning assistance while in bed.

3. Side Support Rails

Provide accommodations for safe positioning and transfer support.

4. Transportation and Storage Adaptor Kit

Storage and delivery management of the ETUDE HC Bed is made simple,

5. Emergency Battery Back-up kit

Add an additional level of safety and comfort. The AH rechargeable battery provides a temporary power supply enabling the user or caregiver to safely re-position the bed in emergency situations.

ETUDE HC	
Dimensions	36" (W) X 82" (L)
Deck-to-Floor Height	7.9"- 26.4"
Weight Capacity	440 lbs.
Warranty	Electronics - 2 Years • Actuators - 3 Years • Frame- 5 Years
Deck	Slat Deck
Mobility	Any Height
Motors	4 DC Motors
Pendant	Ergonomically designed, rubberized, sealed and water-resistant lockable pendant
Rails	4 rail options available
Bed Ends	1 style and 1 finish
Compliance	UL 606 • IEC 60601-2-52
Assembly	Tool-less assembly
Locking Mechanism	Locking Casters

NOTE: Whenever bed rails are used as part of any hospital or home care bed system, the potential for bed rail entrapment exists. To reduce this risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment please visit the US FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment."

CARROLL CS SERIES BEDS

An industry standard that is anything but.

Since its introduction, the Carroll CS Series line of beds has set a new standard in performance for Post-Acute Care beds. From hallmark design elegance to clinical and facility-friendly features, these beds offer benefits for caregivers and residents alike.

- Slat deck design allows for simple disinfection
- Pure vertical rise
- Tool-free assembly
- Wide variety of bed end styles and finishes
- Manufactured in the USA




CS3™ BED



CS7™ BED

CS SERIES BED RAILS



1. Carroll Assist Rail® Side Rails

- Half rail that mounts to the midsection of the bed
- No Guard 2-way positioning option available (Assist and Transfer only)
- Available for all 36", 39" and 42" CS Series

3. Carroll Assist Bar

- Designed for assisting residents in and out of bed
- Easy pull knob release, quickly folds out of the way with two way positioning
- Available for all 36", 39" and 42" CS Series beds

2. Invacare® ThinkSoft® Positioning Device

- Arthritic-friendly, blow molded device provides positioning support at head end
- Easy for residents to grasp and hold for positioning
- Available for all 36", 39" and 42" CS Series beds

4. Partial Rails

- Traditional rail that extends 3/4 of the length of the bed
- Leaves open area at foot end for access to care or transfers
- Available for Carroll CS Series™ CS3™ Bed, and CS7™ (36" width only)

CARROLL CS SERIES		
	CS3	CS7
Dimensions	36", 39", 42" (W) X 76"- 84" (L)	
Deck-to-Floor Height	8.5"- 27"	7"- 30"
Weight Capacity	450 lbs.	500 lbs.
Warranty	Electronics - 3 Years • Mechanical - 5 Years • Frame - 10 Years	
Deck	Slat Deck	
Mobility	Low Height	Any Height
Motors	4 DC Motors	
Pendant	Ergonomically designed, rubberized, sealed and water-resistant	
Rails	4 rail options available	
Bed Ends	8 styles and 8 finishes available in laminate	
Compliance	CAN/CSA C22.2 No. 601.1 • UL 60601-1 • IEC 60601-2-38	
Assembly	Tool-less assembly	
Locking Mechanism	Locking Casters	One-touch floor lock

NOTE: Whenever bed rails are used as part of any hospital or home care bed system, the potential for bed rail entrapment exists. To reduce this risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment please visit the US FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment."

INVACARE DLX SERIES BEDS

The durable workhorse of Post-Acute Care.

The workhorse of Post-Acute Care beds comes by its durable reputation honestly through years of reliable service and dependable performance.

- Durable grid deck
- Built in hand crank in case of emergencies
- Mobility at low height
- Wide variety of bed end styles and finishes
- Manufactured in the USA



820 DLX BED



SC 900 DLX BED

6-Function Hand Control Pendant

The electric Hand Control features bright, universal symbols and allows residents and staff to easily operate the comfort controls of the bed.



Adjustable Mattress Deck

With adjustable bed-end brackets and corner mattress keepers, the SC900DLX beds are expandable from 76" to 80".



Built-in hand crank


Bed can be raised or lowered in the event of a power outage.

Accessories

- Vinyl bumpers
- Trapeze/mounting frame
- Pan Deck



INVACARE DLX SERIES BED RAILS



1. Carroll Assist Rail® Side Rails

- Half rail that mounts to the mid section of the bed
- Easy pull knob for 3-way positioning: Guard, Assist and Transfer
- No Guard 2-way positioning option available (Assist and Transfer only)
- Low profile for resident comfort and ease of use

2. Invacare® ThinkSoft® Positioning Device

- Arthritic-friendly, blow molded device provides positioning support at head end
- Easy for residents to grasp and hold for positioning
- Reliable, pinch-free latching mechanism

3. Carroll Assist Bar

- Designed for assisting residents in and out of bed
- Easy pull knob release, quickly folds out of the way with two way positioning: Storage and Assist

4. Partial Rails

- Traditional rail that extends 3/4 of the length of the bed
- Leaves open area at foot end for access to care or transfers
- Quick release for easy transfers and resident care

5. Full, Half and Quarter Length Rails

- Provide safety and comfort for the resident
- One button release for caregiver access or transfers
- Half length rails can be used on both head and foot sections as a split rail design

DLX SERIES		
	DLX 820 Bed	SC900DLX
Dimensions	36"(W) X 76" - 80" (L)	
Deck-to-Floor Height	13.75" - 27"	9" - 21"
Weight Capacity	450 lbs.	
Warranty	Electronics - 1 Year • Mechanical - 2 Years • Frame - 10 Years	
Deck	Grid or Pan	
Mobility	Mobile at low height	
Motors	3 DC Motors	
Pendant	Ergonomically designed, rubberized, water-resistant	
Rails	7 rail options available	
Bed Ends	6 styles and 8 finishes available laminate	
Compliance	CAN/CSA C22.2 No. 601-1 • UL 60601-1 • IEC60601-2-38	
Assembly	N/A	
Locking Mechanism	Locking Casters	

NOTE: Whenever bed rails are used as part of any hospital or home care bed system, the potential for bed rail entrapment exists. To reduce this risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment please visit the US FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment."

INVACARE BEDS FOR POST-ACUTE CARE

IH820DLX-KDLSO-QS



820DLX Bed

Item Number	Description
IH820-3MDLX-QS	820DLX Full Electric Bed, Grid Deck - Quick Ship
IH820DLX-KDLSO-QS	820DLX Full Electric Bed, Assist Rail, Kendal Bed Ends in Solar Oak
IH820DLX-KDLAW-QS	820DLX Full Electric Bed, Assist Rail, Kendal Bed Ends in African Walnut
IH820DLX-KDLBC-QS	820DLX Full Electric Bed, Assist Rail, Kendal Bed Ends in Biltmore Cherry
IH820DLX-AMSSO-QS	820DLX Full Electric Bed, Assist Rail, Amherst Bed Ends in Solar Oak
IH820DLX-AMSAW-QS	820DLX Full Electric Bed, Assist Rail, Amherst Bed Ends in African Walnut
IH820DLX-AMSBC-QS	820DLX Full Electric Bed, Assist Rail, Amherst Bed Ends in Biltmore Cherry



SC900DLX Bed

Item Number	Description
SC900-80WDLX-QSP	SC900DLX Full Electric Bed, Grid Deck
SC900-IHKDLSO-QS	SC900DLX Full Electric Bed, Grid Deck, Assist Rail, Kendal Bed Ends in Solar Oak
SC900-IHKDLAW-QS	SC900DLX Full Electric Bed, Grid Deck, Assist Rail, Kendal Bed Ends in African Walnut
SC900-IHKDLBC-QS	SC900DLX Full Electric Bed, Grid Deck, Assist Rail, Kendal Bed Ends in Biltmore Cherry
SC900-IHAMSSO-QS	SC900DLX Full Electric Bed, Grid Deck, Assist Rail, Amherst Bed Ends in Solar Oak
SC900-IHAMSAW-QS	SC900DLX Full Electric Bed, Grid Deck, Assist Rail, Amherst Bed Ends in African Walnut
SC900-IHAMSBC-QS	SC900DLX Full Electric Bed, Grid Deck, Assist Rail, Amherst Bed Ends in Biltmore Cherry



ETUDE HC Bed

Item Number	Description
ETUDE-HC	ETUDE HC Homecare Bed Package

CS3 Beds

NEW
NEW



Item Number	Description
IHCS38-QSP	CS3 bed only
IHCS3W39-QSP	CS3 39" wide, bed only
IHCS3W42-QSP	CS3 42" wide, bed only
IHCS3AMAC-QSP	Package: CS3 Bed, Amherst Bed Ends in Amber Cherry
IHCS3AMBC-QSP	Package: CS3 Bed, Amherst Bed Ends in Biltmore Cherry
IHCS3AMS0-QSP	Package: CS3 Bed, Amherst Bed Ends in Solar Oak
IHCS3DRAC-QSP	Package: CS3 Bed, Drake Bed Ends in Amber Cherry
IHCS3DRBC-QSP	Package: CS3 Bed, Drake Bed Ends in Biltmore Cherry
IHCS3DRSO-QSP	Package: CS3 Bed, Drake Bed Ends in Solar Oak
IHCS3AMACBAAS-QSP	Package: CS3 Bed, Amherst Bed Ends in Amber Cherry, Assist Bar
IHCS3AMBCBAAS-QSP	Package: CS3 Bed, Amherst Bed Ends in Biltmore Cherry, Assist Bar
IHCS3AMSOBAAS-QSP	Package: CS3 Bed, Amherst Bed Ends in Solar Oak, Assist Bar
IHCS3AMWCBAAAS-QSP	Package: CS3 Bed, Amherst Bed Ends in Williamsburg Cherry, Assist Bar
IHCS3DRACBAAS-QSP	Package: CS3 Bed, Drake Bed Ends in Amber Cherry, Assist Bar
IHCS3DRBCBAAS-QSP	Package: CS3 Bed, Drake Bed Ends in Biltmore Cherry, Assist Bar
IHCS3DRSOBAAS-QSP	Package: CS3 Bed, Drake Bed Ends in Solar Oak, Assist Bar
IHCS3AMACRLAS-QSP	Package: CS3 Bed, Amherst Bed Ends in Amber Cherry, Assist Rail
IHCS3AMBCRLAS-QSP	Package: CS3 Bed, Amherst Bed Ends in Biltmore Cherry, Assist Rail
IHCS3AMSORLAS-QSP	Package: CS3 Bed, Amherst Bed Ends in Solar Oak, Assist Rail
IHCS3AMWCRLAS-QSP	Package: CS3 Bed, Amherst Bed Ends in Williamsburg Cherry, Assist Rail

CS3 Beds continued next page

INVACARE BEDS FOR POST-ACUTE CARE



CS3 Beds

Item Number	Description
IHCS3DRBCRLAS-QSP	Package: CS3 Bed, Drake Bed Ends in Biltmore Cherry, Assist Rail
IHCS3DRSORLAS-QSP	Package: CS3 Bed, Drake Bed Ends in Solar Oak, Assist Rail
IHCS3AMACTSPD-QSP	Package: CS3 Bed, Amherst Bed Ends in Amber Cherry, Thinksoft Positioning Device
IHCS3AMBCTSPD-QSP	Package: CS3 Bed, Amherst Bed Ends in Biltmore Cherry, Thinksoft Positioning Device
IHCS3AMSOTSPD-QSP	Package: CS3 Bed, Amherst Bed Ends in Solar Oak, Thinksoft Positioning Device
IHCS3AMWCTSPD-QSP	Package: CS3 Bed, Amherst Bed Ends in Williamsburg Cherry, Thinksoft Positioning Device
IHCS3DRACTSPD-QSP	Package: CS3 Bed, Drake Bed Ends in Amber Cherry, Thinksoft Positioning Device
IHCS3DRBCTSPD-QSP	Package: CS3 Bed, Drake Bed Ends in Biltmore Cherry, Thinksoft Positioning Device
IHCS3DRSOTSPD-QSP	Package: CS3 Bed, Drake Bed Ends in Solar Oak, Thinksoft Positioning Device

CS7 Beds

NEW
NEW
NEW



Item Number	Description
HCS78-QSP	CS7 36", bed only
IHCS7W39-QSP	CS7 39" wide, bed only
IHCS7W42-QSP	CS7 42" wide, bed only
IHCS7ACP-BB	CS7 36" wide, ACP, Battery backup bed only
IHCS7ACP-C	CS7 Bed, Attendant Control Panel setup
IHCS7AMAC-QSP	Package: CS7 Bed, Amherst Bed Ends in Amber Cherry
IHCS7AMBC-QSP	Package: CS7 Bed, Amherst Bed Ends in Biltmore Cherry
IHCS7AMS0-QSP	Package: CS7 Bed, Amherst Bed Ends in Solar Oak
IHCS7DRAC-QSP	Package: CS7 Bed, Drake Bed Ends in Amber Cherry
IHCS7DRBC-QSP	Package: CS7 Bed, Drake Bed Ends in Biltmore Cherry
IHCS7DRSO-QSP	Package: CS7 Bed, Drake Bed Ends in Solar Oak
IHCS7AMACBAAS-QSP	Package: CS7 Bed, Amherst Bed Ends in Amber Cherry, Assist Bar
IHCS7AMAWBAAS-QSP	Package: CS7 Bed, Amherst Bed Ends in African Walnut, Assist Bar
IHCS7AMBCBAAS-QSP	Package: CS7 Bed, Amherst Bed Ends in Biltmore Cherry, Assist Bar
IHCS7AMS0BAAS-QSP	Package: CS7 Bed, Amherst Bed Ends in Solar Oak, Assist Bar
IHCS7AMWCBAAAS-QSP	Package: CS7 Bed, Amherst Bed Ends in Williamsburg Cherry, Assist Bar
IHCS7DRACBAAS-QSP	Package: CS7 Bed, Drake Bed Ends in Amber Cherry, Assist Bar
IHCS7DRBCBAAS-QSP	Package: CS7 Bed, Drake Bed Ends in Biltmore Cherry, Assist Bar
IHCS7DRSOBAAS-QSP	Package: CS7 Bed, Drake Bed Ends in Solar Oak, Assist Bar

CS7 Beds continued next page

INVACARE BEDS FOR POST-ACUTE CARE



CS7 Beds

Item Number	Description
IHCS7AMACRLAS-QSP	Package: CS7 Bed, Amherst Bed Ends in Amber Cherry, Assist Rail
IHCS7AMAWRLAS-QSP	Package: CS7 Bed, Amherst Bed Ends in African Walnut, Assist Rail
IHCS7AMBCRLAS-QSP	Package: CS7 Bed, Amherst Bed Ends in Biltmore Cherry, Assist Rail
IHCS7AMSORLAS-QSP	Package: CS7 Bed, Amherst Bed Ends in Solar Oak, Assist Rail
IHCS7AMWCRLAS-QSP	Package: CS7 Bed, Amherst Bed Ends in Williamsburg Cherry, Assist Rail
IHCS7DRBCRLAS-QSP	Package: CS7 Bed, Drake Bed Ends in Biltmore Cherry, Assist Rail
IHCS7DRSORLAS-QSP	Package: CS7 Bed, Drake Bed Ends in Solar Oak, Assist Rail
IHCS7AMACTSPD-QSP	Package: CS7 Bed, Amherst Bed Ends in Amber Cherry, Thinksoft Positioning Device
IHCS7AMAWTSPD-QSP	Package: CS7 Bed, Amherst Bed Ends in African Walnut, Thinksoft Positioning Device
IHCS7AMBCTSPD-QSP	Package: CS7 Bed, Amherst Bed Ends in Biltmore Cherry, Thinksoft Positioning Device
IHCS7AMSOTSPD-QSP	Package: CS7 Bed, Amherst Bed Ends in Solar Oak, Thinksoft Positioning Device
IHCS7AMWCTSPD-QSP	Package: CS7 Bed, Amherst Bed Ends in Williamsburg Cherry, Thinksoft Positioning Device
IHCS7DRACTSPD-QSP	Package: CS7 Bed, Drake Bed Ends in Amber Cherry, Thinksoft Positioning Device
IHCS7DRBCTSPD-QSP	Package: CS7 Bed, Drake Bed Ends in Biltmore Cherry, Thinksoft Positioning Device
IHCS7DRSOTSPD-QSP	Package: CS7 Bed, Drake Bed Ends in Solar Oak, Thinksoft Positioning Device
IHCS7ACPAMBCBAS-QS	Package: CS7 Bed, Attendant Control Panel, Amherst Bed Ends in Biltmore Cherry, Assist Bar
IHCS7ACPAMSOBAS-QS	Package: CS7 Bed, Attendant Control Panel, Amherst Bed Ends in Solar Oak, Assist Bar
IHCS7ACPAMBCLAS-QS	Package: CS7 Bed, Attendant Control Panel, Amherst Bed Ends in Biltmore Cherry, Assist Rail
IHCS7ACPAMSOLAS-QS	Package: CS7 Bed, Attendant Control Panel, Amherst Bed Ends in Solar Oak, Assist Rail
IHCS7ACPAMBCSPD-QS	Package: CS7 Bed, Attendant Control Panel, Amherst Bed Ends in Biltmore Cherry, Thinksoft Positioning Device
IHCS7ACPAMSOSPD-QS	Package: CS7 Bed, Attendant Control Panel, Amherst Bed Ends in Solar Oak, Thinksoft Positioning Device

INVACARE RAILS AND BED ENDS

CS Series Rails



Item Number	Description	Compatible With
IHCSRLAS-QSP	Assist Rail, 36" Wide Deck	CS3, CS5, CS7 36" Wide
IHCSRLAS-39	Assist Rail, 39" Wide Deck	CS3, CS7 39" Wide
IHCSRLAS-42	Assist Rail, 42" Wide Deck	CS3, CS5, CS7 42" Wide
IHCSBAAS-QSP	Assist Bar, 36" Wide Deck	CS3, CS5, CS7 36" Wide
IHCSBAAS-39	Assist Bar, 39" Wide Deck	CS3, CS7 39" Wide
IHCSBAAS-42	Assist Bar, 42" Wide Deck	CS3, CS5, CS7 42" Wide
IHCSTSPD-QSP	ThinkSoft Positioning Device	CS3, CS5, CS7, CS9 All Sizes
IHCSRLPA	Partial Rail, 36" Wide Deck	CS3, CS5, CS7 36" Wide

DLX Series Rails



Item Number	Description	Compatible With
IH6801A	Full Length Rail	820DLX, SC900DLX
IH6803A	1/4 Length Positioning Rail	820DLX, SC900DLX
IH6808A	3/4 Length Rail	820DLX, SC900DLX
IH6815A	1/2 Rail Head Rail	820DLX, SC900DLX
IH6816A	1/2 Rail Foot Rail	820DLX, SC900DLX
IHRAILAB-DLX	Assist Bar	820DLX, SC900DLX
IHRAILAE-DLX	Assist Rail	820DLX, SC900DLX
IHRAILTS-DLX	ThinkSoft Positioning Device	820DLX, SC900DLX

ETUDE Rails and Accessories



Item Number	Description	Compatible With
EVSR-1823	ETUDE HC Verso Side Rails (PAIR)	ETUDE-HC
ESR-2477	ETUDE HC Side Rail 16" x 12"	ETUDE-HC
ESR-2478	ETUDE HC Side Rail 16" x 16"	ETUDE-HC
ESSHL-2479	ETUDE HC Swivel Support Left	ETUDE-HC
ESSHR-2480	ETUDE HC Swivel Support Right	ETUDE-HC
ELSP-2093	ETUDE HC Lifting Support Pole Kit	ETUDE-HC
ESP-2716	ETUDE HC Wood Side Panels	ETUDE-HC
ETA-3779	ETUDE HC Transport Adaptor Kit	ETUDE-HC
EK-3993	ETUDE HC Deck Extension Kit 4"	ETUDE-HC
EK-3994	ETUDE HC Deck Extension Kit 8"	ETUDE-HC
EBBK-4186	ETUDE HC Battery Back Up Kit	ETUDE-HC

INVACARE RAILS AND BED ENDS



CS Series Bed Ends*

	Item Number	Description	Bed Frame Compatible With
	IHBMTSO-QSP	Belmont Bed Ends in Solar Oak 36" Wide	CS3, CS5, CS7 - Any width frame
	IHCSAMSAC-QSP	Amherst Bed Ends in Amber Cherry 36" Wide	CS3, CS5, CS7 - Any width frame
	IHCSAMSAW-QSP	Amherst Bed Ends in African Walnut 36" Wide	CS3, CS5, CS7 - Any width frame
	IHCSAMSBC-QSP	Amherst Bed Ends in Biltmore Cherry 36" Wide	CS3, CS5, CS7 - Any width frame
	IHCSAMSSO-QSP	Amherst Bed Ends in Solar Oak 36" Wide	CS3, CS5, CS7 - Any width frame
	IHCSAMSWC-QSP	Amherst Bed Ends in Williamsburg Cherry 36" Wide	CS3, CS5, CS7 - Any width frame
	IHDRKAC-QSP	Drake Bed Ends in Amber Cherry 36" Wide	CS3, CS5, CS7 - Any width frame
	IHDRKBC-QSP	Drake Bed Ends in Biltmore Cherry 36" Wide	CS3, CS5, CS7 - Any width frame
	IHDRKSO-QSP	Drake Bed Ends in Solar Oak 36" Wide	CS3, CS5, CS7 - Any width frame
	IHCSAMSAW-ACP	Amherst Bed Ends in African Walnut ACP 36" Wide	IHCS7ACP-C, IHCS7ACP-BB, IHCS990-ACP
	IHCSAMSBC-ACP	Amherst Bed Ends in Biltmore ACP 36" Wide	IHCS7ACP-C, IHCS7ACP-BB, IHCS990-ACP
	IHCSAMSSO-ACP	Amherst Bed Ends in Solar Oak ACP 36" Wide	IHCS7ACP-C, IHCS7ACP-BB, IHCS990-ACP
	IHCSAMSB42-QSP	Amherst Bed Ends in Biltmore Cherry, 42" Wide	CS3, CS5, CS7 - Any width frame
	IHCSAMSS042-QSP	Amherst Bed Ends in Solar Oak, 42" Wide	CS3, CS5, CS7 - Any width frame
	IHCSKDS042-QSP	Kendall Bed Ends in Solar Oak, 42" Wide	CS3, CS5, CS7 - Any width frame
NEW	IHCSAMSCA-42	Amherst Bed Ends in Cafelle, 42" Wide	CS3, CS5, CS7 - Any width frame
NEW	IHCSAMSKM-QSP	Amherst Bed Ends in Kensington Maple, 36" Wide	CS3, CS5, CS7 - Any width frame
NEW	IHCSBMTSO-ACP	Belmont Bed Ends in Solar Oak ACP 36" Wide	IHCS7ACP-C, IHCS7ACP-BB, IHCS990-ACP
NEW	IHCSDRKSO-ACPT	Drake Bed Ends in Solar Oak ACP TREND 36" Wide	IHCS7ACP-C, IHCS7ACP-BB, IHCS990-ACP
NEW	IHCSDRKBC-ACPT	Drake Bed Ends in Biltmore Cherry ACP TREND 36" Wide	IHCS7ACP-C, IHCS7ACP-BB, IHCS990-ACP
NEW	IHCSDRKCA-ACPT	Drake Bed Ends in Cafelle ACP TREND 36" Wide	IHCS7ACP-C, IHCS7ACP-BB, IHCS990-ACP
NEW	IHCS9DRKSO-ACPT	Drake Bed Ends in Solar Oak ACP TREND, Adjustable Width	IHCS990ACP-BB
NEW	IHCS9DRKBC-ACPT	Drake Bed Ends in Biltmore Cherry ACP TREND, Adjustable Width	IHCS990ACP-BB
NEW	IHCS9DRKCA-ACPT	Drake Bed Ends in Cafelle ACP TREND, Adjustable Width	IHCS990ACP-BB
NEW	IHCSDRKCA-QSP	Drake Bed Ends in Cafelle 36" Wide	CS3, CS5, CS7 - Any width frame
NEW	IHCS9DRKSO-ADJ	Drake Bed Ends in Solar Oak Adjustable Width	Replacement bed ends for the CS9 beds with adjustable bed end brackets
NEW	IHCS9DRKBC-ADJ	Drake Bed Ends in Biltmore Cherry Adjustable Width	Replacement bed ends for the CS9 beds with adjustable bed end brackets



DLX Series Bed Ends*

Item Number	Description	Bed Frame Compatible With
IHAMSAW-QSP	Amherst Bed Ends in African Walnut 36" Wide	820DLX, SC900DLX
IHAMSB42-QSP	Amherst Bed Ends in Biltmore Cherry 36" Wide	820DLX, SC900DLX
IHAMSSO-QSP	Amherst Bed Ends in Solar Oak 36" Wide	820DLX, SC900DLX
IHKDLAW-QSP	Kendall Bed Ends in African Walnut 36" Wide	820DLX, SC900DLX
IHKDLBC-QSP	Kendall Bed Ends in Biltmore Cherry 36" Wide	820DLX, SC900DLX
IHKDLSO-QSP	Kendall Bed Ends in Solar Oak 36" Wide	820DLX, SC900DLX



Etude Bed Ends

Item Number	Description	Bed Frame Compatible With
ETUDEHFB	ETUDE HC Head and Footboard	ETUDE-HC

*Builder Bed Ends are available with a 4 week lead time.

THERAPEUTIC SUPPORT SURFACES

Create an advanced bed system by combining any Invacare bed for Post-Acute Care with a complete range of Invacare Therapeutic Support Surfaces. From the innovative Invacare® Glissando™ Gliding Mattress to the new microAIR® powered mattresses, you can rest easy knowing residents are resting in clinical comfort.



INVACARE® GLISSANDO® GLIDING MATTRESS


Designed to help prevent friction and shear that can contribute to pressure ulcer development, while providing comfort for the resident.

Dual Layers of High Density Foam

Two-part construction with a patented "gliding" layer enable top and bottom sections to move independently as mattress is articulated.

Castellations

Break surface tension of foam allowing maximized body contact, optimizing pressure reduction.



INVACARE® MICROAIR® THERAPEUTIC SUPPORT SURFACES

Designed to be clinically effective yet durable, with easily removable cells for maintenance and a unique detachable bladder sleeve design for easy cleaning.

Designed for Clinical Efficacy

- HeelSense™ Technology on microAIR 600, microAIR 800, and microAIR 1000 Therapeutic Support Surfaces
- Nylon non-slip base within covers
- Soft nylon urethane backed cells
- Modular air cell design for easy replacement and cleaning
- Unique bladder sleeve design provides stable cell support and is detachable from the mattress base for ease of cleaning



Yes, you can:

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324
©2020

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols™ and®.
All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. Etude is a registered trademark in the EU, United Kingdom, Switzerland, Norway, Australia and New Zealand. Specifications are subject to change without notification
Form No. 20-002 Rev. 10/21
210825

