


INVACARE®  
**BED RAILS**  
FOR INVACARE® DLX AND  
CARROLL CS SERIES™ BEDS



From a company that strives to ensure individual safety, you will have the peace of mind knowing that **Carroll CS Series™ and Invacare® DLX Series Rails** are classified to the international standard, IEC 60601-2-38 when used as part of an Invacare bed system.\*

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit [www.invacare.com/bedsafety](http://www.invacare.com/bedsafety) for additional safety and entrapment information.



### 1. Assist Rail

Half rail that mounts to the midsection of the bed  
 Easy pull knob for 3-way positioning: Guard, Assist and Transfer  
 No Guard 2-way positioning option available (Assist and Transfer only)  
 Low profile for resident comfort and ease of use  
 Available for all 36", 39" and 42" CS Series and DLX Series Beds

### 2. Invacare® ThinkSoft® Positioning Device

Arthritic-friendly, blow molded device provides positioning support at head end  
 Easy for residents to grasp and hold for positioning  
 Reliable, pinch-free latching mechanism  
 Available for all 36", 39" and 42" CS Series beds and DLX Series Beds

### 3. Carroll® Assist Bar

Designed for assisting residents in and out of bed  
 Easy pull knob release, quickly folds out of the way with two way positioning: Storage and Assist  
 Available for all 36", 39" and 42" CS Series beds and DLX Series Beds

### 4. Partial Rails

Traditional rail that extends 3/4 of the length of the bed  
 Leaves open area at foot end for access to care or transfers  
 Quick release for easy transfers and resident care  
 Available for Carroll CS Series™ CS3™, CS5™, and CS7™ (36" width only) and DLX series beds

### 5. Full, Half and Quarter Length Rails

Provide safety and comfort for the resident  
 One button release for caregiver access or transfers  
 Half length rails can be used on both head and foot sections as a split rail design  
 Available for DLX Series only



Model #	Rail Dimensions L x W	Limited Warranty	Description
<b>CS Series</b>			
IHCSRLAS	26.5" x 20"	1 year	Assist Rail, 36" Wide Deck
IHCSRLAS-39	26.5" x 20"	1 year	Assist Rail, 39" Wide Deck
IHCSRLAS-42	26.5" x 20"	1 year	Assist Rail, 42" Wide Deck
IHCSRLAS-NG	26.5" x 20"	1 year	Assist Rail, No Guard, 36" Wide Deck
IHCSRLAS-39NG	26.5" x 20"	1 year	Assist Rail, No Guard, 39" Wide Deck
IHCSRLAS-42NG	26.5" x 20"	1 year	Assist Rail, No Guard, 42" Wide Deck
IHCSBAAS	5" x 23.5"	1 year	Assist Bar, 36" Wide Deck
IHCSBAAS-39	5" x 23.5"	1 year	Assist Bar, 39" Wide Deck
IHCSBAAS-42	5" x 23.5"	1 year	Assist Bar, 42" Wide Deck
IHCSTSPD	24" x 15"	1 year	ThinkSoft Positioning Device
IHCSRLPA	68" x 13"	1 year	Partial Rail, 36" Wide Deck
<b>DLX Series</b>			
IH6801A-DLX	67" x 13"	1 year	Full Length Rail
IH6808A-DLX	55" x 13"	1 year	Partial Rail
IH6803A-DLX	22" x 13"	1 year	Quarter Length Positioning Rail
IH6815A-DLX	26.5" x 13"	1 year	Half Length Head Rail
IH6816A-DLX	26.5" x 13"	1 year	Half Length Foot Rail
IHRAILAB-DLX	5" x 23.5"	1 year	Assist Bar
IHRAILAE-DLX	26.5" x 20"	1 year	Assist Rail
IHRAILTS-DLX	24" x 15"	1 year	ThinkSoft Positioning Device

**Invacare Corporation**  
[www.invacare-cc.com](http://www.invacare-cc.com)

**USA**  
 One Invacare Way  
 Elyria, Ohio  
 44035-4190  
 (800) 668-2337

**Canada**  
 570 Matheson Blvd. E.,  
 Unit 8  
 Mississauga, Ontario  
 L4Z 4G4 Canada  
 (800) 668-5324

©2017  
 Invacare Corporation.  
 All rights reserved.  
 Trademarks are  
 identified by the  
 symbols™ and®.  
 All trademarks are  
 owned by or licensed  
 to Invacare Corporation  
 unless otherwise  
 noted. Specifications  
 are subject to change  
 without notification.  
 Form No. 12-619  
 Rev. 10/17  
 171052



**Yes, you can.®**

\*Testing conducted to select portions of UL/IEC Standards.