

CS3 bed is featured with Assist Bar, Drake Bed Ends and Invacare* Solace* Mattress

CARROLL CS SERIES™ CS3™ BED



The **Carroll CS Series™ CS3™ Bed** was developed with safety in mind for those who require a bed that is stable throughout the travel range and mobile in the lowest position only.

Invacare® beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the recommended FDA Guidelines for bed rail entrapment. Visit **www.invacare.com/bedsafety** for additional safety and entrapment information.

Features

Mobility at Low Height – The Maxi-Lock system secures the bed to the floor in all elevated positions, providing reassurance that the bed is locked and secured for maximum resident safety.

Exclusive Slat Deck Design - Provides a comfortable sleep surface, allows for improved infection control and helps reduce wear and tear on the mattress.

Ergonomic Hand Pendant – Water-resistant¹, sealed hand pendant designed to fit comfortably in the resident's hand.

Dual Pendant Connections – Pendant connects to either side of the bed, providing customized convenience for the resident.

Tool-less Assembly - The CS3 bed features tool-less assembly, minimizing set-up time.

Customizable Options -

- 39" or 42" Wide Deck
- 4" Bed Extender (84" Length)
- Attendant Control Panel
- Portable Battery Backup
- Hi-Lo Lockout



2. Maxi-Lock System 4. S

Model#

IHCS3

Height Range

8 5" - 27" H

Sleen

Surface

76" - 80" I

36", 39", 42" W



450 lb







Compliance

†Classified to:
 UL 60601-1
IEC 60601-2-38
CAN/CSA C22.2

No. 601.1

Regulation

Warranty

10 Years

Frame
5 Years

Limited

5 Years Mechanical 3 Years Electronics

NOTE: Whenever bed rails are used as part of any hospital or homecare bed system, the potential for bed rail entrapment exists. To reduce this risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment please visit the US FDA's website at http://www.fda.gov and search for "Bed Rail Entrapment."

1. Pendant is water-resistant to IPX6 per IEC 60529.

†Testing conducted to select portions of UL/IEC Standards.



Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

©2020

Invacare Corporation.

All rights reserved.

Trademarks are identified by the symbols ™ and

®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject to change without notification.

Form No. 20-030 200202

^{*} Maximum safe working load.