




PUMP REPLACEMENT Diaphragms

#letloveflo
f i y
@evenflofeeding

Missing or worn out diaphragms can reduce your pump's suction and overall performance. These extra diaphragms work with your Evenflo Advanced Electric Pumps.

TO USE: Sanitize all parts prior to use. Install diaphragm (A) and diaphragm cap (B) to flange assembly (C). Reference your Evenflo pump manual or go to www.evenflofeeding.com for more detailed instructions.

WARNING: FOR ADULT USE ONLY – Keep out of reach of children. This product contains small parts which could present a choking hazard.



Les membranes manquant ou usé peut réduire l'aspiration et les performances de votre tire-lait. Ces membranes supplémentaires fonctionnent avec les tire-lait à la fine pointe de la technologie Evenflo.

UTILISATION : Nettoyer toutes les pièces avant l'utilisation. Monter le clapet (A) et le bouchon du clapet (B) sur l'ensemble de l'entonnoir (C). Se reporter au manuel du tire-lait Evenflo ou se rendre sur www.evenflofeeding.com pour davantage d'instructions détaillées.

AVERTISSEMENT : POUR ADULTES UNIQUEMENT – Garder hors de la portée des enfants. Ce produit contient des petites pièces pouvant présenter un risque d'étouffement.

Parent's Ally™

USA / Canada 1-855-334-2229

www.evenflofeeding.com

EVENFLO FEEDING, INC.
WEST CHESTER, OH 45069

©2019 EVENFLO FEEDING, INC.

Made in China/Fabriqué en Chine

2951-558 30MAY2019 5146112