



Every Drop Counts®



neonatal
solutions

Medela's human milk management solutions

Medela's human milk management solutions give NICU professionals and other healthcare providers easy, convenient solutions for handling, storing, and transporting breast milk for premature babies. Whether preparing bulk breast milk feeds, initiating and collecting colostrum or breast milk from mothers, or transferring contents to ENFit® or oral syringes, our solutions are designed with healthcare professionals and mothers in mind.



Colostrum Protector

101046344

Features & Benefits:

- Minimizes loss of colostrum by catching it prior to reaching the pump parts.
- Easy to use and seamlessly connects to the Symphony® Breast Pump, without impacting pump performance.
- More effective than traditional methods¹
- 95% of nurses indicated it should define the standard of practice¹
- 89% of nurses recommend it¹



Breast Milk Transfer Lid

90003S-100 or ENFTL

250 mL Breast Milk Bulk Preparation System

90001S-100

Features & Benefits:

- Helps prevent spills or loss of breast milk.
- Easier ENFit or oral syringe milk transfers.
- Sterile and ready-to-use.*
- Non-contact breast milk transfers from both Medela and most commonly used containers.

For Enteral Feeding Only. Not for I.V. use.

*in undamaged, unopened package



35 mL Collection & Storage Container

87110S-50

- Panel for writing name, date and time.
- Storage capacity up to 35 mL.

**35mL container has mL measurements only.*



80 mL Breast Milk Collection Container

101036902

- The curved bottom helps minimize loss of breast milk.
- Storage capacity up to 80 mL.



250 mL Breast Milk Container

87150S-100

- Larger capacity for bulk feed preparation.
- Patterned volume marking for improved visibility.

Collection Container Features & Benefits:

- Screw-on lid provides leak-proof seal with Medela containers.
- Easy-to-read ounce and milliliter measurement.*
- Wide base helps eliminate tip-overs.

Caution: Plastic bottles and component parts become brittle when frozen and may break when dropped. Take appropriate care in handling bottles and components. Do not use the breast milk if bottles or components become damaged.

Product Specifications:

Product Name	Legacy Item #	SAP Item #	Qty/Case	Shipping Carton Weight (lbs)	Shipping Carton Dimensions
Colostrum Protector 2 Pack (Sterile)	N/A	101046344	40 Packs of 2	2.40	16 x 10.875 x 9.25
250 mL Breast Milk Bulk Preparation System (Sterile)	90001S-100	101028973	100 per case	9.66	21.91 x 14.51 x 18.91
Breast Milk Transfer Lid, oral (Sterile)	90003S-100	101028971	100 per case	2.60	8.61 x 15.81 x 11.41
Breast Milk Transfer Lid, ENFit® (Sterile)	ENFTL	101033082	100 per case	2.87	16.10 x 11.40 x 8.70
35 mL Container Bulk (Sterile)	87110S-50	101028972	50 Packs of 2	3.56	9.01 x 16.11 x 10.31
80 mL Breast Milk Collection Container Bulk (Ready-to-Use)	N/A	101036902	100 per case	5.35	17.13 x 5.25 x 19.00
250 mL Breast Milk Container Bulk (Sterile)	87150S-100	101028974	100 per case	9.49	21.51 x 18.71 x 14.41

To order or for more information, contact your Medela Sales Consultant.

1. Medela 2022 Colostrum Protector Evaluation Survey

2. Commercial Discharge Packs and Breast-Feeding Counseling: Effects on Infant-Feeding Practices in a Randomized Trial – Deborah A. Frank, Stephen J. Wirtz, James R. Sorenson and Timothy Heeren Pediatrics 1987;80:845-85

Medela wordmark and logo, Every Drop Counts are registered in the U.S. Patent and Trademark Office. Quick Clean and Micro-Steam are trademarks of Medela. ENFit is a federally registered trademark of GEDSA and is used with permission.