

Enteral Feeding: Critical Information for Informed Choices

Color is a Preference, Not a Mandate

There is no required color standard for enteral feeding products established by the U.S. Food and Drug Administration (FDA), The International Organization for Standardization (ISO), The Joint Commission (TJC), The Institute for Safe Medication Practices (ISMP), or any other regulatory or accrediting body.

ENFit® or TwistLok®: The Choice is Yours

- ENFit is a global standard developed by the Global Enteral Device Supplier Association (GEDSA), GEDSA consists of medical device manufacturers and distributors. **ENFit adoption is optional.**
- TwistLok connectors remain a trusted choice in many hospitals—valued for their familiar handling, precise performance, and easy maintenance.
- TwistLok is compatible with oral syringes, enhancing flexibility for both oral and enteral administration - enabling safe and accurate oral delivery.



ENFit or TwistLok: Your connection, Your choice

Medela offers enteral feeding products with both connectors for safe, compatible care.

Contact your NICU & Maternity Sales Consultant to learn more.

TwistLok[®] Connectors

Designed for neonatal precision, safety and usability

Safety: *Prevents Misconnections* -

Like ENFit[®], TwistLok is designed to be incompatible with IV systems.

Precision: *Designed for Small Volume Dosing* -

TwistLok is designed to reduce over- and under-dosing, a concern with ENFit connectors.

Cleanliness: *Smaller/No Moat Area* -

Reducing residual milk or medication build-up and complex cleaning protocols.

Efficiency: *Simpler Transition* -

TwistLok is easy to adopt with minimal training.



ENFit or TwistLok: Your connection, Your choice

Medela offers enteral feeding products with both connectors for safe, compatible care.

Contact your NICU & Maternity Sales Consultant to learn more.