

Value Care



Semi-electric Bed

Rails should be installed per the Instructions for Use (IFU) provided with/on the Value Care Bed.



The *Invacare Value Care Semi-Electric Bed* utilizes a single motor design making the bed lightweight and easy to assemble. The advanced pendant is designed to allow the head and knee sections to move either simultaneously or independently at the push of a button.

Features

- ▶ Manual crank provides convenient protection in case of power outages
- ▶ Pendant is designed to allow the head and knee section to operate simultaneously or independently
- ▶ 450 lb.*: Safe Working load
- ▶ Metal frame construction
- ▶ Hand pendant with easy to read buttons
- ▶ Electric head and knee elevation
- ▶ Manual crank height elevation
- ▶ **HCPCS Code: E0295****



▲ Use the detachable hand crank to raise and lower the height of the bed. Perfect for power outages. When not in use, store the hand crank on the bracket at the end of the bed.



▲ The dual motor system is easy to assemble/disassemble and comes with battery backup in case of power outage (battery sold separately).



▲ The Carbon steel material and the double support structure gives the Invacare Value Care Semi-Electric Bed high loading capacity and stable deck adjustment.



▲ Mattress retainer at the foot end helps stop the mattress from sliding during articulation.

Technical Specifications ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Dimensions ▼	Sleep Surface ▼	Weight Capacity ▼	Shipping Weight ▼	Carton Dimensions ▼	Limited Warranty ▼
5310VC	87.5"L 15"-22.5"H 36"W	80"L 36"W	Safe Working Load: 450 lb.* Max. Patient Weight: 350 lb.	N/A	N/A	Frame - 3 Years Motor - 1 Year All Other Parts - 3 Years
5310VCDECK				114 lb.	45" L 9.5" W 38" D	Frame - 3 Years Motor - 1 Year All Other Parts - 3 Years
5310VCHFB Footboard				55 lb.	39" L 6" W 27" D	Frame - 3 Years Motor - 1 Year All Other Parts - 3 Years

*Includes patient, mattress and all accessories.

NOTE: Whenever bed rails are used as part of any hospital or homecare bed system, the potential for bed rail entrapment exists. To reduce this risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment please visit the US FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment."

**All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

Tel: (800) 333-6900
www.invacare.com

Invacare Corporation
One Invacare Way
Elyria, Ohio
44035-4190