

Redefining standards in care



A new adjustable electric bed

Yes, you can:



Promoting comfort and rest

Introducing the new generation of CS Series beds from Invacare. Our new CS600™ adjustable electric bed promotes comfort and rest for patients at home or in long-term care environments. With a focus on safety and on delivering our highest quality standards to ensure all of our beds support patients as they rest.



Introducing the new CS600™ adjustable electric bed, the next generation of our popular CS Series of beds.

Built from a design foundation of exceptional comfort, safety and versatility, the CS600™ design is focused on improving standards in patient care.

Safety first. No compromise.

When it comes to patient and caregiver safety there is no compromise. The CS600[™] bed meets the highest international safety and performance standard.*

*Medical Bed Equipment Standard for safety, quality and performance - IEC 60601-2-52.





Safe

Reassurance in quality is at the core of what we do. The CS600™ is tested to the highest safety and performance standard* for peace of mind.





The versatile width design of the CS600™ and up to 600lb safe working load meets a varied profile of patients and care environments.



upright position helps to improve circulation

and general healing.



Efficient

Tool-free assembly reduces

setup time and and complex

installation. Slat deck design

and wipe down electronics

make cleaning a breeze.

helps to promote improvements in infection control.

Choice of bed ends

Available in 3 finishes

One touch floor lock

Ease and efficiency

Accessories available:



Quarter length and **Assist rail** options





set to 36", 39" or 42"



secures in 3 positions to accommodate patient limitations. 200 lb safe working load

Safe working

load 600 lbs



for easy transportation

Backlit hand control for improved visibility



Battery backup for operation in case of power outage



for improved floor visibility



Height adjustable

The CS600™ lowers to

the possibility of patient falls. For caregiver convenience, the bed also reaches a height of 30".

the floor. A lower bed height may help reduce

Blue: Indicates Item Number Technical data 🔻













Bed Deck Length ▼

Bed Deck Width ▼

Overall Length (incl. head/foot boards) ▼

Maximum Safe Working Load ▼

Product Weight ▼

Minimum/Maximum Bed Deck Height ▼

(CS600-36																				

80"

36'

88'

600 lbs

217 lbs

7" - 30"

CS600-EXP

80"

36", 39", 42"

88"

600 lbs

243 lbs

7" - 30"













Maximum Back Angle ▼

Maximum Knee Angle ▼

Maximum Reverse Trendelenburg v

Length Extension ▼

Width Extension ▼

ACP (Attendant Control Panel) ▼

60°+5°

30°+5°

80" - 84" w/accessory Optional

80" - 84"

w/accessory

Optional

CS600-EXP 60°± 5°

30°± 5°

17°± 5°

36"

w/accessory

Integrated

Optional

Frame

finish

Bed end finishes

CS600-36



Solar Oak

36" CS600BE36S0

42" CS600BE42S0

African

Walnut

CS600BE36AF

CS600BE42AF

Biltmore

Cherry

CS600BE36BC

CS600BE42BC



Light Gray

Bed end style



Amherst

Please note: colors may vary slightly from those displayed above

Quarter length rail options

Standard width, CS600-36

Head section, Qty: 2 CS600HALFRAILH

Foot section, Qty: 2 CS600HALFRAILF

Expandable width, CS600-EXP

Head section, Qty: 2 CS600HALFRAILHEXP

Foot section, Qty: 2 CS600HALFRAILFEXP Assist rail options

Standard width, CS600-36

Assist rail, Qty: 1 CS600ASTRAIL36 Expandable width, CS600-EXP

Assist rail, Qty: 1 CS600ASTRAILEXP

Other

Length extension

Trolley CS600TROLLEY **Battery backup** CS600BATBCK Width extension CS600WIDTHEXP **ACP (Attendant Control Panel)** CS600ACPPAD

accessories

CS600LENGTHEXP

Trapeze CS600TRAPEZE

Underbed light CS600UNDLIGHT

Invacare Corporation Tel: (800) 333-6900 Worldwide Headquarters www.invacare.com

1 Invacare Way Elyria, OH 44035-4190





©2024 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and *. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Please refer to the User Manual for instructions for use. Form: 23-020 Job: 240135 US-04/2024