




DynaGinate™ Calcium Alginate Dressings



ITEM 3010	4" x 4"	10 / bx	12 bxs / cs	NEW SIZE
ITEM 3026	2" x 2"	10 / bx	12 bxs / cs	
ITEM 3027	4¼" x 4¼"	10 / bx	12 bxs / cs	
ITEM 3028	4" x 8"	5 / bx	12 bxs / cs	
ITEM 3029	2 g x 30 cm (Rope)	5 / bx	12 bxs / cs	

Dürma+ DynaGinate Calcium Alginate Dressings form a soft, fibrous gel upon contact with wound exudate, helping to create a moist environment while contouring to the wound bed.

- Calcium alginate dressing transforms into a fibrous gel upon contact with wound exudate
- Soft gel contours to the wound bed and does not adhere to tissue during dressing changes
- Highly absorbent; ideal for heavily exuding wounds
- Maintains a moist wound environment to support the healing process
- Suitable as a primary dressing for partial- and full-thickness wounds at stages 3 or 4



1.888 DYNAREX | DYNAREX.COM |

155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |