Pegasus Metro Feel the Freedom

Quality, comfort and safety

for independence, confidence and style



Yes, you can:



Feel the difference

Our scooters have been designed together with healthcare experts, to meet our highest standards in quality, safety and comfort.

Through continuous dedicated research and our extensive experience in designing for mobility, at Invacare we pride ourselves in scooters that help improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From ergonomic steering controls to the adaptable cushioned seats and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.

Quality



High-quality build with exceptional attention to detail.

Comfort



Adjustable design to help provide lasting comfort.



Feel in control with our Pegasus Metro scooter

Introducing our Pegasus Metro, a high-quality scooter designed to enhance safety, comfort and added driving performance. The Pegasus Metro offers everything you need from a mobility scooter.

Compact and easy to drive, this scooter handles a wide variety of surfaces. Designed with a powerful motor, the Pegasus Metro drives across uneven ground, around corners and over most obstacles without

compromising safety or comfort. Incredibly simple to steer and start, the Pegasus Metro puts you in total control.

up to

12 km/h

max speed



136 kg max user weight



Advanced driving control

- Ergonomic steering for added control
- Designed for all-terrain performance*
- · Anti-splash guards

Designed to help **provide traction**

around corners

performance*

and safe, all-terrain



Safety in mind

- Automatic speed reduction helps reduce speed around corners
- Two-step disengaging lever to help prevent free-wheeling
- Brake lights and indicators
- · Headlights for improved visibility



Built around comfort

- · High quality cushioned seat
- · Adjustable steering column
- · Adjustable seat height
- Flip-up armrests and swivel seat makes getting on and off easy





Swivel seat for easy transfer



High quality lighting




Two-step disengaging lever to help prevent free-wheeling



Add to your scooter

Tailor the Pegasus Metro to your taste and needs. Turn over to explore the optional extras available.



Add to your scooter

We have a wide range of optional extras available, so you can customise your Pegasus Metro scooter. Ask your retailer to discover the full range.



Rear basket



view mirrors



bumper





Crutch holder



Seatbelt






shroud panels

Technical data

For more information about this product, including the user manual, please visit our website www.pro.invacare.ca





Seat depth v

Seat to floor height ▼







Driving unit width ▼

Total length V

470 mm

Min. 440 mm Max. 510 mm

Comfort: 500 mm Deluxe: 490 mm

660 mm

1270 mm



Empty weight ▼

The heaviest part ▼



Max. user weight ▼



Battery capacity v



Turning diameter v



Max. climbable obstacle height 🔻

110 kg

57 kg

136 kg

2 x 12 V/50 Ah (C20) AGM

2600 mm

80 mm

Ground clearance v



Max. safe slope ▼



Motor capacity ▼



Max. driving distance ▼



Max. driving speed ▼

80 mm

10 °

6 km/h: 240 W / 500 W¹ 10 km/h: 240 W / 600 W¹ 12 km/h: 550 W / 1500 W ¹

54 km

6 / 10 / 12 km/h

Colours*







Silver Grey





Upholstery



Black Vinyl

Tel: (800) 668-5324 canadaorders@invacare.com www.pro.invacare.ca



Yes, you can:

Peak performance

^{*} Scooter colours may vary slightly. Please check with your retailer.