

Invacare Leo

Feel the Freedom



Compact, safe with an exceptional drive

Bringing independence to those moments that matter



Yes, you can.

Feel the difference

Our scooters have been designed together with healthcare experts, to meet our highest standards in quality, safety and comfort.

Through continuous dedicated research and our extensive experience in designing for mobility, at Invacare we pride ourselves in scooters that help improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From transportable designs, cushioned seats and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.

Quality



High-quality build with exceptional attention to detail.

Comfort



Adaptive and adjustable design to help provide lasting comfort.



Feel in control with our Invacare Leo scooter

The Invacare Leo is a compact scooter that gives a smooth, safe drive. Simply dismantle it for easy transportation and take it with you for lots of exciting days out and those moments that matter.

Leo is packed with many safety and performance features. Its intelligent lighting package includes indicators that automatically switch off and a rear brake light for safe driving. A visual warning for the battery

level display, alerts with an audible warning (beep) when the battery needs recharging. Incredibly easy to use the Leo puts you in control with confidence.

up to
8km/h
max speed

up to
38km
max distance
on a single
charge

136 kg
max user
weight



Easy to use and transport

- Dismantles into five parts for easy transportation
- Ergonomic steering for added control
- Simple to steer and drive



Safety in mind

- Two-step disengaging lever to help prevent free-wheeling
- Brake lights, headlights and indicators
- Additional optional extras available



Built around comfort

- High quality cushioned seat
- Adjustable seat height
- Flip up armrests and swivel seat makes getting on and off easy



Adjustable steering for maximum comfort



Easy battery charge

with off-board charger kit



Anti-splash guards help protect the scooter from dirt and water



Swivel seat for easy transfer



Easy to transport in your car



Two-step disengaging lever to help prevent free-wheeling



Add to your scooter

Tailor the Invacare Leo to your taste and needs. Turn over to explore the optional extras available.

Add to your scooter

We have a wide range of optional extras available, so you can customize your Invacare Leo scooter. Ask your retailer to discover the full range.



Large mirror



Crutch holder



Seatbelt



Lockable front box

Technical data

For more information about this product, including the user manual, please visit our website www.pro.invacare.ca



Backrest width ▼

470 mm



Seat depth ▼

410 mm



Seat to floor height ▼

460 - 560 mm



Backrest height ▼

475 mm



Driving unit width ▼

590 mm



Total length ▼

1220 mm



Empty weight ▼

59 kg (Excl. batteries)
83 kg (Incl. batteries)



Max. user weight ▼

136 kg



Battery capacity ▼

2 x 36 Ah



Turning diameter ▼

2620 mm



Max. climbable obstacle height ▼

60 mm



Max. safe slope ▼

10°



Motor capacity ▼

240 W



Max. driving distance ▼

38 km



Max. driving speed ▼

8 km/h

Colours*



Onyx
Blue



Ruby
Red

Upholstery



Black Vinyl

* Scooter colours may vary slightly.
Please check with your retailer.

Tel: (800) 668-5324
canadaorders@invacare.com
www.pro.invacare.ca

Invacare Canada LP
84 Citation Drive, Unit 11
Concord, ON, L4K 3C1
Canada



Yes, you can.®