



New Product Announcement

Oxivir® Tb RTU – 5 Gal./18.927 L Pail



In an effort to supply greater amounts of disinfectant to critical areas, Diversey has launched a new Oxivir® Tb product in a 5 gallon / 18.927 L pail.

Technical data	Oxivir® Tb RTU
Diversey SKU	101104054
Format	5 Gallon / 18.927 L Pail
Yield	20 - 32 oz. / 946 mL bottles per pail
Dimensions	12" x 12" x 15.6" (L x W x H)
Weight	45.18 lbs. / 20.493 kg each

Chemical Transfer Options: Due to the unique characteristics of products sold in bulk, operators will need to implement a solution that will allow them to transfer product from a pail to smaller containers/ bottles. While Diversey cannot provide a transfer solution for every situation, there are simple options that will work for most.

Preferred Option: Top Mounted Manual Pump

The preferred option for dispensing chemical from a 5 gallon / 18.927 L pail is a top-mounted mechanical pump that uses a suction tube to pull/pump out the solution. These pumps attach to the 3/4" NPT threaded knock-out in the center of the 70mm cap found on the top of the pail. **For transferring Accelerated Hydrogen Peroxide AHP® products, only a corrosion resistant (non-metal) pump should be used.** Basic tools may be required for installation.

Diversey is able to provide a limited number of [Tolco 160150 - Model TDP-8 pumps](#) for purchase though [Plan B](#) or they can also be sourced through the manufacturer or other suppliers. Diversey Representatives are not able to purchase or provide these pumps on behalf of distributors or customers.

It is the responsibility of the operator to ensure the safety of users. Users should wear all proper Personal Protective Equipment (PPE), and follow all guidelines for safe handling of chemicals.

Product Labels: Aftermarket labels should always be used on any container that a chemical is transferred to.

Due to high demand, Diversey may not be able to supply pre-printed aftermarket labels. Distributors and End-Users should [click here to request self-print .pdf aftermarket label files](#). Please note that it is unlawful to re-sell any dispensed chemical.

Questions: If you have questions please contact your Diversey sales manager or call Diversey Customer Technical Support (TCS) at (800) 558-2332.



70mm cap with 3/4" NPT Threaded Knock-Out



Tolco Model TDP-8