



**SparClean** 

**60**

## Coffee & Tea Carafe Destainer

SparClean® Coffee & Tea Carafe Destainer quickly and easily removes coffee, tea, and other tough stains from beverage serving equipment. Designed for use in conjunction with Spartan's 3 sink wash-rinse-sanitize system, Coffee & Tea Carafe Destainer ensures that coffee pots, urns, tea and beverage pitchers, etc. look clean and stain-free.



Providing Quality Maintenance and Chemical Specialty Solutions with Unparalleled Customer Service



Spartan Chemical Company, Inc.  
 1110 Spartan Drive  
 Maumee, OH 43537  
 1-800-537-8990  
 www.spartanchemical.com

**FEATURES AND BENEFITS:**

- Works fast – quickly removes stains from coffee pots, urns, carafes, tea and beverage pitchers, and more
- No waste – Easy-to-use pump dispenses just enough solution for the job
- Delights customers – Leaves coffee pots sparkling clean

**DIRECTIONS FOR USE**

Some individuals may be sensitive to ingredients in this product. Before using, read product label and Safety Data Sheet. If questions remain, consult your employer or a physician.

1. Remove any remaining coffee from pot.
2. Pump 2 squirts of SparClean Coffee & Tea Carafe Destainer directly into the coffee pot.
3. Fill coffee pot to top with warm/hot water.
4. Allow to set for at least 30 seconds.
5. Empty solution from coffee pot and wash in soapy water.
6. Rinse thoroughly.

**NOTE: CONTAINS BLEACH. Contact with fabric, clothing or carpet will cause discoloration.** Prolonged contact may discolor metal trim and laminate surfaces.

**SPECIFICATION DATA:**

**Specific Gravity:** 1.036 g/ml

**Density:** 8.64 lbs./gal.

**Viscosity:** Water thin

**Flash Point (Cleveland Open Cup):** None

**Storage Stability:**

- a. Shelf @ 24°C/75°F – 1 year minimum
- b. Accelerated @ 49°C/120°F – 60 days minimum
- c. Freeze/thaw cycle: Product will withstand a minimum of three freeze/thaw cycle.

**Miscibility:** Completely soluble in all proportions with hot or cold water

**Rinsability:** Rinses freely with hot or cold water and leaves no film.

Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Material Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting [www.spartanchemical.com](http://www.spartanchemical.com).



Distributed by:

PACKAGING	
766003	1-quart (12 per case - 2 pumps included)



**GUARANTEE:** Spartan's modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of the date of manufacture. Use product as directed and read all precautionary statements.

For institutional and industrial use only.

© SCC 07/14 L7660