

HD LONG-TERM CARE BED DB300



CREATING REORDER NUMBERS FOR YOUR BED COMBINATIONS

STEP 1: START WITH THE BED REORDER NUMBER (12013)


STEP 2: ADD THE 2 LETTER CODE FOR THE BOARDS OF YOUR CHOICE




MH - MAHOGANY



CH - AMBER CHERRY



OA - LIGHT OAK



CM - COMPOSITE

STEP 3: ADD THE 2 LETTER CODE FOR THE RAIL/BAR OF YOUR CHOICE



MR - METAL SWING RAIL



CR - COMPOSITE
SWING RAIL



AB - ROTATING
ASSIST BAR



PB - PIVOT ASSIST BAR

EXAMPLE: 12013-MHMR