

flex™

Flexible Modular Seating

FLEXSTEEL®

HOW TO:

STORE YOUR SLEEP KIT

INCLUDED PIECES

- (1) Memory Foam Topper
(in a bamboo zippered cover)
- (1) Fitted Sheet
- (1) Top Sheet
- (1) Pillowcase
- (1) Zippered Storage Bag



COOL COMFORT





BAMBOO COVER



ANTI SLIP
DESIGN





MACHINE
WASHABLE




1.

Remove and fold fitted sheet, top sheet and pillowcase from your Flex™ cushions and set aside. Zip open the storage bag and also set aside.





3.

Fold the opposite end over the previously folded side and slip both elastic bands over the folded memory foam topper.



2.

Slip off memory foam topper's elastic bands from the corners of each cushion. Begin to fold one end in thirds.



4.

Place memory foam topper, fitted sheet, top sheet and pillowcase inside the open storage bag. Zip the bag shut, applying pressure with your hand to push down the sleep kit.

Note: Store your sleep kit in your Flex™ Storage Ottoman or in your Flex Seat Box.