



New Product Announcement

Oxivir® Tb RTU – 55 gal. / 208 L Drum



In an effort to supply greater amounts of disinfectant to critical areas, Diversey has launched a new Oxivir $^{\circ}$ Tb product in a 55 gal. / 208 L drum.

Technical data	Oxivir [®] Tb RTU
Diversey SKU	101104052
Format	55 gal. / 208 L Drum
Yield	220 - 32 oz. / 946 mL bottles per drum
Dimensions	23.3" x 23.3" x 34.8" (L x W x H)
Weight	485.55 lbs. / 220.24 kg each


Chemical Transfer Options: Due to the unique characteristics of products sold in bulk, operators will need to implement a solution that will allow them to transfer product to smaller containers/bottles. While Diversey cannot provide a transfer solution for every situation, there are simple options that will work for most.

Preferred Option: Top Mounted Manual Pump

The preferred option for dispensing chemical from a 55 gal. / 208 L drum is a top-mounted mechanical pump that uses a suction tube to pull/pump out the solution. These pumps attach to the 2" threaded cap fitting found on the top of the drum. For transferring Accelerated Hydrogen Peroxide® (AHP®) products, only a corrosion resistant (non-metal) pump should be used. Basic tools may be required for installation.

Diversey is working to be able to provide a <u>Tolco 160114 pump</u> for this purpose and will have a limited supply available for purchase though <u>Plan B</u> using a credit card. They can also be sourced elsewhere. Contact your Diversey Sales Representative if you have questions regarding this option.

It is the responsibility of the operator to ensure the safety of users. Users should wear all proper Personal Protective Equipment (PPE), and follow all guidelines for safe handling of chemicals.



Due to high demand, Diversey may not be able to supply pre-printed aftermarket labels. Distributors and End-Users should <u>click here to request self-print .pdf aftermarket label files</u>. Please note that it is unlawful to re-sell any dispensed chemical.

Questions: If you have questions please contact your Diversey sales manager or call Diversey Customer Technical Support (TCS) at (800) 558-2332.



2" Threaded Cap



Tolco 160114 Pump