

LYSOL®

Laundry Sanitizer Additive

Lysol $^\circ$ Laundry Additive kills 99.9% of bacteria like Staphylococcus Aureus & Klebsiella pneumoniae. Simply add to rinse cycle. Contains 0% bleach. Gentle on most fabrics including whites & colors.*

Features & Benefits

- Suitable on most fabric types. Use on towels, bedding, work uniforms, clothing for all ages, athletic wear
- Use in cold or hot water
- Works in standard and HE machines



* When used as directed.







Use instructions

- Dosing: Standard machines: Fill 2 capfuls to line 1 of the dosing cup. HE machines: Fill 2 capfuls to line 2 of the dosing cup.
- Pour into Machine: Add two capfuls to the fabric softener compartment at the start of laundry cycle.
- If machine does not have a fabric softener dispenser, wait until rinse cycle and add directly to drum of machine.
- Leave in rinse cycle for 16 minutes.

To clean surfaces:

- Spray soiled area until thoroughly wet, then wipe clean with clean cloth or sponge.
- Repeat for stubborn stains.
- No rinsing required.

To Sanitize/Disinfect:

- Pre-clean surface.
- Apply to surface until thoroughly wet.

To sanitize:

• Leave for 30 seconds before wiping.

To disinfect:

• Leave for 10 minutes before wiping.

Rinse all food contact surfaces with water after use.

Technical data	LYSOL® Laundry Sanitizer Additive		
Certifications	DIN 02455803		
Color/Form	Blue, Liquid		
рН	9.4 to 10.4		
Shelf Life	3 years		

Product	Pack size	Dilution	Product code	
LYSOL® Laundry Sanitizer Additive - Crisp Linen	4 x 2.7 L Bottle	Ready-to-use	CB986877	I+I

Safe handling

Please make sure your employees read and understand the product label and Safety Data Sheet before using this product. The label contains directions for use; and both the label and SDS contain hazard warnings, precautionary statements and first aid procedures. SDS are available online at www.diversey.com or by calling 888.352.2249. Improper use or dilution may result in damage to surfaces and may result in health and physical hazards that match those of the concentrate.

LYSOL® is a registered trademark of Reckitt Benckiser and is used with their permission.