




DynaDerm™

Hydrocolloid Dressings



BORDERED

ITEM 3014	2" × 2" Thin	20 / bx	12 bxs / cs
ITEM 3015	4" × 4" Thin	10 / bx	12 bxs / cs
ITEM 3016	6" × 6" Thin	5 / bx	12 bxs / cs
ITEM 3017	6" × 7" Sacral	5 / bx	12 bxs / cs

NON-BORDERED

ITEM 3018	4" × 4" X-Thin	10 / bx	12 bxs / cs
ITEM 3019	6" × 6" X-Thin	5 / bx	12 bxs / cs
ITEM 3020	4" × 4" Foam	5 / bx	12 bxs / cs

- Thermally insulated to reliably protect a variety of non-infected wounds from contamination
- Semi-occlusive formula helps to maintain moisture and promote autolytic debridement for optimal healing
- Premium water-resistant materials help prevent bandage peeling for secure protection
- Flexible, conformable dressings with tapered edges provide a comfortable, anatomical fit
- Engineered for safe, long-lasting wear, depending on wound condition and exudate levels
- Caution: Avoid using on diabetic foot wounds without consulting physician for assessment
- Sterile, not made with natural rubber latex



| 1.888 DYNAREX | DYNAREX.COM |

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |

