



# DynaCol AG

## Collagen Wound Dressings



<b>ITEM 3095</b>	<b>2" x 2"</b>	<b>10 /bx</b>	<b>10 bxs / cs</b>
<b>ITEM 3100</b>	<b>4" x 4"</b>	<b>10 /bx</b>	<b>10 bxs / cs</b>
<b>ITEM 3101</b>	<b>6" x 6"</b>	<b>10 /bx</b>	<b>10 bxs / cs</b>
<b>ITEM 3105</b>	<b>7" x 7"</b>	<b>5 /bx</b>	<b>10 bxs / cs</b>

DynaCol AG Collagen Wound Dressings are made from purified bovine collagen with silver chloride, designed to maintain a moist wound environment and to help protect against microbial contamination in partial and full-thickness wounds.

- Bovine Type I collagen wound dressings combined with silver chloride designed for partial and full-thickness wounds
- Soft gel forms upon contact with wound exudate to help maintain a moist healing environment
- Collagen matrix supports new tissue growth and natural healing
- Silver chloride helps protect the dressing from microbial contamination
- Pliable, absorbent material conforms to irregular wound surfaces for effective coverage
- Sterile, single-use dressings individually packaged for hygiene and convenience
- Suitable for use on ulcers, donor sites, and surgical wounds
- Available in multiple sizes to accommodate various wound types



1.888 DYNAREX | DYNAREX.COM |

155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |