

Quick Clean™ Micro-steam™ Bags

Steam clean with confidence

Steam clean pump parts and breastfeeding accessories in the microwave and **eliminate 99.9% of most common bacteria and germs.**

- Works in **90 seconds**
- Re-usable: each bag can be used **20 times**



CDC Guidance



U.S. Department of
Health and Human Services
Centers for Disease
Control and Prevention

Sanitize kit parts for extra protection

*"For extra germ removal, sanitize pump parts at least once daily. Sanitizing is especially important for newborns less than 2 months old, born prematurely, or has a weakened immune system due to illness or medical treatment ."**

Ordering Information

Legacy #	New SAP #	Description	Quantity
N/A	101049153	Quick Clean™ Micro-Steam™ Bags (Bulk Singles)	100/cs
87024NA	101029211	Quick Clean™ Micro-Steam™ Bags (5ct)	5/pk

For more information contact your Medela Sales Consultant.