

G1 SCBA Cleaning Options

Providing you with more options for cleaning



Easy Facepiece Cleaning

- No electronic components are included on the facepiece to worry about during facepiece cleaning
- Easily clean the entire facepiece without removing components

Removable Soft Goods

- The G1 SCBA Soft Goods can be easily removed for enhanced cleaning protocols
- The G1 Soft Goods Kit provides you with a second set of G1 Soft Goods to reduce out-of-service time and extend the life of your SCBA's soft good components



Clean in a SCBA Washing Machine!

- Clean your entire G1 SCBA (including cylinder and facepiece) in a qualified SCBA Washing Machine*
- Solo Rescue and Meiko are qualified manufacturers to clean the G1 SCBA

* Internal testing determined no degradation of the G1 SCBA after 180 wash cycles.

PART NO.	DESCRIPTION
10181021	G1 Soft Goods Kit (serviceable tunnel straps, basic lumbar, 4-pt harness)
10181022	G1 Soft Goods Kit (serviceable tunnel straps, swiveling lumbar, 4-pt harness)
10181023	G1 Soft Goods Kit (serviceable tunnel straps, basic lumbar, 5-pt harness)
10181024	G1 Soft Goods Kit (serviceable tunnel straps, swiveling lumbar, 5-pt harness)
10207584	G1 SCBA Washing Machine Kit