

Enteral Feeding Syringes

Mitigate Risk & Costs in Your NICU



8/16/24 Updated FDA Recommendation:

"Use syringes not manufactured in China, if possible. If you only have syringes manufactured in China, continue to use until you can use alternatives. Closely monitor for leaks, breakage, and other problems." ([Learn more](#))



Section 301 China Tariff Update 9/13/24:

Enteral feeding syringe tariffs increase from 0% to 100% effective 1/1/2026. ([Learn more](#))

Syringes designed exclusively for neonatal needs

Medela Syringes have a multi-functional design and can be used for pump feeds, closed gravity feeding, and venting. Available in ENFit® and TwistLok™ sizes 1-60 mL. Compatible with most syringe pumps. All syringes individually wrapped and sterile.



Bold, highly visible gradient markings helpful in low-lit environments.



Oral tip in 1 mL and 3 mL syringes for precise, gentle medication delivery and oral care.



Gravifeed Syringe Lid for venting with built-in hydrophobic filter. No gauze or hanging necessary.



Medela's ENFit® and TwistLok® syringes are produced in carefully selected locations known for their high-quality manufacturing standards.

Medela LLC
1101 Corporate Drive, McHenry, IL 60050
www.medela.us

© 2025 Medela LLC. Not responsible for errors. Content appearance may vary from pictures. All rights reserved. Medela wordmark and logo are registered in the U.S. Patent and Trademark Office. ENFit® is a federally registered trademark of GEDSA and is used with permission.

medela
THE SCIENCE OF CARE™