

Contact™ Nipple Shields

Support when latching is difficult

Understanding Latch Challenges

Mothers may experience challenges such as difficulty latching, sore, cracked, or flat/inverted nipples.^{1,2}

Provide temporary support for breastfeeding when latching is difficult or painful, while also promoting healing.



Contact™ Nipple Shields

- Unique shape maximizes skin-to-skin contact
- Ultra-thin, flexible silicone maximizes nipple stimulation
- Ideal for flat or inverted nipples
- Pre-term and term infant sizing may vary
- Multiple sizes of sterile and non-sterile versions available

Ordering Information

Legacy #	New SAP #	Description	Quantity
67251S	101028955	Contact™ Nipple Shield - 16 mm (Sterile)	20/cs
67218S	101028958	Contact™ Nipple Shield - 20 mm (Sterile)	20/cs
67203S	101028960	Contact™ Nipple Shield - 24 mm (Sterile)	20/cs
n/a	101034109	Contact™ Nipple Shields, 16 mm, with Storage Case (Retail Packaging, Non-Sterile)	6 pairs/pk
n/a	101034108	Contact™ Nipple Shields, 20mm, with Storage Case (Retail Packaging, Non-Sterile)	6 pairs/pk
n/a	101034107	Contact™ Nipple Shields, 24mm, with Storage Case Retail Packaging, Non-Sterile)	6 pairs/pk

For more information, contact your Medela Sales Consultant

1. Kronborg H et al. Matern Child Nutr. 2017; 13(1). 2.Coentro VS et al. J Obstet Gynecol Neonatal Nurs. 2022; 51(1):73-82.