


Resp-CO₂

Equal Length EtCO₂ Sampling Line & O₂ Tubing

EtCO₂ NASAL SAMPLING CANNULAS

Dynarex Resp-O₂™ EtCO₂ Nasal Sampling Cannulas with Cushion Nasal Tips are designed to help monitor exhaled carbon dioxide during respiratory treatment. Each soft, cushioned cannula features a permanent barrier between the nasal prongs so that EtCO₂ can be sampled from one nostril while oxygen continues to flow to the other nostril.

- Permanent barrier between nasal prongs allows for EtCO₂ sampling with uninterrupted oxygen delivery
- Constructed from soft materials for added patient comfort
- Crush-resistant tubing helps prevent occlusion and kinking
- Features EtCO₂ sampling line and O₂ tubing of equal length
- Transparent design allows for clear visualization
- Not made with natural rubber latex or DEHP



<p>7' (2.1 m) Tubing</p> <p>MALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33180</p> <p>INFANT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>7' (2.1 m) Tubing</p> <p>MALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33181</p> <p>PEDIATRIC</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>7' (2.1 m) Tubing</p> <p>MALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33182</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>7' (2.1 m) Tubing</p> <p>FEMALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33183</p> <p>INFANT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>7' (2.1 m) Tubing</p> <p>FEMALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33184</p> <p>PEDIATRIC</p> <p>PACKAGE QUANTITY 25/Case</p>
--	---	---	--	---

<p>7' (2.1 m) Tubing</p> <p>FEMALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33185</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>10' (3.0 m) Tubing</p> <p>MALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33186</p> <p>PEDIATRIC</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>10' (3.0 m) Tubing</p> <p>MALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33187</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>14' (4.2 m) Tubing</p> <p>MALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33188</p> <p>PEDIATRIC</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>14' (4.2 m) Tubing</p> <p>MALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33189</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>10' (3.0 m) Tubing</p> <p>FEMALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33190</p> <p>PEDIATRIC</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>10' (3.0 m) Tubing</p> <p>FEMALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33191</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>14' (4.2 m) Tubing</p> <p>FEMALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33192</p> <p>PEDIATRIC</p> <p>PACKAGE QUANTITY 25/Case</p>
---	--	--	--	--	--	--	--

<p>14' (4.2 m) Tubing</p> <p>FEMALE LUER CONNECTOR STANDARD O₂ TUBING CONNECTOR</p> <p>ITEM 33193</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>7' (2.1 m) Tubing</p> <p>MALE LUER CONNECTOR UNIVERSAL O₂ TUBING CONNECTOR</p> <p>ITEM 33194</p> <p>PEDIATRIC</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>7' (2.1 m) Tubing</p> <p>MALE LUER CONNECTOR UNIVERSAL O₂ TUBING CONNECTOR</p> <p>ITEM 33195</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>7' (2.1 m) Tubing</p> <p>FEMALE LUER CONNECTOR UNIVERSAL O₂ TUBING CONNECTOR</p> <p>ITEM 33196</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>10' (3.0 m) Tubing</p> <p>FEMALE LUER CONNECTOR UNIVERSAL O₂ TUBING CONNECTOR</p> <p>ITEM 33197</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>10' (3.0 m) Tubing</p> <p>MALE LUER CONNECTOR UNIVERSAL O₂ TUBING CONNECTOR</p> <p>ITEM 33198</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>14' (4.2 m) Tubing</p> <p>MALE LUER CONNECTOR UNIVERSAL O₂ TUBING CONNECTOR</p> <p>ITEM 33199-M</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>	<p>14' (4.2 m) Tubing</p> <p>FEMALE LUER CONNECTOR UNIVERSAL O₂ TUBING CONNECTOR</p> <p>ITEM 33199-F</p> <p>ADULT</p> <p>PACKAGE QUANTITY 25/Case</p>
--	--	--	--	---	---	---	---

