



# DynaGinateAG™

## Calcium Alginate Dressings



<b>ITEM 3030</b>	<b>2" × 2"</b>	<b>10 / bx</b>	<b>10 bxs / cs</b>
<b>ITEM 3031</b>	<b>4" × 5"</b>	<b>10 / bx</b>	<b>10 bxs / cs</b>
<b>ITEM 3032</b>	<b>4" × 8"</b>	<b>5 / bx</b>	<b>10 bxs / cs</b>
<b>ITEM 3033</b>	<b>¼" × 12" (Rope)</b>	<b>5 / bx</b>	<b>10 bxs / cs</b>

Calcium Alginate Dressings are formulated to conform to wound contours, transforming into a gel-like covering upon contact with moisture, promoting an optimal healing environment with an antimicrobial barrier.

- Formulated with a combination of alginate and silver ions that change to a fibrous gel when in contact with wound exudate
- Creates a moist wound environment and provides an antimicrobial barrier
- Suitable for infected wounds or where there is a risk of infection
- High absorbency sterile dressing used for heavily exuding wounds
- Primary dressing ideal for the management of partial and full thickness wounds at stages 3 or 4



| 1.888 DYNAREX | DYNAREX.COM |

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |