

5L OXYGEN CONCENTRATOR

REORDER NO. 33950

CONTINUOUS OXYGEN FLOW 5 LPM

OXYGEN CONCENTRATION 90 to 96%

OXYGEN OUTLET PRESSURE 6-10 Max PSI

OPERATING ALTITUDE

(FEET ABOVE SEA LEVEL) 0 to 7,500 or 13,123 as a setting

DECIBEL LEVEL 40

WEIGHT 37 lb.

37 ID.

PACKAGE QUANTITY 1/Case The Dynarex Resp-O₂[™] Oxygen Concentrator is designed to provide supplemental oxygen to help relieve breathing difficulties in respiratory patients. The 5-liter Oxygen Concentrator includes audible alarms, oxygen purity indicator (OPI), and user-friendly flow-meter knob to provide prescribed liters-per-minute. The 37-pound concentrator is easily transportable with accessible top handle and four rolling casters.

- Features a flow range of 0.5-5 LPM
- Audible and visual alarms include Low Flow-Rate, Low/High Pressure, and Power Failure
- Manual reset button helps ensure uninterrupted oxygen delivery
- Designed for quiet operation with low noise output
- Easy-to-read control panel for optimal visibility

ALSO AVAILABLE:

Dry Bubble Humidifier - 400 ML · 6 PSI (HOME USE)

REORDER NO. 34221

Oxygen Humidifier Elbow Connectors with 12" Tubing

REORDER NO. 34505

Oxygen Humidifier Straight Connectors with 18" Tubing

REORDER NO. 34506

dynarex

Sterile Water for Inhalation, USP - PREFILLED BUBBLE HUMIDIFIER

REORDER NO. 34300 – SIZE 350 ML

REORDER NO. 34301 – SIZE 550 ML





R240905



Sterile Water for Inhalation

Dynarex Resp-O₂[™] Sterile Water supplies are disposable, singleuse respiratory therapy bottles filled with water designed for patient inhalation. The water is sterilized to ensure patient safety and to help avoid cross-contamination.

- Reservoir bottles feature built-in oxygen tubing connectors
- Sterile, transparent adaptor body provides optimal visualization
- Audible alarm alerts clinician to air-flow constrictions

OXYGEN CONCENTRATOR ACCESSORIES

Bubble Humidifier


The Dynarex Resp-O₂[™] Disposable Bubble Humidifier is designed to humidify supplemental oxygen for patients receiving oxygen therapy at home. Ideal for use with all types of oxygen concentrators with either face masks or cannulas, the Bubble Humidifier is intended for oxygen low flow rates up to 6 PSI. The transparent humidifier features clearly defined minimum and maximum lines for easy monitoring.


- Features pressure-relief valve with audible alarm for optimal safety
- Includes PVC diffuser to help minimize noise
- Disposable, single-patient use for optimal hygiene
- Not made with natural rubber latex or DEHP

Oxygen Humidifier Connectors

Dynarex Resp-O2[™] Oxygen Humidifier Connectors are designed to attach dry humidifier bottles or pre-filled bubble humidifiers to continuous-flow oxygen concentrators. They come with either sturdy oxygen tubing on one end and a male threaded grip elbow on the other end, or with star lumen oxygen tubing on one end and molded threaded straight connector on the other end. Connectors are universal and fit most humidifier bottles with a DISS connection.

- Transparent oxygen tubing allows for easy visualization
- Universal fit connects with most humidifier bottles and concentrators
- Ribbed connector provides a secure fit without damaging connection points
- Designed to attach to barb-style concentrator outlets
- Not made with natural rubber latex or DEHP









R240905