



# Xeroform™

## Petrolatum Gauze Dressings



<b>ITEM 3042</b>	<b>3" × 18"</b>	12 / bx	6 bxs / cs
<b>ITEM 3053</b>	<b>4" × 4"</b>	25 / bx	4 bxs / cs
<b>ITEM 3054</b>	<b>5" × 9"</b>	50 / bx	4 bxs / cs

A non-adherent knitted gauze impregnated with a blend of petroleum, 3% Bismuth, mineral, and paraffin oils. These highly conformable, non-occlusive dressings allow the wounds to breathe and are indicated for use in deep to shallow wounds with minimal exudate.

- Helps maintain a moist wound environment
- Non-adherent and non-occlusive
- Conformable for use with deep and shallow wounds
- Knitted gauze, impregnated wound dressing
- Designed to allow the wound to breathe, minimizing skin and tissue maceration
- Sterile primary dressing for low-exuding wounds at stages 2, 3, or 4
- Active ingredient 3% Bismuth reduces wound odor



| 1.888 DYNAREX | DYNAREX.COM |

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |